


Go!

Temp

**Collect
temperature
data on your
computer!**



Quick-Start Guide

Get Going in 3 Easy Steps!

1. Install Logger Lite software.
2. Plug in your Go!Temp.
3. Collect data!

For detailed installation instructions, open the Installation Guide found on the enclosed CD.

Tutorials

A great way to get started using your Go!Temp is to work through the tutorials. Simply start Logger Lite, choose Open from the File menu, then open the Tutorials folder. Open Tutorial 1 and follow the on-screen instructions.

Logger Lite Toolbar Buttons



Create a new file.



Autoscale the graph.



Open a file from *Let's Go! Investigating Temperature* or one of the Tutorials.



Examine your data point by point.



Save your file.



Calculate statistics on your data.



Print the screen.
More print options are found in the File menu.



Predict what will happen before you collect.



Want to collect a second run?
Store them as you go.



Switch between °C and °F.



Switch between pages.



Choose from three styles of meters!



Start collecting data!






Double-click the column heading to edit

Double-click to add graph title

Digital meter provides live readouts



Thermometer can show °C and °F

Grab axes to stretch

Go!Temp Specifications

Range: -20°C - 110°C

Maximum temperature that the sensor can tolerate without damage: 130°C

Resolution: 0.07°C

Accuracy: +/- 0.5°C


Response time: 4s (to 90% of full reading in water)

Go!™

Logger Lite™
version 1.2



Go!



Go!Temp

Go!Temp is an easy-to-use temperature probe that connects directly to your computer's USB port.

Using **Logger Lite software** (included), you can collect, analyze, print, and save your Go!Temp data right on your computer.

Logger Lite offers simple functionality for those new to data collection. If you are looking for more advanced features for data collection and analysis with your Go!Temp, visit www.vernier.com/loggerpro to learn more about our award-winning software, *Logger Pro*®.




New!

Use Go!Temp in so many ways!

- Science fair projects
- Monitor the temperature of your fish tank
- Monitor the temperature of your greenhouse
- Environmental monitoring
- Check your refrigerator and freezer temperatures
- Check the weather
- Monitor your hot tub's temperature
- Measure the efficiency of a thermos
- Measure the temperature of the air coming from your air conditioner or furnace
- Monitor the soil temperature of your outside plants

Windows system requirements: Windows 98, 2000, Me, or XP on a Pentium processor or equivalent
Macintosh system requirements: Mac OS X (10.2 or newer)

Vernier Software & Technology
13979 S.W. Millikan Way • Beaverton, OR 97005-2886
Toll free 1-888-837-6437 • FAX 1-503-277-2440
info@vernier.com • www.vernier.com



Vernier Software & Technology is the leader in educational data-collection technology. For other innovative products from Vernier, such as Go! Link, Vernier LabPro, and *Logger Pro*, visit our web site at www.vernier.com.



Go! and Logger Lite are trademarks and Vernier LabPro and Logger Pro are registered trademarks of Vernier Software & Technology in the United States of America. Windows is a registered trademark of Microsoft Corporation. Macintosh is a registered trademark of Apple Computer, Inc.

Note: This product is not appropriate for industrial, medical, research, or commercial applications.