



Engage the Scientists
of Tomorrow



Vernier Software started in 1981 as a two-person, part-time company run by a physics teacher and a social worker. Today we have 110 employees, and we sell in over 150 countries. We are proud to celebrate our 40th year in 2021!

Last year in this catalog letter, we made a comment about “living in interesting times.” We were just talking about the complications of tariffs on our pricing. Now, in 2021, we are all really living through interesting times! With all the disruption of our business, we are pleased to say that we have been able to retain all of our employees (while working mostly from home). We have changed the way we do a lot of things, but like you, we are adapting.

When schools shut down in March, we quickly posted free experiment data so instructors would have something to share with students as they improvised ways to teach science remotely. And we literally pivoted—with Pivot Interactives—one of the best and most popular products for remote learning.

We now have a great collection of software tools for teaching science remotely, and they all work on Chromebooks, as well as computers and tablets. Our new Vernier Graphical Analysis Pro app includes data from many of the experiments in our lab books with videos taken of the procedure during data collection. With the Vernier Video Analysis app, students can take videos with their cell phones and analyze their motion data. If you have not tried out these programs, please do. All are available for a free 30-day trial.

And as we do every year, we have introduced some new Go Direct sensors. This year we added the Go Direct Weather System, Go Direct Thermocouple, Go Direct Static Charge, and Go Direct Platinum-Cell Conductivity.

And finally, we are excited to celebrate our 40th year with the introduction of LabQuest 3! It is a major upgrade to our LabQuest line of handheld data-collection tools with a large screen and advanced touch-screen abilities.

We know school budgets will be affected by COVID-19 this year, and we have decided not to raise prices on Vernier products. We encourage you to give our products a try on a 30-day (or longer) preview basis. Feel free to contact any of us personally at any time.

Stay positive and test negative!

John Wheeler
CEO
jwheeler@vernier.com

Dave and Christine Vernier
Co-Presidents
dvernier@vernier.com and cvernier@vernier.com



About Vernier Software & Technology

Vernier Software & Technology was co-founded in 1981 by Dave and Christine Vernier. Dave's background as a physics teacher and Christine's knack for business combined to form a company with a deep commitment to education.

Forty years later, the company is still owned by Christine and Dave, along with nine employee owners who have backgrounds in science and math education, as well as business.

Vernier is proud to be recognized for its philanthropic commitment, environmental policies, steady growth, and as one of the Best 100 Companies to Work For in Oregon for 20 years.



2020 Best Companies to Work For in Oregon





2020 Healthiest Employers of Oregon



2020 Best Green Companies in Oregon



2020 Corporate Philanthropy Award



On the Cover

Monitoring ecosystem
abiotic factors

Why Vernier?

Instill a Love of Learning in All Students

Your passion and dedication, along with the implementation of high-quality sensors, experiments, and resources in your classroom, enable your students to explore science in new ways.

Our mission is to provide you with the tools you need to encourage scientific curiosity in all students—see what partnering with us can do.



10 Tips for Writing Your Best Grant Proposal

We understand that grants are essential for you to get the supplies, tools, and resources necessary to address the many needs of your students.

This year, with school budgets in such a precarious place and remote learning still playing such a large role, securing grant funding means you and your students can have the support needed to thrive, no matter where learning takes place.

We have created an infographic with 10 tips for grant writing to help you perfect your proposal with newfound confidence.

www.vernier.com/grants

Contents

What's New

PAGE 2

Elementary School

PAGE 4

Middle School

PAGE 16

High School

GETTING STARTED PAGE 30

BIOLOGY PAGE 42

ENVIRONMENTAL SCIENCE PAGE 58

EARTH SCIENCE PAGE 70

CHEMISTRY PAGE 74

PHYSICAL SCIENCE PAGE 92

PHYSICS PAGE 96

ENGINEERING, CODING, AND ROBOTICS
PAGE 122

TEXAS INSTRUMENTS PAGE 132

COLLEGE www.vernier.com/college



Sensors & Accessories

PAGE 134

Index

PAGE 138

What's New?



Vernier Graphical Analysis Pro

We are enhancing our award-winning Vernier Graphical Analysis™ app with advanced features supporting remote learning and more advanced analysis of experiment data.

Learn more at www.vernier.com/graphical-analysis-pro

LabQuest 3



LabQuest 3 is a powerful, advanced, easy-to-navigate, and versatile data-logging solution for STEM students.

The all-new LabQuest® 3 reimagines data collection by providing students with an innovative, easy-to-use interface. A larger screen and advanced touch-screen abilities make it easier for students to collect, graph, and analyze data wherever they are—the classroom, at home, or in the field. Challenge your students to gain a deeper understanding of science through data with the accessible, groundbreaking LabQuest 3.

Learn more on pp. 31–33.

Remote Learning



Keep students engaged in STEM with our remote learning solutions, including remote alternatives to hands-on experiments, coding activities, and more.


Learn more at www.vernier.com/remote-learning

Investigations



Food Chemistry Experiments

This new lab book is filled with experiments that use food as a means to explore crucial chemistry concepts. Students are more likely to engage with science when they see concepts applied to the real world. Learn more on page 83.



Climate and Meteorology Experiments

Challenge students to use data-collection technology to explore storm systems and other important weather-related topics. Learn more on page 23.



Vernier Video Analysis: Motion and Sports

Expand students' learning opportunities and further connect the study of motion to their daily lives with these investigations using Vernier Video Analysis™. Learn more on page 119.



Sensor Cart Physics

Explore introductory AP*-level concepts in kinematics, dynamics, and conservation of energy and momentum using the Go Direct® Sensor Cart. Learn more on page 102.



Vernier Coding Activities with Arduino®: Analog Sensors

Integrate Vernier sensor technology with Arduino and connect the physical world to the computer-centric activity of learning to code. Learn more on page 126.



Human Physiology Experiments: Volume 2

This lab book contains 15 experiments designed to encourage students to explore the physiology of various human organ systems. Learn more on page 50.



OpenSciEd

Our partnership with OpenSciEd gives you access to free, field-tested units that support the three-dimensional learning approach. Learn more on pp. 19–21.

Sensors



Go Direct Weather System

Easily monitor a wide variety of environmental factors with just one sensor. Learn more on page 51.



Go Direct Static Charge

Unlike most electroscopes, this new sensor provides digital measurements of positive and negative charges of objects that would not be possible in a traditional lab. Learn more on page 108.



Go Direct Thermocouple

With this affordable and easy-to-use sensor, students can collect reliable data during experiments that involve extreme temperatures. Learn more on page 85.

Elementary School

www.vernier.com/elementary-school

Why Vernier?

Technology engages young students. Our carefully designed hands-on data-collection technology helps elementary school teachers introduce young learners to science and STEM. We've created easy-to-use resources to help you educate and inspire your students.

EASY

Simple for students and teachers to use

AFFORDABLE

Priced to fit school budgets

VERSATILE

Compatible with a variety of devices



I can't even imagine all of the amazing things I'll be able to do with the kids with your products. I'm just beyond grateful for companies like yours who give back and help teachers inspire tomorrow's science leaders.

*Covey Denton,
Greenfield School*



Topics

Explore a sampling of our featured experiments by topic to learn how Vernier technology helps your students deepen their understanding of key STEM concepts.

Temperature

PAGE 7

Gas Pressure

PAGE 8

Motion

PAGE 8

Force

PAGE 9

Light

PAGE 9

Magnetism

PAGE 10

Voltage

PAGE 10

Wind Energy

PAGE 12

Solar Energy

PAGE 12

Robotics

PAGE 13

Coding

PAGE 13



Instill a Lifelong Love of Learning

Young minds are naturally curious; engage your students with fun, interactive lessons that encourage investigation of their world and instill a lifelong love of learning.



New Lessons? They're Now a Breeze

From bubbling bread and baking soda reactions to reflectivity of light and simple motion, we offer a variety of student-ready, easy-to-implement investigations designed to help excite and engage your young learners.




Educational Standards

Helping students meet standards is an important aspect of teaching. Vernier technology helps teachers as they prepare students to meet the NGSS and state standards through investigations that support three-dimensional learning.

www.vernier.com/standards

Getting Started



What You Need to Get Started

A Go Direct Sensor

These versatile sensors connect to your device via Bluetooth® wireless technology or USB.

B Device

Go Direct® sensors connect to a wide variety of devices commonly used in classrooms, including Chromebooks, computers, compatible mobile devices, and LabQuest® 3.

C Vernier Graphical Analysis App

Our data-collection app facilitates student understanding with real-time graphs of experimental data.

D Lab Book

Step-by-step instructions at your fingertips save valuable time when integrating probeware into your curriculum. Most of our lab books for elementary school provide support for Go Direct sensors and the Graphical Analysis™ app.

Our lab books come with a generous site license—purchase once and share files school wide.

Next Generation Science Standards

Hands-on learning has been at the core of Vernier's mission for 40 years, and as we create new products—whether it is hardware, software, or written investigations—we work to align to the NGSS, making it easy for teachers and science supervisors to help students meet these standards.

Vernier Book	NGSS DCI Topics			
	Physical Science	Life Science	Earth and Space Science	Engineering Design
<i>Investigating Temperature</i>	●			●
<i>Investigating Gas Pressure</i>	●	●		
<i>Investigating Motion</i>	●	●		
<i>Investigating Force</i>	●			
<i>Investigating Light</i>	●		●	
<i>Investigating Magnetism</i>	●			
<i>Investigating Voltage</i>	●			
<i>Elementary Science with Vernier</i>	●	●	●	●
<i>Investigating Wind Energy</i>	●			●
<i>Investigating Solar Energy</i>	●			●
<i>Coding with Codey Rocky™: Mission to Mars</i>	●		●	●

Temperature

Investigating Temperature

Watch a video



Download only
ELB-TEMP-E \$20

Download + print
ELB-TEMP \$25

In this book, students investigate topics related to temperature, including melting and freezing of water, insulation design, and chemical reactions.

10 Experiments Included

Physical Science

STRUCTURE AND PROPERTIES OF MATTER

- I'm Melting! Water Changes States
- Solid, Liquid, Gas: Water Can Do It All

ENERGY

- Are We Cool or What?
- Why Do We Need Thermometers?
- Celsius or Fahrenheit: What's the Difference?

- Getting it Just Right! Adjusting Water Temperature
- The Temperature Probe Spends the Night
- Hold Everything! Comparing Insulators
- Keeping it Cool! Design Your Own Thermos
- Cool Reaction! The Reaction of Baking Soda and Vinegar (*shown above*)

Sensor Used

Go Direct Temperature

Students use this rugged, general-purpose sensor to monitor temperature.

GDY-TMP \$69

Teacher pack also available
(includes 8 Go Direct Temperature Probes and a Charge Station)
GDY-TMP-TP \$599



Learn more at www.vernier.com/elb-temp

Gas Pressure

Motion

Investigating Gas Pressure



Download only
ELB-GP-E \$10

Students investigate the behavior of gas pressure when more gas is added or the volume of the container changes.

4 Experiments Included in E-book

- Learning to Use a Pressure Sensor

STRUCTURE, FUNCTION, AND
INFORMATION PROCESSING

Life Science

MATTER AND ENERGY IN ORGANISMS
AND ECOSYSTEMS

- Bubbles in Your Bread

- Get a Grip! (*shown above*)

Physical Science

FORCES AND INTERACTIONS

- Under Pressure

Products Used



Go Direct® Gas Pressure

Measure the change in gas pressure as variables such as temperature and volume change.

GDX-GP \$89



Gas Pressure Sensor Bulb

GPS-BULB1 \$6

Learn more at www.vernier.com/elb-gp-e

Investigating Motion



Download only
ELB-MD-E \$10

The motion of a bouncing ball and a toy car are just two examples of the investigations about motion that students conduct using this e-book.

7 Experiments Included in E-book

- Learning to Use a Motion Detector

ENERGY

Physical Science

FORCES AND INTERACTIONS

- e-Motion!
- Spring into Action
- Air Ball! (*shown above*) also uses Go Direct Gas Pressure.

- Driving with Energy

- Weigh Station—All Trucks Stop!

Life Science

STRUCTURE, FUNCTION, AND
INFORMATION PROCESSING

- Batty About Science

Sensor Used

Go Direct Motion

Monitor the position of a moving object using ultrasound.

GDX-MD \$99



Learn more at www.vernier.com/elb-md-e

Force

Light

Investigating Force



Download only
ELB-FOR-E \$10

Everyday forces, such as the frictional force on a shoe, are investigated in this e-book.

4 Experiments Included in E-book

- Learning to Use a Force Sensor

Physical Science

FORCES AND INTERACTIONS

- Lift the Load!
- What a Drag! (*shown above*)
- Oh! My Aching Back! How Ramps Make Lifting Easier

Sensor Used

Go Direct Force and Acceleration

Use this force sensor to measure the force of pushes and pulls in the classroom and outdoors. This sensor can also measure acceleration.

GDX-FOR \$99



Learn more at www.vernier.com/elb-for-e

Investigating Light



Download only
ELB-LC-E \$10

Students investigate light properties including how light changes with distance, reflects off different colors, and varies with the seasons.

5 Experiments Included in E-book

- Learning to Use a Light Sensor

Physical Science

WAVES: LIGHT AND SOUND

- Sunshine on My Shoulders

Earth and Space Science

EARTH'S SYSTEMS

- Summer and Winter
- Reflectivity of Light (*shown above*)

SPACE SYSTEMS: STARS AND THE SOLAR SYSTEM

- Distance From the Sun

Sensor Used

Go Direct Light and Color

Students use this sensor to measure the brightness of a light bulb or the reflectance of light off of various objects. They can also measure UV light and relative amounts of red, blue, and green light.

GDX-LC \$79




Learn more at www.vernier.com/elb-lc-e

Magnetism

Voltage

Investigating Magnetism



Download only
ELB-3MG-E \$10

In this e-book, students investigate the magnetic field of magnets and electromagnets.

4 Experiments Included in E-book

- Learning to Use a Magnetic Field Sensor

Physical Science

FORCES AND INTERACTIONS

- Exploring the Poles (*shown above*)
- Making Magnets
- Electromagnets

Sensor Used

Go Direct® 3-Axis Magnetic Field


Use this sensor to explore properties of magnets, electromagnets, and the Earth's magnetic field.

GDX-3MG \$69



Learn more at www.vernier.com/elb-3mg-e

Investigating Voltage



Download only
ELB-VOLT-E \$10

Do C-cell batteries provide a higher voltage than AA batteries? Students investigate this type of question in this e-book focused on voltage.

4 Experiments Included in E-book

- Learning to Use a Voltage Probe

Physical Science

ENERGY

- Are All Batteries the Same? (*shown above*)
- Stacked Batteries
- All Worn Out

Sensor Used

Go Direct Voltage

This sensor is an excellent choice for investigating batteries, circuits, and electromagnets.

GDX-VOLT \$69



Learn more at www.vernier.com/elb-volt-e


Elementary Science with Vernier



This collection of experiments for elementary students includes the topics of temperature, motion, force, magnetism, light, electricity, and gas pressure.

Includes Experiments from These E-books

- Investigating Temperature
- Investigating Gas Pressure
- Investigating Motion
- Investigating Force
- Investigating Light
- Investigating Magnetism
- Investigating Voltage



INCLUDES
43
EXPERIMENTS

Download only

EWV-E \$40

Printed book + download

EWV \$48

Learn more at www.vernier.com/ewv

Elementary Go Direct Package

8 Products • GDP-EL-DX • \$579

Buy 8 or more packages at \$562 and save \$136



This package includes

Go Direct
Temperature

Go Direct Light
and Color

Go Direct
Motion

Go Direct
3-Axis
Magnetic Field

Go Direct
Gas Pressure

Go Direct
Voltage

Go Direct
Force and
Acceleration

Gas Pressure
Sensor Bulb

All sensors work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3.

Learn more at www.vernier.com/gdp-el-dx

Wind Energy

Investigating Wind Energy



Download only
ELB-WIND-E \$20

Download + print
ELB-WIND \$25

Students investigate wind energy to learn about energy transfer, basic electric circuits, and blade design.

11 Experiments Included

- Introduction to Wind Turbines
- Exploring Wind Energy
- Introduction to the Energy Sensor
- Wind Turbine Output: The Effect of Load (*shown above*)
- Exploring Wind Turbine Blades
- Blade Design: Pitch
- Blade Design: Area
- Blade Design: Quantity
- Blade Design: Mass
- Blade Design: Material
- Project: Power Up! (Engineering Design)

Package Available

Investigating Wind Energy Package

Contains the following products

- Go Direct® Energy
- Vernier Resistor Board
- KidWind MINI Wind Turbine with Blade Design



GDP-EL-WE \$172
Buy 8 or more at \$167
and save \$40

Learn more at www.vernier.com/elb-wind

Solar Energy

Investigating Solar Energy



Download only
ELB-SOLAR-E \$20

Download + print
ELB-SOLAR \$25

Solar energy provides a real-world example where students investigate energy transfer, series and parallel circuits, and other factors that affect solar panel output.

11 Experiments Included

- Introduction to Solar Panels
- Exploring Solar Energy
- Introduction to the Energy Sensor
- Making Connections: Circuits
- Solar Panel Output: Effect of Load
- Solar Panel Output: Effect of Shade
- Solar Panel Output: Effect of Angle (*shown above*)
- Pumping Water with Solar Energy
- Exploring Surface Temperature
- Project: Solar Homes (Engineering Design)
- Project: What's Cookin'? (Engineering Design)

Package Available

Investigating Solar Energy Package

Contains the following products

- Go Direct Energy
- Go Direct Surface Temperature
- Solar Energy Exploration Kit
- Vernier Resistor Board

GDP-EL-SE \$265
Buy 8 or more at \$257 and save \$64




Learn more at www.vernier.com/elb-solar

Robotics

Coding

Coding with Codey Rocky: Mission to Mars



Download only
MBCR-M2M-E
\$20*

*Free with purchase
of Codey Rocky
from Vernier

Students program their Codey Rocky robot to explore, learn about, and survive on Mars.

6 Experiments Included in E-book

- Houston, This is Codey
- Dance of the Martians
- Surviving the Desert of Mars
- Wild, Wild Mars
- Daily Life on Mars
- Surveying Mars

Product Used

Codey Rocky™ by Makeblock®

Easy-to-use robotic hardware combined with block-based programming provides students with the ideal introduction to coding.

MB-CR \$104.99



Learn more at www.vernier.com/mbcr-m2m-e

Coding with Scratch



Integrate Go Direct Force and Acceleration into your classroom activities with Scratch. Your students can learn coding by applying their skills to fun, collaborative, hands-on coding projects.

We've designed a free module of Vernier Scratch activities—including a teacher's guide—that helps students sharpen coding skills and gain valuable experience with data-collection technology.

Example Projects

- Storytelling in Scratch: Use block-based coding to tell the story of Newton's "year of wonders."
- Interactive Art: Write code in Scratch to create a parallax effect.
- Ideal Gas Laws: Combine coding and an exploration of the ideal gas laws.

Product Used

Go Direct Force and Acceleration

With Go Direct Force and Acceleration, your students can make a sprite move in response to spinning, tilting, falling, or applying a force to the sensor.

GDX-FOR \$99






Learn more at www.vernier.com/scratch

Featured Products

Go Direct Sensors

Sensor		Order Code	Price
Go Direct® 3-Axis Magnetic Field		GDX-3MG	\$69
Go Direct Energy		GDX-NRG	\$89
Go Direct Force and Acceleration		GDX-FOR	\$99
Go Direct Gas Pressure		GDX-GP	\$89
Go Direct Light and Color		GDX-LC	\$79
Go Direct Motion		GDX-MD	\$99

Go Direct Sound		GDX-SND	\$89
Go Direct Surface Temperature		GDX-ST	\$79
Go Direct Temperature		GDX-TMP	\$69
Go Direct Voltage		GDX-VOLT	\$69
Go Direct Weather		GDX-WTHR	\$99

Go Direct Charge Station







Accessory		Order Code	Price
Go Direct Charge Station		GDX-CRG	\$69

See all our products for elementary school science at www.vernier.com/elementary-school

Additional Products

Product		Order Code	Price
Davis® Weather Stations		www.vernier.com/weather	
Gas Pressure Sensor Bulb		GPS-BULB1	\$6
KidWind MINI Wind Turbine with Blade Design		KW-MWTBD	\$65
Solar Energy Exploration Kit		KW-SEEK	\$79
USB Digital Microscope		BD-EDU-100	\$119
Vernier Resistor Board		VES-RB	\$18

Coding and Robotics

Product		Order Code	Price
Codey Rocky™		MB-CR	\$104.99
mBot™ (blue)		MBOT-B	\$73.49
mBot (pink)		MBOT-P	\$73.49
mBot Explorer		MBOT-S	\$83.99
Neuron Inventor Kit by Makeblock®		MB-NEURON	\$125.99
Go Direct Force and Acceleration (for use with Scratch)		GDX-FOR	\$99

Lab Books

Title	Order Code	Price
<i>Elementary Science with Vernier</i>	Download only: EWW-E Download + print: EWW	\$40 \$48
<i>Investigating Temperature*</i>	Download only: ELB-TEMP-E Download + print: ELB-TEMP	\$20 \$25
<i>Investigating Motion*</i>	Download only: ELB-MD-E	\$10
<i>Investigating Light*</i>	Download only: ELB-LC-E	\$10
<i>Investigating Magnetism*</i>	Download only: ELB-3MG-E	\$10
<i>Investigating Gas Pressure*</i>	Download only: ELB-GP-E	\$10
<i>Investigating Force*</i>	Download only: ELB-FOR-E	\$10
<i>Investigating Voltage*</i>	Download only: ELB-VOLT-E	\$10
<i>Investigating Solar Energy</i>	Download only: ELB-SOLAR-E Download + print: ELB-SOLAR	\$20 \$25
<i>Investigating Wind Energy</i>	Download only: ELB-WIND-E Download + print: ELB-WIND	\$20 \$25
<i>Coding with Codey Rocky: Mission to Mars</i> (Included with purchase of Codey Rocky from Vernier)	Download only: MBCR-M2M-E	\$20
<i>Coding with mBot: Self-Driving Vehicles</i> (Included with purchase of mBot from Vernier)	Download only: MBOT-MSDV-E	\$20

* All experiments from this e-book are included in Elementary Science with Vernier.

See all our products for elementary school science at www.vernier.com/elementary-school

Middle School

www.vernier.com/middle-school

Why Vernier?

Hands-on learning with technology is ideal for middle school students. Enhance their discovery and understanding of the world around them with the use of Vernier technology. Using our versatile, cutting-edge products and ready-to-go experiments correlated to the NGSS and state standards, you can encourage your students' curiosity and prepare them for high school—and the world beyond.

EASY

Simple for students and teachers to use

AFFORDABLE

Priced to fit school budgets

VERSATILE

Supports a variety of devices and investigations



The technology's ease of use and accessibility allows students to really take charge of the learning process as they acquire data; the technology has been a game changer.

*Susan Foster,
Manlius Pebble Hill School*



Contents

Explore our offerings for middle school and learn how Vernier technology helps your students deepen their understanding of key STEM concepts.

Getting Started


PAGE 18

Three-Dimensional Learning Approach

PAGE 19

Classic Approach

PAGE 19



Next Generation Science Standards

Hands-on learning has been at the core of our mission for over 40 years, and as we create new products—whether it is hardware, software, or written investigations—we work to align them to the NGSS, making it easy for you to help students meet these standards.

www.vernier.com/ngss-correlations



Coding and Robotics

Set up your middle school students for success with cutting-edge products and partnerships that encourage curiosity, develop computational thinking skills, and enhance their understanding of the world around them.

www.vernier.com/middle-school/engineering




Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training

Getting Started



What You Need to Get Started

A Go Direct® Sensor

These versatile sensors connect to your device via Bluetooth® wireless technology or USB.

B Device

Go Direct sensors connect to a wide variety of devices commonly used in classrooms, including Chromebooks, computers, compatible mobile devices, and LabQuest® 3.

C Vernier Graphical Analysis™ App

Our free data-collection app facilitates student understanding with real-time graphs of experimental data.

The new Graphical Analysis Pro offers additional exclusive features, such as the ability to perform live experiments and share the data over the internet in real time. Sign up for a free 30-day trial today!

Learn more at www.vernier.com/graphical-analysis-pro

D Lab Book

Step-by-step instructions at your fingertips save valuable time when integrating probeware into your curriculum. Most of our lab books for middle school provide support for Go Direct sensors and the Graphical Analysis app.

Our lab books come with a generous site license—purchase once and share files school wide.

Three-Dimensional Learning Approach



Vernier and OpenSciEd

Vernier knows that science education is not static. Your students need to understand critical scientific concepts, use these concepts to solve problems, and understand how they connect to the real world. These objectives are incorporated into the main pillars of the three-dimensional learning framework developed by the National Research Council. Vernier provides downloadable e-books, shown on the next page, that incorporate the three-dimensional learning approach.

We are proud to partner with OpenSciEd, a provider of high-quality, open-source science instructional materials. Our partnership gives you access to free, field-tested and EQuIP-approved units that support the three-dimensional learning approach. Vernier provides free downloadable supplements that integrate data-collection technology into these units. When Vernier technology is paired with OpenSciEd's classroom-tested curriculum, your students establish a deep understanding of critical scientific concepts through data collection.

Learn more at www.vernier.com/openscienced

Classic Approach



Vernier Lab Books

While the three-dimensional learning approach is valuable, sometimes a more classic approach to instruction is a better fit for your students, teaching style, and resources. In a classic approach, students follow detailed directions to conduct an experiment or investigate a specific science concept, topic, or law.

Vernier supports this more classic approach by providing a robust library of lab books covering most science disciplines. Our lab books provide teacher-created, step-by-step experiments that help your students work toward meeting the NGSS performance expectations and guide students through conducting hands-on experiments in a more structured way.

Learn more at www.vernier.com/lab-books


Vernier Supplements to OpenSciEd

GRADE 6

Thermal Energy



18 Lessons




Free Download
OSE-62TE-E

Students plan and carry out investigations to systematically test cup systems, tracking the flow of matter and energy into or out of the system as they develop a model of thermal energy.


Sensor Used

Go Direct®
Temperature
GDX-TMP \$69

Weather, Climate, and Water Cycling



22 Lessons



Free Download
OSE-63WC-E

In this Earth science unit, students use data-collection technology to explain small-scale storms, mesoscale weather systems, and global-level patterns of precipitation. In the culminating lesson, students explain how climate varies in different parts of the world.

Sensors Used

Go Direct
Temperature
GDX-TMP \$69

Go Direct
Light & Color
GDX-LC \$79

Go Direct Weather
GDX-WTHR \$99

GRADE 7

Metabolic Reactions



14 Lessons




Free Download
OSE-73MR-E

In this unit on metabolic reactions, students use a real case study of a middle school student to develop models to explain how the body uses food and how the body's subsystems work together.


Sensor Used

Go Direct
CO₂ Gas
GDX-CO₂ \$199

Matter Cycling and Photosynthesis



15 Lessons



Free Download
OSE-74MC-E


In this series of interactive experiments, students study the relationship between the food they eat and photosynthesis. Students investigate why plants need light, how they can survive without it, and so much more.

Sensor Used

Go Direct CO₂ Gas
GDX-CO₂ \$199

GRADE 8

Sound Waves



14 Lessons




Free Download
OSE-82SW-E

Students engage in model-based reasoning, argumentation, and computational and mathematical reasoning to develop models to explain what makes a sound, how sound moves through air, and how it makes something move.

Sensor Used

Go Direct Motion
GDX-MD \$99

Forces at a Distance



12 Lessons



Free Download
OSE-83FD-E

Students explore magnetism using Go Direct Sensor Carts to answer complex scientific questions such as how distance affects the strength of force pairs in a magnetic field.

Products Used

Go Direct Sensor Cart (Green) GDX-CART-G \$169
Go Direct Sensor Cart (Yellow) GDX-CART-Y \$169
Go Direct Sensor Cart Accessory Kit GDX-CART-AK \$58

Learn more at www.vernier.com/opensciEd

Life Science

Physical Science

Exploring Life Science



Download only
MSB-LS-E \$10

From yeast to humans, this e-book provides opportunities for students to learn about life science.

5 Experiments Included in E-book

Structure, Function, and Information Processing

- Get a Grip (*shown above*)
- Heart Rate and Body Position
- Heart Rate and Exercise

Matter and Energy in Organisms and Ecosystems

- Diffusion: How Fast?

Growth, Development, and Reproduction of Organisms

- Yeast Beasts in Action

Package Available

Exploring Life Science Go Direct Package

This package contains the following:
Go Direct Gas Pressure, Go Wireless® Heart Rate,
Go Direct Conductivity, Gas Pressure Sensor Bulb

GDP-MS-LS \$283

Buy 8 or more at \$275
and save \$64



Learn more at www.vernier.com/msb-ls-e

Middle School Explorations: Chemical Reactions



Download only
MSB-CR-E \$20

In the six experiments in this book, students gain an understanding of various types of chemical reactions as they build a model to explain what goes on at the molecular level during a chemical reaction.

6 Experiments Included in E-book

Students investigate endothermic and exothermic reactions, precipitate formation, conservation of mass, and other reactions.

Sensor Used



Go Direct Temperature

This is a rugged, general-purpose sensor that students can use to monitor temperature.

GDX-TMP \$69

Teacher pack also available

(includes 8 Go Direct Temperature Probes and a Charge Station)

GDX-TMP-TP \$599

Learn more at www.vernier.com/msb-cr-e

Physical Science

Exploring Physical Science



Watch
a video



Download only
MSB-PS-E \$20

From matter and energy to motion and forces, students explore a wide variety of topics in basic chemistry and physics in this e-book.

22 Experiments Included in E-book

Structure and Properties of Matter

- Fun with Pressure

Chemical Reactions

- Boiling Temperature of Water
- Freezing Temperature of Water
- How Low Can You Go? Freezer Bag Ice Cream

PLUS 2 MORE

Forces and Interactions

- Friction

- First Class Levers
 - Pulleys (*shown above*)
- PLUS 7 MORE

Energy

- A Hot Hand
- A Good Sock
- Lemon "Juice"

Waves and Electromagnetic Radiation

- Reflectivity of Light
- Mapping a Magnetic Field
- Electromagnets

Package Available Exploring Physical Science Go Direct® Package GDP-MS-PS \$642



This package contains the following Go Direct sensors: Temperature (2), Gas Pressure, Force and Acceleration, Motion Detector, Voltage, 3-Axis Magnetic Field, Light and Color

Buy 8 or more at \$623 and save \$152



Learn more at www.vernier.com/msb-ps-e

Exploring Motion and Force with Go Direct Sensor Cart



Download only
MSB-CART-E
\$20

In this e-book, students explore the force of friction, aspects of motion, and simple machines such as the lever, ramp, and pulley.

7 Experiments Included in E-book

- Investigating Friction
- Levers as Machines
- Pulleys as Machines (*shown above*)
- Ramps as Machines
- Getting Faster
- Crash Test
- Newton's Second Law

Package Available Exploring Motion and Force with Go Direct Sensor Cart Package

This package contains the following Go Direct sensors: Sensor Cart (Green) and Sensor Cart (Yellow)

GDP-MS-SC \$338

Buy 8 or more at \$328 and save \$80



Learn more at www.vernier.com/msb-cart-e

Earth and Space Science

Exploring Earth and Space Science



Download only
MSB-ESS-E \$20

Weather, soil, and water quality are a few of the Earth science topics students explore in this e-book.

12 Experiments Included in E-book

Earth's Systems

- Soil Study
- Ocean Floor Mapping
- Water Hardness Study
- A Water Field Study

Weather and Climate

- Heating of Land and Water
- The Greenhouse Effect
- Relative Humidity
- Absorption of Radiant Energy
- Reflectivity of Light
- Schoolyard Study
- What Causes the Seasons? (*shown above*)
- Solar Homes (Engineering Design)

Package Available

Exploring Earth and Space Science Go Direct Package

This package contains the following Go Direct sensors: Temperature (2), Light and Color, Motion Detector, Conductivity, pH


GDP-MS-ESS \$504

Buy 8 or more at \$489
and save \$120



Learn more at www.vernier.com/msb-ess-e

Climate and Meteorology Experiments



Download only
HSB-CM-E \$20

This new lab book is packed with interactive investigations that challenge students to use data-collection technology to explore weather, climate, and other important weather-related topics.

11 Experiments Included in E-book

Weather and Climate

- Modeling Solar Insolation
- What Causes Land and Sea Breezes?
- Investigating Albedo
- Exploring the Greenhouse Effect
- Effect of Air Temperature on Humidity
- What is Dew Point?
- Measuring Wind Chill
- Changes in Barometric Pressure
- Formation of Clouds
- Measuring Wind Direction
- Studying Microclimates: Urban Heat Islands

Package Available

Climate and Meteorology Experiments Go Direct Package

This package contains the following Go Direct sensors: Surface Temperature (2), Light and Color, Weather System

GDP-CM \$365

Buy 8 or more at \$354
for a savings of \$88



Learn more at www.vernier.com/hsb-cm-e

Engineering, Technology, and Coding

Coding with Scratch

Engage your students with scientific and computational concepts with our free activities for the popular coding platform Scratch.

Go Direct® Force and Acceleration brings real-world data into your Scratch project. With this integration, students can learn coding by purposefully and successfully applying their skills to fun, hands-on coding projects. This helps students make natural connections between the digital and physical worlds.



The LEGO® MINDSTORMS® Education EV3 Core Set with Charger is a hands-on cross-curricular STEM solution that engages students by providing the resources to design, build, and code their creations. Expand the possibilities of Scratch using this robotics kit.



LEGO-EV3-CORE
\$439.90

Learn more at www.vernier.com/scratch

Sam Labs and Vernier Technology

Bring STEAM, data collection, and coding to life for students with SAM Labs and Vernier sensor kits.



Teaching programming with SAM Blocks and Vernier sensors shows students how real experiment data interact with code to bring SAM Labs output blocks to life. When educators select from one of our custom packages based on science topics such as force or motion, they will get access to free, engaging activities to help students start coding right away.

Available Packages (includes sensor and output block)

Force	SL-FOR-PKG	\$167 (pictured above)
Magnetism	SL-3MG-PKG	\$145
Motion	SL-MD-PKG	\$167
Sound	SL-SND-PKG	\$153
Temperature	SL-TMP-PKG	\$100

Learn more at www.vernier.com/sam-labs

ACTIVITY 6

Driving Outside the Lines

Students write, run, and troubleshoot mBot code to navigate their mBot using dead reckoning; in other words, using time to measure and predict distance traveled and degrees turned.



Products Used in This Activity



Can also be done with

mBot Explorer

MBOT-S \$83.99

mBot™ by Makeblock®

mBot provides students with a fun and tactile way to learn entry-level coding with simple Scratch-based software. Included with your purchase is our *Coding with mBot: Self-Driving Vehicles* e-book.

MBOT-P (pink) or MBOT-B (blue) \$69.99 each

Experiment Source



**Coding with mBot:
Self-Driving Vehicles**

FREE DOWNLOAD*

*Free with purchase
of mBot from Vernier

Learn more at www.vernier.com/mbot-msdv-e


Wind Energy

Wind Energy Explorations

Students gain an understanding of energy, circuits, and loads, as well as practice engineering design as they use this e-book to explore wind energy.

Experiments Included in E-book

- Energy Transformation
- Measuring Wind Energy
- Exploring Wind Turbines
- Wind Turbines: Effect of Load
- Blade Variable: Pitch
- Blade Variable: Quantity
- Blade Variable: Area
- Blade Variable: Shape
- Project: Max Power (Engineering Design)



INCLUDES
9
EXPERIMENTS

Download only
MSB-WIND-E \$20

Wind Energy Explorations Go Direct Packages

Single Station Package *(shown below)*

This package includes

- Go Direct® Energy (1)
- Vernier Resistor Board (1)
- KidWind Basic Wind Experiment Kit (1)

GDP-MS-WE \$231



Classroom Package

This package includes

- Go Direct Energy Sensors (3)
- Vernier Resistor Boards (3)
- KidWind Basic Wind Experiment Classroom Pack (includes materials for 6 to 10 groups of 2 to 4 students each) (1)

GDP-MS-WEC \$630

Learn more at www.vernier.com/msb-wind-e


Solar Energy

Solar Energy Explorations

Solar energy provides a relevant topic for students to explore energy, temperature, and electrical circuits, culminating in an engineering design project.

Experiments Included in E-book

- Renewable Energy
- Introduction to Solar Panels and Solar Energy
- Measuring Energy
- Making Connections: Circuits
- Solar Panel Output: Effect of Load
- Solar Panel Output: Effect of Shade
- Solar Panel Output: Effect of Angle
- Solar Panel Output: Effect of Temperature
- Project: Build a Solar Car (Engineering Design)



INCLUDES
9
EXPERIMENTS

Download only
MSB-SOLAR-E \$20


Solar Energy Explorations Go Direct Package

This package includes two sensors that both work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3. It also includes an experiment kit and a resistor board.

- Go Direct Energy
- Solar Energy Exploration Kit
- Go Direct Surface Temperature
- Vernier Resistor Board

GDP-MS-SE \$265

Buy 8 or more packages at \$257 and save \$64














Learn more at www.vernier.com/msb-solar-e

Featured Products

Go Direct Sensors

Sensor		Order Code	Price
Go Direct® 3-Axis Magnetic Field		GDX-3MG	\$69
Carts and Tracks			
Dynamics Cart and Track System with Go Direct Sensor Carts		DTS-GDX	\$535
Go Direct Sensor Cart (Green)		GDX-CART-G	\$169
Go Direct Sensor Cart (Yellow)		GDX-CART-Y	\$169
Go Direct Conductivity		GDX-CON	\$99
Go Direct Current		GDX-CUR	\$79
Go Direct Energy		GDX-NRG	\$89
Go Direct Force and Acceleration		GDX-FOR	\$99
Go Direct Gas Pressure		GDX-GP	\$89
Go Wireless® Heart Rate		GW-HR	\$89

Go Direct Light and Color		GDX-LC	\$79
Go Direct Motion		GDX-MD	\$99
Go Direct Optical Dissolved Oxygen		GDX-ODO	\$298
pH Sensors			
Go Direct pH		GDX-PH	\$89
Go Direct Tris-Compatible Flat pH		GDX-FPH	\$115
Go Direct Sound		GDX-SND	\$89
Go Direct Structures & Materials Tester		GDX-VSMT	\$999
Temperature Probes			
Go Direct Surface Temperature		GDX-ST	\$79
Go Direct Temperature		GDX-TMP	\$69
Go Direct Voltage		GDX-VOLT	\$69
Go Direct Weather System		GDX-WTVA	\$128

See all our products for middle school science at www.vernier.com/middle-school

Looking for Replacement Parts?

Visit www.vernier.com/replacements







Go Direct Charge Station

Accessory	Order Code	Price
Go Direct Charge Station 	GDX-CRG	\$69

LabQuest 3 Interface and Sensors

Learn more about LabQuest® 3 and sensors at www.vernier.com/labq3

Additional Products

Products	Order Code	Price
Davis® Weather Stations 	www.vernier.com/weather	
pH Storage Solution 	PH-SS	\$20
KidWind Basic Wind Experiment Kit 	KW-BWX	\$124
OHAUS® Balances 	www.vernier.com/ohaus	
Solar Energy Exploration Kit 	KW-SEEK	\$79
Vernier Resistor Board 	VES-RB	\$18

Coding and Robotics

Products	Order Code	Price
Go Direct Force and Acceleration (for use with Scratch)	GDX-FOR	\$99
LEGO® MINDSTORMS® Education EV3 Core Set with Charger	LEGO-EV3-CORE	\$439.90
mBot™ (blue)	MBOT-B	\$73.49
mBot (pink)	MBOT-P	\$73.49
mBot Explorer	MBOT-S	\$83.99
Vernier Coding with SAM Labs Packages	varies	\$100–\$167

Lab Books

Title	Order Code	Price
<i>Middle School Science with Vernier</i>	Download + print: MSV Download only: MSV-E	\$48 \$40
<i>Exploring Motion and Force with Go Direct Sensor Cart</i>	MSB-CART-E	\$20
<i>Exploring Physical Science*</i>	MSB-PS-E	\$20
<i>Exploring Life Science*</i>	MSB-LS-E	\$10
<i>Exploring Earth and Space Science*</i>	MSB-ESS-E	\$20
<i>Solar Energy Explorations</i>	MSB-SOLAR-E	\$20
<i>Wind Energy Explorations</i>	MSB-WIND-E	\$20
<i>Coding with mBot: Self-Driving Vehicles</i>	MBOT-MSDV-E	\$20
<i>Earth Science with Vernier</i>	Download + print: ESV Download only: ESV-E	\$48 \$40
<i>Climate and Meteorology Experiments</i>	HSB-CM-E	\$20

See all our products for middle school science at www.vernier.com/middle-school

* All experiments from this e-book are included in Middle School Science with Vernier.

High School

www.vernier.com/high-school

Encourage your students and build their confidence in pursuing a STEM career path with hands-on experience using data-collection technology from Vernier. Our technology supports you as you set up students for success for standardized testing, as well as prepare them to meet the NGSS and state standards through experiments that support three-dimensional learning.



Contents

Explore a sampling of our featured experiments by topic to learn how Vernier technology helps your students deepen their understanding of key STEM concepts.

Lab Books & Investigations

PAGE 29

A Guide to Vernier Data Collection

PAGE 30

LabQuest® 3

PAGE 32

Interfaces

PAGE 35

Software and Digital Curriculum

PAGE 36

Subjects

BIOLOGY
PAGE 42

ENVIRONMENTAL
SCIENCE
PAGE 58

EARTH SCIENCE
PAGE 70

CHEMISTRY
PAGE 74

PHYSICAL SCIENCE
PAGE 92

PHYSICS
PAGE 96


ENGINEERING, CODING,
AND ROBOTICS
PAGE 122

TEXAS INSTRUMENTS
PAGE 132

College

www.vernier.com/college

Lab Books & Investigations



E-books and Printed Books—the Choice is Yours

Many of our popular, award-winning lab books are available in both e-version and printed formats. When you purchase a printed book, you also receive the electronic version. When you purchase either format, you receive

- Anytime access to the most up-to-date versions of experiments on all supported Vernier software (free Vernier web account required)
- Editable student files and complete teacher information files, including sample data and supplies lists
- A generous site license—purchase once and share files with other teachers in your school

Helping You Meet Standards and Learning Objectives

Vernier understands that helping students meet standards is an important part of teaching. As standards change, we are committed to providing you with the most current information. You will find the following alignments and correlations for Vernier lab books at

www.vernier.com/standards

- NGSS (Next Generation Science Standards)
- CSTA (Computer Science Teachers Association)
- AP* (Advanced Placement Program)
- IB† (International Baccalaureate Diploma Program)

Ideas for Your Science Classroom

If you are looking for experiments that can help you excite your students about STEM, check out our extensive library of experiments. We make it easy to find ideas from fellow educators and Vernier professionals.

Visit www.vernier.com/ideas

NGSS Aligned


To learn about the Next Generation Science Standards and Vernier, visit www.vernier.com/ngss

Learn more at www.vernier.com/lab-books

* AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

† The IB Diploma Program is an official program of the International Baccalaureate Organization (IBO) which authorizes schools to offer it. The material available here has been developed independently of the IBO and is not endorsed by it.

Getting Started with Go Direct Sensors



Why Choose Go Direct Sensors?

With over 50 sensors to choose from, our Go Direct® family of sensors offers an affordable solution that includes free software. Go Direct sensors are easy to use—just connect and start collecting data with your device.

What You Need to Get Started

A Go Direct Sensor

These versatile sensors connect to your device via Bluetooth® wireless technology or USB.

B Device

Go Direct sensors connect to a wide variety of devices commonly used in classrooms, including Chromebooks, computers, compatible mobile devices, and LabQuest® 3.

C Vernier Graphical Analysis™ App

Our data-collection app facilitates student understanding with real-time graphs of experimental data.

The new Graphical Analysis Pro app offers additional exclusive features, such as the ability to perform live experiments and share the data over the internet in real time. Sign up for a free 30-day trial today!


Learn more at
www.vernier.com/graphical-analysis-pro

D Lab Book

Step-by-step instructions at your fingertips save valuable time when integrating probeware into your curriculum. Many of our lab books provide support for Go Direct sensors and the Graphical Analysis app.

Our lab books come with a generous site license. Purchase once and share files school wide.

Getting Started with LabQuest 3



Why Choose LabQuest 3?

LabQuest 3 is a powerful, easy-to-use, and versatile data-logging solution for STEM students. A full-featured data-collection platform, LabQuest 3 is an excellent choice for laboratories, classrooms, or in-the-field investigations.

What You Need to Get Started

A LabQuest 3

With its large, high-resolution screen, LabQuest can be easily navigated using gestures. It also offers fast data collection, wireless connectivity with Wi-Fi and Bluetooth wireless technology, and a rechargeable, high-capacity battery.

B Sensors

Compatible with all Vernier sensors, LabQuest 3 connects wirelessly to the family of Go Direct sensors and connects easily with our wired LabQuest sensors.

C Software

LabQuest 3 has built-in software, LabQuest App, that gives your students real-time graphing and analysis capabilities in one handheld device. LabQuest 3 offers built-in apps, such as a Periodic Table, Sound Recorder, and more, and includes student instructions for over 75 of our most popular experiments.

D Lab Book

Looking for even more lab ideas? Our popular, award-winning lab books provide hundreds of well-tested, customizable experiments.

Our lab books come with a generous site license. Purchase once and share files school wide.

NEW LabQuest 3



LabQuest 3 is a powerful, easy-to-navigate, and versatile data-logging solution for STEM students.

The all-new LabQuest® 3 reimagines data collection by providing students with an innovative, easy-to-use interface. A larger screen and advanced touch screen abilities makes it easier for students to collect, graph, and analyze data wherever they are—the classroom, at home, or in the field. Challenge your students to gain a deeper understanding of science through data with the accessible, groundbreaking LabQuest 3.

- Connects wirelessly to the family of Go Direct® sensors
- Easy-to-use platform enables students to generate graphs and analyze results
- An excellent choice for laboratories, classrooms, or in-the-field investigations

LABQ3 \$339

LabQuest 3 purchase includes: LabQuest 3 unit, Rechargeable battery (in unit), AC power adapter, Micro-USB computer connection cable, and Quick-Start Guide



Full-Featured Data-Collection Platform

The most engaging and effective approach to science is interactive, with students collecting and analyzing data to understand and apply core concepts. Graphing and analyzing data is an essential component of the inquiry and learning process. LabQuest 3, with its built-in data-collection and analysis app that works with all Vernier sensors, supports hands-on data collection in the classroom, in the lab, and in the field.

- Is a Chromebook™ not available? No problem. LabQuest 3 can do it all—data collection, data analysis, and data sharing.
- Keep your expensive computers safe from spills, drops, and crashes—use LabQuest 3 in the chemistry lab, at the watershed, or next to your bridge tester. LabQuest 3 does not need another device for data collection or analysis.
- With a portable design, LabQuest 3 lets your students take it anywhere they go.
- LabQuest 3 works with both LabQuest and Go Direct sensors.



Connectivity to Other Platforms

One-to-Many Data Sharing

Students can share real-time data with multiple devices for a truly hands-on, collaborative learning experience. Use LabQuest 3 to transfer data wirelessly to computers, Chromebooks, or mobile devices running Vernier Graphical Analysis™.

USB Sensor Interface

If you want to use your own computer or Chromebook™ to collect data, use LabQuest 3 as a conduit between our wired LabQuest sensors and these devices. LabQuest 3 works as a USB sensor interface with our Logger Pro® software, Vernier Graphical Analysis™ app, or Vernier Graphical Analysis Pro.


Learn more at www.vernier.com/labq3

LabQuest 3

LabQuest App


LabQuest 3 has built-in software that gives your students real-time graphing capabilities in a handheld device. It's powerful, yet beautifully simple.

- Collect data and view them in a Data Table, Meter, and Graph.
 - Perform curve fits.
 - Use built-in sensors—GPS, accelerometers, and more.
 - Draw a prediction before collecting data.
 - Display two graphs at once.
 - Display a tangent line or use the Integral function tool.
 - Calculate statistics for your data.
- Learn more about built-in applications and other great features at www.vernier.com/labq3



One-Touch Simplicity

Your students can collect data and view them in a Meter, Graph, or Data Table.



Meter



Graph

Site	Latitude (°)	Longitude (°)	Altitude (m)	DO Concentration (mg/L)
new stream	45.50782	-122.85773	66	10.07
pond	45.50790	-122.85690	61	10.71
entrance	45.50841	-122.85613	51	13.77

Data Table

Learn more at www.vernier.com/labq3

Accessories and Replacement Parts

Product	Order Code	Price
LabQuest Charge Station	LQ3-CRG	\$129
LabQuest 3 Stand	LQ3-STN	\$5
LabQuest Power Supply*	LQ3-PS	\$14
LabQuest Lanyard	LQ3-LAN	\$9
LabQuest 3 Battery	LQ3-BAT	\$34
LabQuest Battery Boost 3	LQ-BOOST3	\$119
Vernier Micro USB Cable*	CB-USB-MICRO	\$5
Vernier Micro USB to USB-C Cable	CB-USB-C-MICRO	\$9

*Included with LabQuest 3

LabQuest Viewer App



LabQuest Viewer®

Teach students how to use LabQuest® by projecting your LabQuest screen. Display live images of all LabQuest units in your lab to monitor student progress or compare group data. This is compatible with both macOS® and Windows® computers.

Computer software includes a site license for every teacher's computer in your school.

CD: LQ-VIEW \$79

Download: LQ-VIEW-E \$79

For more information, visit www.vernier.com/lq-view


LabQuest Viewer for iPad®

Use LabQuest Viewer app for iPad on your classroom iPad to wirelessly view and control LabQuest. When your iPad is used with a projector, you can easily display any LabQuest screen for the entire class to see.

For more information, visit www.vernier.com/lq-view-ipad



Download on the
App Store




LabQuest Mini



LabQuest Mini

LabQuest Mini brings the power of our award-winning LabQuest technology to you when you don't need the versatility of a standalone device. The perfect solution for educators collecting data with a computer or Chromebook,[™] LabQuest Mini interfaces with Vernier Graphical Analysis,[™] Vernier Graphical Analysis Pro, and Logger Pro[®] software.

LQ-MINI \$169



Learn more at www.vernier.com/lq-mini


LabQuest Stream



LabQuest Stream[®]


LabQuest Stream brings data collection with LabQuest sensors to even more platforms—computers, Chromebooks, smartphones, and tablets. LabQuest Stream makes a one-to-one connection to your technology either via USB or wirelessly via Bluetooth[®] wireless technology without the need to connect to your school's network. LabQuest Stream is our recommended interface for BYOD classrooms using LabQuest sensors.

LQ-STREAM \$229




Learn more at www.vernier.com/lq-stream

Vernier Graphical Analysis



View a graph, table, and meter simultaneously.



Use analysis tools, including text annotations and statistics.

Collect, share, and analyze sensor data with our free software for Chrome OS™, iOS, iPadOS®, Android™, Windows®, and macOS®.

Using the Vernier Graphical Analysis™ app, you can collect data from Go Direct® sensors or LabQuest® sensors connected to a compatible interface.

Enter data manually, copy data saved on your clipboard, or receive data from a Data Sharing source (LabQuest 3 or Logger Pro® 3) using Wi-Fi.

Free Download

Learn more and download Vernier Graphical Analysis for free at www.vernier.com/graphical-analysis

Key Features


Data Collection

- Collect data from multiple sensors simultaneously.
- Select time-based or event-based data collection, including events with entry.
- Adjust data-collection rate and duration as needed.
- Enter data manually or using the clipboard.
- Draw predictions before data collection.
- Perform graph matching exercises with a motion detector.


Data Analysis

- View data in a meter, on a graph, in a table, or all three at once.
- Display one, two, or three graphs as needed.
- Easily select what columns and data sets are plotted on each graph.
- Calculate descriptive statistics and fit lines and curves to some or all of your data.

NEW Vernier Graphical Analysis Pro



Videos synchronized with sensor data help students understand experiment phenomena such as simple harmonic motion.



With the included sample experiments, students can experience an experiment like Boyle's Law even when lab equipment is unavailable.

Make Virtual Science Real with Our Reimagined Graphical Analysis Pro App

Our award-winning Vernier Graphical Analysis™ app went Pro with a new, enhanced version that enables students to engage in real-time experimentation and analysis of data—either at home or in the classroom. With Vernier Graphical Analysis Pro, educators can create experiments and share the data with students in real time. Seeing data collected right before their eyes gives students the ability to connect abstract concepts to real-world applications. Plus, the enhanced features of Graphical Analysis Pro help students experience three-dimensional learning. As always, we strive to make your job as an educator easier. This is why we created a wide variety of videos with sample experiments synced with data and complete instructions that cover common topics in biology, chemistry, and physics.

Free Trial for Educators


Try out Graphical Analysis Pro for free for 30 days. Access the sample experiment library and enhanced analysis tools to use with your students and make virtual science real!

Get a free trial and learn about site license options at www.vernier.com/graphical-analysis-pro


Key Features

- With Graphical Analysis Pro, educators can create their own live experiments using their Vernier sensors and share the data in real time to students. Educators have the power to create their own videos—synced with actual data—to distribute to students studying remotely, and students with access to Vernier sensors can create their own videos to share with others.
- Educators can draw from our sample experiments with synced data to share with students.
- Graphical Analysis Pro is compatible with most Vernier sensors, so educators can continue to use them, even as students learn remotely.
- Our new app is compatible with multiple computer operating systems and mobile device platforms—including Chromebooks, which provides flexibility and cost savings, as students can use their own device for analysis.
- Graphical Analysis Pro offers an intuitive interface making it easy to use and get started.
- As always, we offer great customer service from the experts and educators on our technical support team.


Logger Pro 3



After you click Collect, Logger Pro 3 draws the graph in real time, and the data table and digital meter update continuously.



Plot position, velocity, and acceleration data from a Motion Encoder Cart.



Collect absorbance data from Vernier spectrometers, including our Go Direct SpectroVis Plus and Vernier UV-VIS Spectrophotometers.

Real-Time Graphing and Powerful Analytical Tools

Logger Pro® 3 is our data-collection and analysis software for LabQuest sensors on Windows® and macOS® computers. With a complete suite of data-collection and analysis tools, Logger Pro 3 is suitable for all students, from beginning to advanced.

One program does it all—for only \$249—for all of your school's computers AND your students' personal computers.

Logger Pro 3 can gather data from a variety of sources, including LabQuest® 3, LabQuest Mini, LabQuest Stream®, Go!Link®, OHAUS® balances, compatible TI graphing calculators, and spectrometers.

Key Features

Logger Pro 3 includes a site license for your entire high school.

- Site license includes home computers of teachers and students

Logger Pro 3 Data Sharing

- Use Logger Pro 3 for lecture demonstrations. Collect data on your computer and Data Share your data to student devices running our free Vernier Graphical Analysis™ app or Graphical Analysis Pro app.

Advanced Features

- Import remotely collected data from LabQuest 3 and TI-84 Plus calculators.
- Lay out graphs, tables, and text across multiple pages to describe your experiment.

- Graph data in a variety of ways, including log graphs, double-Y graphs, strip charts, and FFT graphs.
- Model data with user-adjustable functions.
- Extract data from movies using frame-by-frame video analysis.
- Capture video from video cameras or import compatible movie files.
- IB* curriculum support—manual curve fits and error bars

Note: Logger Pro 3 cannot be used to collect data with our Go Direct® sensors (other than Go Direct SpectroVis® Plus).

*The IB Diploma Program is an official program of the International Baccalaureate Organization (IBO) which authorizes schools to offer it. The material available here has been developed independently of the IBO and is not endorsed by it.

Logger Pro 3

with manual, CD, and download

LP \$249


download only

LP-E \$249


Windows® and macOS® computers only

Learn more at www.vernier.com/logger-pro

Vernier Spectral Analysis



Absorbance spectra of green food coloring at different concentrations



Wavelength selection screen for Beer's Law and kinetics experiments

Collect, share, and analyze spectrometer data with our free software for Chrome OS,[™] Windows,[®] macOS,[®] iOS, iPadOS,[®] and Android.[™]

Benefits

The free Vernier Spectral Analysis[®] app makes it easy to incorporate spectroscopy into your biology, chemistry, and physics experiments. Using the app, students can collect a full spectrum and explore topics such as Beer's law, enzyme kinetics, and plant pigments.


The user-friendly software includes analysis features such as curve fitting and data interpolation.

Features

- Follow on-screen instructions for simplified Beer's law or kinetics data collection.
- Collect full absorbance spectrum or % transmittance data in less than one second.
- Analyze data with built-in analysis tools, including data interpolation and curve fittings.
- Determine the order of kinetics reaction with the calculated columns function.
- Understand color transmission using the color strip shown on full spectrum graphs.
- View a full spectrum of your sample while collecting data for Beer's law or kinetic experiments.
- View spectral lines by collecting intensity vs. wavelength data.

Learn more at www.vernier.com/spectral-analysis

Vernier Instrumental Analysis



The separation of cyclohexane and toluene

Compatible with Chrome OS, iOS, iPadOS, Android, Windows, and macOS


Our free Vernier Instrumental Analysis[™] app makes it easy to incorporate instrumentation into your chemistry curriculum. With this app, students can collect and analyze data from Mini GC, Mini GC Plus, Go Direct Mini GC,[™] and Go Direct Polarimeter using computers, Chromebooks, or other mobile devices.

Learn more at www.vernier.com/instrumental-analysis

Vernier Video Analysis



WORKS ON
CHROMEBOOKS!



Investigate projectile motion

Study Motion Everywhere

The Vernier Video Analysis™ app brings video analysis to your students in an easy-to-use, streamlined application.

Benefits

Vernier Video Analysis makes studying motion easy and accessible. Students can design their own scientific investigations, record videos, and then analyze the motion. This app brings video analysis to all your students regardless of device—it even works with Chromebooks!

Free 30-Day Trial

Get a 30-day free trial and learn about site license options at www.vernier.com/video-analysis

Features


- Vernier Video Analysis app is compatible with multiple devices and platforms: macOS®, iPadOS®, iOS, Windows® 10, Chrome OS™, and Android™.
- Students can use prepared videos, found videos, or collect their own videos for analysis.
- The app makes it possible to do experiments that cannot be done with sensors, such as analyzing the motion of a basketball in flight—objects can be tracked automatically by the app.
- Analysis is easy with multiple graphing options, so students are able to think critically about the collected data—they can even analyze the motion of multiple objects in a single video.
- With this app, you can apply vectors and vector components over the video after tracking a moving object, illuminating changes in position, velocity, and acceleration.
- When multiple objects have been marked, just enter their masses and the app can automatically calculate and display the center of mass location.
- Annual site-licensing makes purchasing and renewing quick and easy.

NEW

Vernier Video Analysis: Motion and Sports

The *Vernier Video Analysis: Motion and Sports* lab book features 12 investigations using Vernier Video Analysis. In addition to traditional physics concepts such as velocity and acceleration, its investigation of sports activities expands learning opportunities and further connects the study of motion to students' daily lives.

Download only
HSB-VVAMS-E
\$25



Pivot Interactives



Deepen Student Understanding with Pivot Interactives

Benefits

Pivot Interactives is a powerful supplement to hands-on experimentation, enabling students to vary experimental parameters one at a time to view results from a set of many recordings of the same experiment. These high-quality videos give your students the opportunity to observe and study hard-to-replicate phenomena. Students make measurements and analyze their data directly within the Pivot Interactives online environment.

Features

- Augment hands-on learning with interactive videos to teach concepts in biology, chemistry, and physics.
- Use Pivot Interactives for formative and summative assessment.
- Assign pre-made activities to students or author new ones.
- Provide feedback to students through Pivot Interactives.
- Pricing for high schools is \$5 per seat (10-seat minimum) with site licensing pricing available.

Free Trial for Educators

Start a free 30-day trial today at www.pivotinteractives.com

Try Pivot Interactives free for 30 days. Browse the entire library of videos, explore the analysis tools, and use it with your students.



Students build cellular respiration equation models by observing and collecting data from a diverse set of organisms.

See Pivot
in Action



Watch
a video



Students change what happens in the videos by varying a parameter such as acid-base combination or indicator, and observe how it affects the outcome.



Students measure the total power output of the sun by comparing the intensity of the sun's light at Earth's surface to the intensity of a known source of light.

Biology

www.vernier.com/biology

Our biology solutions include high-quality sensors, easy-to-use software, and exceptional technical support to set up you and your students for classroom success.

Topics

Explore our featured experiments by topic to learn how Vernier technology helps your students engage with data-collection technology and deepens their understanding of key biological concepts.

Biology

PAGE 44

Human Physiology

PAGE 48

Agricultural Science

PAGE 51

Spectroscopy

PAGE 52


Biotechnology

PAGE 54



Bring Your Biology Lessons to Life

From cellular biology to ecology to human physiology, get your students excited about biology using Vernier technology. Our sensors, software, and investigations help biology students explore phenomena, develop their understanding of living organisms, and encourage their scientific curiosity. Work with our team to implement high-quality sensors, experiments, and technology solutions in your classroom and set your students up for success in science and beyond.



Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training

EXPERIMENT 11

Cell Respiration

Students measure cellular respiration in germinating peas and determine what effect temperature has on respiration rate.



Sensor Used



Go Direct® CO₂ Gas

Use Go Direct CO₂ Gas to measure CO₂ gas levels, air temperature, and relative humidity. It's an excellent sensor for measuring fermentation, cell respiration, and photosynthesis.

GDx-CO2 \$199

Experiment Source



Biology with Vernier

Download only: BWV-E \$40

Printed book + download: BWV \$48

Learn more at www.vernier.com/bwv-11b

EXPERIMENT 6

Enzyme Action

Students measure the activity of the enzyme catalase and analyze how different factors (e.g., enzyme concentration, pH, and temperature) influence enzyme activity.

Watch a video



Sensor Used



Go Direct Gas Pressure

Use Go Direct Gas Pressure to monitor gas pressure in a variety of experiments. Easily change the displayed units to any one of seven options. This sensor includes a syringe, tubing, and stoppers to ease experiment setup.

GDx-GP \$89

Experiment Source



Biology with Vernier

Download only: BWV-E \$40

Printed book + download: BWV \$48

Learn more at www.vernier.com/bwv-6b

INCLUDES
31
EXPERIMENTS

EXPERIMENT 1

Energy in Food

Students determine and compare the energy content of different foods using calorimetry.



Sensor Used



Go Direct Temperature

This rugged probe measures the temperature of a variety of substances including air, soil, and water.

GDX-TMP \$69

Experiment Source



Biology with Vernier

Download only: BWV-E \$40

Printed book + download: BWV \$48

Learn more at www.vernier.com/bwv-1

Biology with Vernier

Biology with Vernier addresses the fundamentals of a high school biology course with 31 experiments that include cell respiration, photosynthesis, membrane diffusion, osmosis, human physiology, transpiration, fermentation, and more.

The instructor information section included for each experiment contains reagent preparation information, sample data, and tips for successful completion.

Learn more at www.vernier.com/bwv



Download only

BWV-E \$40

Printed book + download

BWV \$48

Biology Go Direct Starter Package

This package includes four sensors, which all work with our free Vernier Graphical Analysis™ 4 app, as well as Graphical Analysis Pro and LabQuest® 3.

- Go Direct Temperature
- Go Wireless® Heart Rate
- Go Direct Gas Pressure
- Go Direct CO₂ Gas

GDP-BIO-ST \$446

Learn more at www.vernier.com/gdp-bio-st

Standard package also available
(see page 49)



EXPERIMENT 25

Primary Productivity

Measuring the effect of light level on net and gross productivity in aquatic ecosystems helps students understand primary productivity.



Sensor Used



Go Direct® Optical Dissolved Oxygen

Use this sensor to measure dissolved oxygen, water temperature, and atmospheric pressure.

GDX-ODO \$298

Accessory Used



Primary Productivity Kit

This kit is an accessory for one of our most popular biology experiments, "Primary Productivity." The kit consists of a box of 7 plastic bottles, 7 rubber stoppers, and a set of screens.

PPK \$44

Experiment Source



Biology with Vernier

Download only: BWV-E \$40

Printed book + download: BWV \$48

Learn more at www.vernier.com/bwv-25

EXPERIMENT 31

Photosynthesis and Respiration (CO₂ & O₂)

Students use a terrestrial plant to measure photosynthesis and cellular respiration.



Sensors Used



Go Direct CO₂ Gas

Measure gaseous carbon dioxide concentration levels, air temperature, and relative humidity using this sensor.


GDX-CO2 \$199



Go Direct O₂ Gas

Use this sensor to measure gaseous oxygen concentration levels and air temperature.

GDX-O2 \$189



BioChamber 2000

BC-2000 ⚡ \$22

Experiment Source



Biology with Vernier

Download only: BWV-E \$40

Printed book + download: BWV \$48

Learn more at www.vernier.com/bwv-31c

Biology Go Direct Standard Package

This package includes 11 sensors that all work with our free Vernier Graphical Analysis™ 4 app, as well as Graphical Analysis Pro and LabQuest® 3. Two sampling chambers are also included.


- Go Direct Temperature
- Go Direct Gas Pressure
- Go Direct O₂ Gas
- Go Direct CO₂ Gas
- Go Direct Colorimeter
- Go Direct Conductivity
- Go Direct EKG
- Go Direct pH
- Go Direct Optical Dissolved Oxygen
- Go Direct Respiration Belt
- Go Wireless® Heart Rate
- BioChamber 250 ⚡
- BioChamber 2000 ⚡

GDP-BIO-ODX ⚡ \$1,528

Buy 8 or more packages at \$1,482 and save \$368

Learn more at
www.vernier.com/gdp-bio-odx

Starter package also available



Pivot Interactives for Biology

Pivot Interactives is a powerful supplement to hands-on experimentation, allowing students to vary experimental parameters one at a time to view results from a set of many recordings of the same experiment.

Start a free 30-day trial today at www.pivotinteractives.com



Biology Lab Books




Biology with Vernier

Download only BWV-E \$40

Printed book + download BWV \$48

31 Experiments




Advanced Biology with Vernier*

Download only BIO-A-E \$40

Printed book + download BIO-A \$48

17 Experiments

* Instructions for Graphical Analysis app not yet available



Investigating Biology through Inquiry

Download only BIO-I-E \$40

Printed book + download BIO-I \$48

22 Investigations

AP[†] AND IB[‡] CORRELATIONS

To see all AP[†] correlations, visit www.vernier.com/ap-correlations

[†] AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

To see all IB[‡] correlations, visit www.vernier.com/ib-correlations

[‡] The IB Diploma Program is an official program of the International Baccalaureate Organization (IBO) which authorizes schools to offer it. The material available here has been developed independently of the IBO and is not endorsed by it.

EXPERIMENT 8

Introduction to Electrocardiography

After obtaining graphical representations of the electrical activity of the heart, students learn to recognize the different waveforms in an EKG and associate them with events in the heart.



Sensor Used



Go Direct® EKG

Go Direct EKG measures electrical activity in the heart and electrical signals produced during muscle contractions.

GDX-EKG \$159

Experiment Source



Human Physiology Experiments : Volume 1

Download only: HSB-HP-E \$30

Printed book + download: HSB-HP \$38

Learn more at www.vernier.com/hsb-hp-8

EXPERIMENT 7

Effect of Exercise on Heart Rate

Observing and measuring how the heart responds to exercise is a fun, hands-on way for students to learn about the cardiovascular system.



Sensor Used



Go Wireless® Heart Rate

This sensor is ideal for continuously monitoring heart rate before, during, and after exercise or while a person is stationary.

CW-HR \$89

Experiment Source



Human Physiology Experiments : Volume 1

Download only: HSB-HP-E \$30

Printed book + download: HSB-HP \$38

Learn more at www.vernier.com/hsb-hp-7

EXPERIMENT 1

Blood Pressure and Autonomic Reflexes

Using a blood pressure sensor, students can compare blood pressures taken before and after exposure to cold. Students obtain graphical representation of blood pressure and observe an example of “fight or flight” response.



Sensor Used



Go Direct Blood Pressure

Designed for versatility, Go Direct Blood Pressure is a non-invasive sensor that measures human blood pressure—systolic, diastolic, and mean arterial pressure—using the oscillometric method.

GDX-BP \$105

Experiment Source



Human Physiology Experiments: Volume 2

Download only: ALB-HP2-E \$30

Printed book + download: ALB-HP2 \$38

Learn more at www.vernier.com/alb-hp2-1

Human Physiology Go Direct Standard Package

This package includes 9 sensors that all work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3. Two useful accessories are also included.

- Go Direct Blood Pressure
- Go Direct EKG
- Go Direct Force and Acceleration
- Go Direct Hand Dynamometer
- Go Direct O₂ Gas
- Go Direct Respiration Belt
- Go Direct Surface Temperature
- Go Direct Spirometer
- Go Wireless Heart Rate
- Reflex Hammer Accessory Kit
- BioChamber 250

GDP-HP-DX \$1,164

Buy 8 or more packages at \$1,129 and save \$280

Learn more at www.vernier.com/gdp-hp-dx

Starter package also available (see page 45)



PLTW Biomedical Science

PLTW Biomedical Science (9–12) inspires students to make an impact on others' lives and empowers them to pursue their life and career goals—whether it's a future in diagnosing, treating, or preventing disease.

Learn more at www.vernier.com/pltw

Learn more about PLTW Engineering

See page 128

Featured Products

NEW Human Physiology Experiments: Volume 2

Human Physiology Experiments: Volume 2 contains 15 experiments designed to encourage students to explore the physiology of various human organ systems. An expansion of our *Human Physiology Experiments: Volume 1* lab book, the setup for these experiments is minimal—students are collecting data within minutes.

Download only ALB-HP2-E \$30

Printed book + download ALB-HP2 \$38

Download a free sample
experiment at
www.vernier.com/alb-hp2



Go Direct Blood Pressure

Go Direct® Blood Pressure is an affordable, non-invasive sensor designed to easily measure human blood pressure. It measures systolic, diastolic, and mean arterial pressure using the oscillometric method. Go Direct Blood Pressure can also report pulse rate and can display both individual pressure pulses and peak-to-peak pulse amplitudes, giving students a few ways to collect data.

GDX-BP \$105

www.vernier.com/gdx-bp



Reflex Hammer Accessory Kit

The Reflex Hammer Accessory Kit converts your Vernier force sensor into a reflex hammer. Use it to capture the strike of the hammer on a tendon. When using the kit with an EKG sensor to record EMGs, students can study reflexes.

RFX-ACC \$29

www.vernier.com/rfx-acc



Go Direct Respiration Belt

Go Direct Respiration Belt uses a force sensor and an adjustable nylon strap to measure human respiration rates before, during, and after exercise.

GDX-RB \$99

www.vernier.com/gdx-rb



Go Direct Spirometer

This is a multi-channel sensor that reports air pressure, flow rate, volume, and respiration rate.

Measuring tidal volumes and other lung function parameters are both simple and easy due to channels that automatically adjust for baseline drift.

GDX-SPR \$199



www.vernier.com/gdx-spr

Agricultural Science

Featured Products

EXPERIMENT 13

Transpiration

Students measure the rate of transpiration from a plant and then investigate how different environmental factors influence water transport in plants.



Sensor Used



Go Direct Gas Pressure

Use Go Direct Gas Pressure to monitor gas pressure in a variety of experiments. Easily change the displayed units to any one of seven options. This sensor includes a syringe, tubing, and stoppers to ease experiment setup.

GDX-GP \$89

Experiment Source



Agricultural Science with Vernier

Download only: AWV-E \$40

Printed book + download: AWV \$48

Learn more at www.vernier.com/awv-13

LabQuest 3

LabQuest® 3 is a powerful, connected, and remarkably versatile data-logging solution.

Why? LabQuest 3 can serve as a standalone data-collection platform that works with all of our sensors. This makes it an excellent choice for teachers and students in the classroom and in the field.

LABQ3 \$339

www.vernier.com/labq3



Go Direct Weather System

Easily monitor a wide variety of environmental factors with just one sensor. The included Go Direct Weather Vane accessory is required to report wind direction.

GDX-WTVA \$128

www.vernier.com/gdx-wthr



CASE
Curriculum for Agricultural
Science Education

Vernier is proud to work with CASE, the Curriculum for Agricultural Science Education. CASE is an ambitious project started by the National Council for Agricultural Education in 2007 and is managed by the National Association of Agricultural Educators. It is committed to the goal of improving educational experiences for agriculture students by empowering agriculture teachers.

Visit the CASE website at www.case4learning.org

INVESTIGATION 14

Plant Pigments

After analyzing the absorbance spectrum of chlorophyll from spinach, students investigate the absorbance spectrum of other pigments commonly found in fruits, vegetables, and other plants.

Free sample experiment available at www.vernier.com/plant-pigments



INVESTIGATION 4

Chemistry of Membranes


After measuring how alcohol damages the cell membranes of beets, students investigate how other compounds can damage cell membranes.



INVESTIGATION 6C

Testing Enzyme Activity

Students measure the enzymatic activity of turnip peroxidase and investigate how different factors (e.g., enzyme concentration, substrate concentration, pH, and temperature) influence enzyme activity.



Investigating Biology through Inquiry


Investigating Biology through Inquiry contains investigations for many fundamental concepts in biology. Each investigation includes a preliminary activity, instructor information, sample researchable questions, and sample data.

Topics covered include

- Cell and molecular biology
- Organismal biology
- Ecology
- Evolution

If you are new to inquiry-based instruction, the extensive Instructor Information sections that accompany each investigation help guide you through the inquiry-based style of biology instruction.

Learn more at www.vernier.com/bio-i



INCLUDES
22
INVESTIGATIONS

Download only

BIO-I-E \$40

Printed book + download

BIO-I \$48


Spectrometers

Go Direct SpectroVis Plus

Introduce your students to spectroscopy with our affordable Go Direct® SpectroVis® Plus. Students can easily collect a full-wavelength spectrum (absorbance, percent transmittance, fluorescence, or intensity), study absorbance vs. concentration (standard curve), or monitor enzyme activity (kinetics). Collect and analyze data using Vernier Spectral Analysis® app, LabQuest® App, or Logger Pro® 3.

GDX-SVISPL \$399

www.vernier.com/gdx-svispl



Plant pigments spectra

Vernier UV-VIS Spectrophotometer

The Vernier UV-VIS Spectrophotometer is a portable ultraviolet and visible light spectrophotometer. It is ideal for measuring the absorbance spectra of various biochemical compounds, including DNA and proteins.

VSP-UV \$2,100

www.vernier.com/vsp-uv



Vernier Fluorescence/UV-VIS Spectrophotometer

This spectrophotometer measures the fluorescence and absorbance spectra of ultraviolet and visible samples such as quinine sulfate, fluorescein, rhodamine, and DAPI.

VSP-FUV \$2,899

www.vernier.com/vsp-fuv




Biotechnology

EXPERIMENT 17

Macromolecules: Experiments with Protein

The protein content of milk and protein drinks are measured and analyzed using the Bradford Assay.



Sensor Used



Go Direct SpectroVis Plus

Use Go Direct® SpectroVis® Plus to collect a full-wavelength spectrum, create standard curves for Bradford and other colorimetric assays, or to monitor enzymatic reactions.

GDX-SVISPL \$399

Download free sample experiments at www.vernier.com/bio-rad-kits

Experiment Source



Advanced Biology with Vernier

Download only: BIO-A-E \$40


Printed book + download: BIO-A \$48

Learn more at www.vernier.com/bio-a-17

EXPERIMENT 6B

Forensic DNA Fingerprinting

Students use prepared DNA samples to determine if any of the five "suspects" from a "crime scene" can be excluded as suspects. Gel electrophoresis, DNA staining, and imaging techniques are used to analyze the samples.



Equipment Used



BlueView™ Transilluminator

This uses super bright blue LEDs to illuminate electrophoresis gels stained with fluorescent dyes (e.g., SYBR® Safe). This combination is a safer alternative to ethidium bromide and a UV transilluminator.

BLUE-VIEW 2 \$419

Download free sample experiments at www.vernier.com/bio-rad-kits

Experiment Source



Advanced Biology with Vernier

Download only: BIO-A-E \$40

Printed book + download: BIO-A \$48

Learn more at www.vernier.com/bio-a-6b

Key Products for Biotech

Go Direct Conductivity

GDX-CON \$99



Go Direct Tris-Compatible Flat pH

GDX-FPH \$115



Go Direct Temperature

GDX-TMP \$69



Go Direct Drop Counter

GDX-DC \$99



Stir Station

STIR \$129



OHAUS® Balances

www.vernier.com/ohaus



BIO-RAD

Vernier and Bio-Rad Laboratories

Bio-Rad® combines high-quality supplies, equipment, and curricula with outstanding customer service and technical support—things we believe are important to teachers. Vernier and Bio-Rad enhance classroom experiences with joint experiments and curricula for biotechnology.

Download free sample experiments at www.vernier.com/bio-rad-kits

Imagers



USB Digital Microscope

This 5 megapixel camera connects to a computer or Chromebook™ via USB. It features 10–300× magnification with manual focus and an adjustable LED light source.

BD-EDU-100 \$119

www.vernier.com/bd-edu-100



Celestron Digital Microscope Imagers

Celestron® Digital Microscope Imagers turn your traditional compound or stereo microscope (not included) into a high-resolution digital imager using a personal computer or Chromebook.

CS-5MP ⚡ \$109

CS-DMI ⚡ \$79

www.vernier.com/cs-dmi

Featured Products

Go Direct Sensors

Sensor		Order Code	Price
Go Direct® Blood Pressure		GDX-BP	\$105
Go Direct CO ₂ Gas		GDX-CO2	\$199
Go Direct Colorimeter		GDX-COL	\$119
Go Direct Conductivity		GDX-CON	\$99
Go Direct EKG		GDX-EKG	\$159
Go Direct Ethanol Vapor		GDX-ETOH	\$149
Go Direct Force and Acceleration (for use with Reflex Hammer Accessory Kit)		GDX-FOR	\$99
Go Direct Gas Pressure		GDX-GP	\$89
Go Direct Hand Dynamometer		GDX-HD	\$109
Heart Rate Monitors			
Go Wireless® Exercise Heart Rate		GW-EHR	\$79
Go Wireless Heart Rate		GW-HR	\$89
Go Direct O ₂ Gas		GDX-O2	\$189

Go Direct Optical Dissolved Oxygen



GDX-ODO

\$298

pH Sensors

Go Direct pH



GDX-PH

\$89

Go Direct Tris-Compatible Flat pH



GDX-FPH

\$115

Go Direct Respiration Belt



GDX-RB

\$99

Go Direct Spirometer



GDX-SPR

\$199

Go Direct SpectroVis® Plus



GDX-SVISPL

\$399

Temperature Probes

Go Direct Surface Temperature



GDX-ST

\$79

Go Direct Temperature



GDX-TMP

\$69

Go Direct Weather System



GDX-WTVA

\$128

Accessories

Accessory	Order Code	Price
Go Direct Charge Station	GDX-CRG	\$69
Reflex Hammer Accessory Kit	RFX-ACC	\$29

See all our products for biology at www.vernier.com/biology

LabQuest Sensors

Sensor	Order Code	Price
25-g Accelerometer	ACC-BTA	\$96
Blood Pressure Sensor	BPS-BTA	\$109
CO ₂ Gas Sensor	CO2-BTA	\$269
Colorimeter	COL-BTA	\$119
Conductivity Probe	CON-BTA	\$99
EKG Sensor	EKG-BTA	\$158
Ethanol Sensor	ETH-BTA	\$119
Gas Pressure Sensor	GPS-BTA	\$89
Goniometer	GNM-BTA ⚡	\$159
Hand Dynamometer	HD-BTA	\$110
Heart Rate Monitors		
Exercise Heart Rate Monitor	EHR-BTA	\$99
Hand-Grip Heart Rate Monitor	HGH-BTA	\$119
O ₂ Gas Sensor	O2-BTA	\$199
Optical DO Probe	ODO-BTA	\$299
PAR Sensor	PAR-BTA	\$229
pH Sensors		
pH Sensor	PH-BTA	\$88
Tris-Compatible Flat pH Sensor	FPH-BTA	\$104
Qubit Sensors		
Qubit EKG/EMG Sensor	Q-S207	\$1099
Qubit GSR Sensor	Q-S222	\$899
Soil Moisture Sensor	SMS-BTA	\$109
Spirometer	SPR-BTA	\$219
Temperature Probes		
Stainless Steel Temperature Probe	TMP-BTA	\$36
Surface Temperature Sensor	STS-BTA	\$25

Spectrophotometers

Equipment	Order Code	Price
Go Direct SpectroVis Plus	GDX-SVISPL	\$399
Vernier Fluorescence/UV-VIS Spectrophotometer	VSP-FUV	\$2,899
Vernier UV-VIS Spectrophotometer	VSP-UV	\$2,100

Digital Microscopes

Equipment	Order Code	Price
Celestron® Digital Microscope Imager	CS-DMI ⚡	\$79
5MP Celestron Digital Microscope	CS-5MP ⚡	\$109
USB Digital Microscope	BD-EDU-100	\$119

Lab Books*

Title	Order Code	Price
<i>Biology with Vernier</i>	BWV	\$48
<i>Investigating Biology through Inquiry</i>	BIO-I	\$48
<i>Advanced Biology with Vernier</i> (LabQuest sensors only)	BIO-A	\$48
<i>Human Physiology Experiments: Volume 1</i> (Go Direct sensors only)	HSB-HP	\$38
<i>Human Physiology Experiments: Volume 2</i> (Go Direct sensors only)	ALB-HP2	\$38
<i>Human Physiology with Vernier</i> (LabQuest® sensors only)	HP-A	\$48
<i>Agricultural Science with Vernier</i> (LabQuest sensors only)	AWV	\$48

* Includes printed book and download; also available as a download only

Looking for Replacement Parts?

Visit www.vernier.com/replacements

See all our products for biology at www.vernier.com/biology

Environmental Science

www.vernier.com/environmental-science

Help your students see that the environmental science concepts discussed in the classroom have serious implications on the world around them. Our hands-on investigations and data-collection technology help students form a better understanding of phenomena.

Topics

Explore a sampling of our featured experiments by topic to learn how Vernier technology helps your students engage with data-collection technology and deepens their understanding of key environmental science concepts.

Environmental Science

page 60

Water Quality

page 62


Renewable Energy

page 64



Show Students How To Investigate Their World

From soil studies to wind energy investigations, the study of environmental science helps students understand how to interact with the natural world. Our easy-to-use sensors support you as you help your students understand key environmental science concepts. Our lab books include ready-to-go investigations to help students establish a deep understanding of key scientific concepts.



Professional Development


We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training

INVESTIGATION 31

The Effect of Acid Deposition on Aquatic Ecosystems

Investigate acid deposition by measuring the magnitude of the change in pH levels in an aquatic environment when dilute acid is introduced dropwise.



Sensors Used



Go Direct® Tris-Compatible Flat pH

The flat glass, double-junction design makes this sensor a good choice for environmental science.

GDX-FPH \$115

Accessories Used



Electrode Support

ESUP \$10



Go Direct Conductivity

Determine the ionic content of an aqueous solution by measuring its electrical conductivity.

GDX-CON \$99

Stir Station

STIR \$129



Investigation Source



Investigating Environmental Science through Inquiry

Download only: ESI-E \$40

Printed book + download: ESI \$48

Learn more at www.vernier.com/esi-31

INVESTIGATION 26

Fossil Fuel Energy

Students calculate the amount of heat transferred from a burning candle to a known volume of water. They also design an experiment to investigate fossil fuels.



Sensor Used



Go Direct Temperature

This is a rugged, general-purpose sensor that students can use to monitor temperature.

GDX-TMP \$69

Investigation Source



Investigating Environmental Science through Inquiry

Download only: ESI-E \$40

Printed book + download: ESI \$48


Learn more at vernier.com/esi-26

INCLUDES
34
INVESTIGATIONS

INVESTIGATION 3

Investigating Dissolved Oxygen

Students analyze the effect temperature has on dissolved oxygen in water by measuring the concentration of dissolved oxygen in different temperatures of water.



Sensor Used



Go Direct Optical Dissolved Oxygen

This optical sensor makes it easy to measure dissolved oxygen in water, atmospheric pressure, and water temperature.

GDx-ODO \$298

Investigation Source



Investigating Environmental Science through Inquiry

Download only: ESI-E \$40

Printed book + download: ESI \$48

Learn more at www.vernier.com/esi-3

Investigating Environmental Science through Inquiry

Investigating Environmental Science through Inquiry contains 34 inquiry-based environmental science investigations.

Topics include

- Earth systems and resources (air, water, and soil)
- The living world
- Global change and population
- Energy resources and consumption
- Pollution

Learn more at www.vernier.com/esi

** Instructions for Vernier Graphical Analysis not yet available*



Download only

ESI-E \$40

Printed book + download

ESI \$48

Environmental Science Go Direct Starter Package

This package includes four sensors that all work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3.

- Go Direct Temperature
- Go Direct Tris-Compatible Flat pH
- Go Direct Conductivity
- Go Direct Optical Dissolved Oxygen

GDP-EV-ST \$581

Learn more at

www.vernier.com/gdp-ev-st



Water Quality

TEST 12

Total Dissolved Solids

Students measure the total dissolved solids of a sample from a local body of fresh water.



Sensor Used



Go Direct® Conductivity

Determine the ionic content of an aqueous solution by measuring its electrical conductivity.

GDX-CON \$99

Accessories Used



Water Quality Bottles

This box of 8 plastic bottles with stoppers is for general water quality use. They could also be used as replacements for the bottles and stoppers in the Primary Productivity Kit. See page 46.

WQ-BOT \$28

Experiment Source



Water Quality with Vernier

Download only: WQV-E \$40

Printed book + download: WQV \$48

Learn more at www.vernier.com/wqv-12

LabQuest 3

LabQuest 3 is a powerful, connected, and remarkably versatile data-logging solution.

Why? LabQuest® 3 can serve as a standalone data-collection platform that works with all of our sensors. This makes it an excellent choice for teachers and students in the classroom and in the field.

LABQ3 \$339

www.vernier.com/labq3



Go Direct Sensor Clamp

The Go Direct Sensor Clamp securely fastens to a wand-style Go Direct sensor, and the included lanyard works as a strap to prevent accidental drops during investigations in the field. Sensors are sold separately.

GDX-CLAMP \$12

Learn more at

www.vernier.com/gdx-clamp



GLOBE[®] & Vernier

The GLOBE Program is an international science and education program that provides students and the public worldwide with the opportunity to participate in data collection and the scientific process as well as contribute meaningfully to our understanding of the Earth system and global environment. Use Vernier sensors to collect GLOBE data.

To learn more about Vernier and GLOBE, see www.vernier.com/globe



Weather

NEW Go Direct Weather System

Easily monitor a wide variety of environmental factors with just one sensor. Go Direct Weather System includes an affordable, wireless handheld sensor used to measure ambient temperature, humidity, wind speed, and more. The included Go Direct Weather Vane accessory is required to report wind direction.

GDX-WTVA (sensor and vane) \$128

Learn more at www.vernier.com/gdx-wtva



Davis Vantage Vue Weather Station

The wireless Vantage Vue Weather Station provides accurate, reliable weather monitoring in a self-contained, easy-to-install system. The sensor suite measures

- Temperature
- Humidity
- Barometric pressure
- Wind speed and direction
- Dew point
- Rainfall



Available Bundles	Stream Live Data on the Internet via Wi-Fi	View Data on Console	Order Code	Price
Davis [®] Vantage Vue [®] Wireless Weather Station (with console)		●	DWVUE	\$395
Davis Vantage Vue + WeatherLink [™] (without console)	●		DWVUE-LWOC	\$465
Davis Vantage Vue + WeatherLink (with console)	●	●	DWVUE-LWC	\$595

For accessories and weather station options, visit www.vernier.com/weather

Renewable Energy



Strengthen students' critical thinking skills by introducing them to alternative energy solutions to real-world problems.

The KidWind Project and Vernier have teamed up to provide the technology, resources, and support you need for your students to investigate renewable energy.

- Engage your students as they watch power output and energy production data develop in real time.
- Inspire creativity as your students build and test prototypes, test solutions to engineering problems, and optimize designs.
- Measure voltage and current, and calculate power, without using a multimeter.
- Set up activities quickly and easily, creating more time for instruction and exploration.

Recommended Classroom Setup for Wind Energy



3 Test Stations



6 to 10 Groups
of 2 to 4 Students

We recommend three test stations for a classroom with 6 to 10 groups of 2 to 4 students.

Each test station should have

- Box fan
- Wind turbine tower with nacelle and generator
- Go Direct® Energy (GDX-NRG)
- Vernier Variable Load (VES-VL)

Each student group needs

- Blade Pitch Protractor
- Wind Turbine Hub
- Blade consumables

KidWind Accessories & Replacement Parts

Part Name	Order Code	Price
Balsa Blade Sheets (100 Pack)	KW-BBS10	\$12
Basic Turbine Building Parts	KW-BTPART	\$16
Blade Design Consumables Classroom Pack	KW-BDC	\$149
Blade Pitch Protractor	KW-BPP	\$4
Chipboard Sheets (50 Pack)	KW-CB50	\$19
Dowels (25 Pack)	KW-D25	\$5
Dowels (100 Pack)	KW-D100	\$10
Gear Set	KW-GEAR	\$9
High Torque Generator with Wires	KW-HIGEN	\$9
KidWind Airfoil Balsa Blade Sheets	KW-ABBS10	\$19
Power Output Board	KW-POBD	\$39
Wind Turbine Generator (10 Pack)	KW-GEN10	\$60
Wind Turbine Hub (3 Pack)	KW-WTH3	\$22

Learn more at www.vernier.com/renewable-energy


Featured Experiments

26 Experiments Available

EXPERIMENT 24

Exploring Solar Collectors

Students measure the temperature change produced when using a solar collector. Students then design an experiment to test the impact a changed variable has on a solar collector.



Sensors Used



Go Direct Surface Temperature

Use this sensor in situations in which low thermal mass or flexibility is required.

GDX-ST \$79



Go Direct Light and Color

Students use this sensor to measure the brightness of a light bulb or the reflectance of light off of various objects.

GDX-LC \$79


Accessory Used



Solar Thermal Exploration Kit

KW-STXK \$59

Experiment Source



Renewable Energy with Vernier

Download only: REV-E \$40

Printed book + download: REV \$48

Learn more at www.vernier.com/rev-24

EXPERIMENT 17

Exploring Solar Panels

Investigate different variables and how they impact electricity production with a solar panel. Students also calculate the efficiency of power production with the solar panel.



Sensors Used



Go Direct Energy

This sensor quantifies the voltage, current, power, and energy output of small wind turbines and solar panels, such as those used in our KidWind Experiment Kits.

GDX-NRG \$89



Go Direct Light and Color

Students use this sensor to measure the brightness of a light bulb or the reflectance of light off of various objects.

GDX-LC \$79

Accessories Used



KidWind 2V/400mA Solar Panel


KW-SP2V \$19



Vernier Variable Load

VES-VL \$64

Experiment Source



Renewable Energy with Vernier

Download only: REV-E \$40

Printed book + download: REV \$48

Learn more at www.vernier.com/rev-17

Featured Experiments

EXPERIMENT 8

Exploring Wind Turbines

Students investigate different variables that affect how a wind turbine moves and produces electricity.



Sensor Used



Go Direct® Energy

This sensor quantifies the voltage, current, power, and energy output of small wind turbines and solar panels, such as those used in our KidWind Experiment Kits.

GDx-NRG \$89

Accessories Used



KidWind Advanced Wind Experiment Kit


KW-AWX \$154

Vernier Variable Load

VES-VL \$64



Experiment Source



Renewable Energy with Vernier

Download only: REV-E \$40

Printed book + download: REV \$48


Learn more at www.vernier.com/rev-8

Renewable Energy with Vernier

The *Renewable Energy with Vernier* lab book features 26 experiments in wind and solar energy. The book contains a combination of explorations, classic experiments, inquiry investigations, engineering projects, and more.

Learn more at www.vernier.com/rev

INCLUDES
26
EXPERIMENTS



Download only

REV-E \$40

Printed book + download

REV \$48

KidWind Competitions—Putting the “E” in STEM

Challenge students to compete in a wind turbine design competition with peers in a supportive environment at local and national events.

To see our recommendations and to get started, visit

www.vernier.com/kidwind-challenges



Featured Products

KidWind Advanced Wind Experiment Kit


Discover advanced concepts of wind turbine technology, including gearboxes and generator construction (with the optional KidWind simpleGEN). Students use the blades they design to generate electricity, lift weights, and pump water. This kit is recommended for use with our lab book *Renewable Energy with Vernier*.

KW-AWX \$154

KidWind Advanced Wind Experiment Kit Classroom Pack

KW-AWXC \$389

Learn more at www.vernier.com/kw-awx



KidWind simpleGEN

The simpleGEN is an easy-to-build AC generator that students can use to demonstrate Faraday's law, light LEDs, and perform experiments that explore how coils, magnets, and rotation affect power generation.

KW-SGEN \$59

Learn more at www.vernier.com/kw-sgen



Solar Energy Exploration Kit

Explore solar energy with this innovative science kit designed to help students investigate energy transformations. Experiment with basic circuits and learn about important factors in photovoltaic systems.

KW-SEEK \$79

Learn more at www.vernier.com/kw-seek



KidWind GENPack

Using the parts in the GENPack, students can construct their own electrical generator and perform experiments with electricity and magnetism. Changing variables in the generator design affects current and voltage output.







KW-GP \$54







Learn more at www.vernier.com/kw-gp



Featured Products

Go Direct Sensors

Sensor		Order Code	Price
Go Direct® CO ₂ Gas		GDX-CO2	\$199
Go Direct Colorimeter		GDX-COL	\$119
Go Direct Conductivity		GDX-CON	\$99
Go Direct Current		GDX-CUR	\$79
Go Direct Energy		GDX-NRG	\$89
Go Direct Ethanol Vapor		GDX-ETOH	\$149
Go Direct Light and Color		GDX-LC	\$79
Ion-Selective Electrodes			
Go Direct Ammonium Ion-Selective Electrode		GDX-NH4	\$249
Go Direct Calcium Ion-Selective Electrode		GDX-CA	\$249
Go Direct Chloride Ion-Selective Electrode		GDX-CL	\$249
Go Direct Nitrate Ion-Selective Electrode		GDX-NO3	\$249









Go Direct O ₂ Gas		GDX-O2	\$189
Go Direct Optical Dissolved Oxygen		GDX-ODO	\$289
pH Sensors			
Go Direct pH		GDX-PH	\$89
Go Direct Tris-Compatible Flat pH		GDX-FPH	\$115
Go Direct SpectroVis® Plus		GDX-SVISPL	\$399
Temperature Probes			
Go Direct Surface Temperature		GDX-ST	\$79
Go Direct Temperature		GDX-TMP	\$69
Go Direct Voltage		GDX-VOLT	\$69
Go Direct Weather System		GDX-WTVA	\$128

See all our products for environmental science at www.vernier.com/environmental-science

Go Direct Accessories

Accessory		Order Code	Price
Go Direct Charge Station		GDX-CRG	\$69
Go Direct Sensor Clamp		GDX-CLAMP	\$12

LabQuest Sensors

Sensor		Order Code	Price
Conductivity Probe		CON-BTA	\$99
Flow Rate Sensor		FLO-BTA	\$129
Optical DO Probe		ODO-BTA	\$299
pH Sensor		PH-BTA	\$88
Tris-Compatible Flat pH		FPH-BTA	\$104
Salinity Sensor		SAL-BTA	\$119
Soil Moisture Sensor		SMS-BTA	\$109
Turbidity Sensor		TRB-BTA	\$112

Digital Microscopes

Equipment	Order Code	Price
Celestron® Digital Microscope Imager	CS-DMI ⚡	\$79
USB Digital Microscope	BD-EDU-100	\$119

Lab Equipment

Equipment	Order Code	Price
KidWind Advanced Wind Energy Kit	KW-AWX	\$154
KidWind Basic Wind Energy Kit	KW-BWX	\$124
Primary Productivity Kit	PPK	\$44
Solar Energy Exploration Kit	KW-SEEK	\$79
Water Depth Sampler	WDS	\$89
Water Quality Bottles	WQ-BOT	\$28

Lab Books

Book Title	Order Code	Price
<i>Investigating Environmental Science through Inquiry</i>	Printed book + download: ESI Download only: ESI-E	\$48 \$40
<i>Water Quality with Vernier</i>	Printed book + download: WQV Download only: WQV-E	\$48 \$40
<i>Renewable Energy with Vernier</i>	Printed book + download: REV Download only: REV-E	\$48 \$40
<i>Climate and Meteorology Experiments</i>	Download only: HSB-CM-E	\$20

Looking for Replacement Parts?

Visit www.vernier.com/replacements

See all our products for environmental science at www.vernier.com/environmental-science

Earth Science

www.vernier.com/earth-science

When you use Vernier technology to teach Earth science you can count on our affordable sensors, intuitive software, and creative solutions to help your students understand key Earth science concepts.



Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training



Earth Science Helps Students Understand Their World


The study of Earth science helps you give students a means to understand the world around them. Your students can explore sea floor spreading, the effect of acid rain on soil, the changing of the seasons, and more with Vernier sensors, software, and experiments.

Weather and Climate

EXPERIMENT 4

Greenhouse Effect

Students use temperature probes to measure temperatures in a model greenhouse, then they analyze collected data to make conclusions about the greenhouse effect.



Sensor Used




Go Direct® Surface Temperature

This sensor has an exposed thermistor that results in an extremely rapid response time, making it perfect for use in air and water.

GDX-ST \$79

Experiment Source



Climate and Meteorology Experiments

Download only: HSB-CM-E \$20

Learn more at www.vernier.com/hsb-cm-e

NEW Climate and Meteorology Experiments

This new lab book is packed with interactive investigations that challenge students to use data-collection technology to explore storm systems and other important weather-related topics.

The experiments in this e-book cover

- The Greenhouse Effect
- Dew Point
- Microclimates

Climate and Meteorology Experiments Go Direct Package

(includes all the sensors needed to do the activities in the book)

- Go Direct Surface Temperature (2)
- Go Direct Light and Color
- Go Direct Weather System


GDP-CM \$365

Buy 8 or more packages at \$354 and save \$88

Learn more at www.vernier.com/hsb-cm-e

11
EXPERIMENTS
INCLUDED IN
E-BOOK

Climate and Meteorology Experiments



Download only

HSB-CM-E \$20

NEW Go Direct Weather System

Easily monitor a wide variety of environmental factors with just one sensor. Go Direct Weather System includes an affordable, wireless handheld sensor used to measure ambient temperature, humidity, wind speed, and more. The included Go Direct Weather Vane accessory is required to report wind direction.

GDX-WTVA (sensor and vane) \$128

Learn more at www.vernier.com/gdx-wtva



Earth Science

EXPERIMENT 29

Seasons and Angle of Insolation

In this experiment, students model how the angle of light from the sun striking various places on Earth is one factor that causes seasons.



Sensor Used



Go Direct® Temperature

This rugged probe measures the temperature of a variety of substances including air, soil, and water.

GDX-TMP \$69

Experiment Source



Earth Science with Vernier

Download only: ESV-E \$40

Printed book + download: ESV \$48

Learn more at [vernier.com/esv-29](http://www.vernier.com/esv-29)

Earth Science with Vernier


In addition to the 33 experiments in *Earth Science with Vernier*, the six projects in this book engage students as they learn about the world around them.

Topics include

- Geology
- Soil analysis
- Water quality tests
- Hydrology/Oceanography
- Meteorology
- Energy

Learn more at [vernier.com/esv](http://www.vernier.com/esv)

INCLUDES
33
EXPERIMENTS



Download only

ESV-E \$40

Printed book + download

ESV \$48

Go Direct 3-Axis Magnetic Field

Useful for topics in geology, this sensor can determine the magnitude and direction of a magnetic field at any point in space.

GDX-3MG \$69

Learn more at www.vernier.com/gdx-3mg



Featured Products

Looking for Replacement Parts?

Visit www.vernier.com/replacements

Go Direct Sensors

Sensor	Order Code	Price
Go Direct 3-Axis Magnetic Field	GDX-3MG	\$69
Go Direct CO ₂ Gas	GDX-CO2	\$199
Go Direct Conductivity	GDX-CON	\$99
Go Direct Current	GDX-CUR	\$79
Go Direct Energy	GDX-NRG	\$89
Go Direct Light and Color	GDX-LC	\$79
Go Direct Motion	GDX-MD	\$99
Go Direct O ₂ Gas	GDX-O2	\$189
Go Direct Optical Dissolved Oxygen	GDX-ODO	\$289
pH Sensors		
Go Direct pH	GDX-PH	\$89
Go Direct Tris-Compatible Flat pH	GDX-FPH	\$115
Temperature Probes		
Go Direct Surface Temperature	GDX-ST	\$79
Go Direct Temperature	GDX-TMP	\$69
Go Direct Voltage	GDX-VOLT	\$69
Go Direct Weather	GDX-WTHR	\$99
Go Direct Weather System	GDX-WTVA	\$128

Go Direct Accessories

Accessory	Order Code	Price
Go Direct Charge Station	GDX-CRG	\$69
Go Direct Sensor Clamp	GDX-CLAMP	\$12

LabQuest Sensors

Sensor	Order Code	Price
Anemometer	ANM-BTA	\$89
Barometer	BAR-BTA	\$71
Flow Rate Sensor	FLO-BTA	\$129
Magnetic Field Sensor	MG-BTA	\$58
Salinity Sensor	SAL-BTA	\$119
Soil Moisture Sensor	SMS-BTA	\$109
Stainless Steel Temperature Probe	TMP-BTA	\$36
Tris-Compatible Flat pH Sensor	FPH-BTA	\$104
Turbidity Sensor	TRB-BTA	\$112

Accessories & Lab Equipment

Product	Order Code	Price
Davis® Weather Stations	vernier.com/weather	
Electrode Support	ESUP	\$10
KidWind 2V/400mA Solar Panel	KW-SP2V	\$19
KidWind Basic Wind Experiment Kit	KW-BWX	\$124
Solar Energy Exploration Kit	KW-SEEK	\$79
Vernier Resistor Board	VES-RB	\$18

Lab Books

Title	Order Code	Price
<i>Earth Science with Vernier</i>	Printed book + download: ESV	\$48
	Download only: ESV-E	\$40
<i>Water Quality with Vernier</i>	Printed book + download: WQV	\$48
	Download only: WQV-E	\$40
<i>Climate and Meteorology Experiments</i>	Download only: HSB-CM-E	\$20

See all our products for Earth science at www.vernier.com/earth-science

Chemistry

www.vernier.com/chemistry

Vernier chemistry resources cover an array of key concepts to help prepare your students for what lies ahead. From gas laws to spectroscopy, our products are backed by an extensive collection of experiments and unparalleled technical support.



Topics

Explore a sampling of our featured experiments by topic to learn how Vernier technology helps your students engage with data-collection technology and deepens their understanding of key chemistry concepts.

General Chemistry

PAGE 76

AP* Chemistry

PAGE 78

Advanced Chemistry

PAGE 80

Inquiry Chemistry

PAGE 82

Food Chemistry

PAGE 83

Organic Chemistry

PAGE 89



Make Your Chemistry Classes More Elemental

Whether you are teaching Beer's law or exploring how humans use food for energy, Vernier technology and investigations help your students better understand important chemistry concepts. Give your students insight into this vital subject with interactive learning opportunities from Vernier.


Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training

*AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

For more information, and to see all our products, visit www.vernier.com



EXPERIMENT 2

Freezing and Melting of Water

Students measure the temperature of water as it changes from a liquid to a solid. The data are analyzed to make predictions about the freezing patterns of other substances.



Sensor Used




Go Direct® Temperature

Students can use this rugged, general-purpose sensor to monitor temperature.

Range: -40 to 125°C

GDX-TMP \$69

Experiment Source



Chemistry with Vernier

Download only: CWV-E \$40


Printed book + download: CWV \$48

Learn more at www.vernier.com/cwv-2

EXPERIMENT 6

Boyle's Law: Pressure-Volume Relationship in Gases

Determine the mathematical relationship between pressure and volume of a gas.



Sensor Used




Go Direct Gas Pressure

Explore pressure changes and gas laws with this sensor that measures the absolute pressure of a gas.

GDX-GP \$89

Experiment Source



Chemistry with Vernier

Download only: CWV-E \$40

Printed book + download: CWV \$48

Learn more at www.vernier.com/cwv-6

EXPERIMENT 21

Household Acids and Bases

Students investigate the pH scale by measuring the pH of household solutions using different methods.



Sensor Used



Go Direct pH

This general-purpose pH sensor is used to monitor pH of aqueous solutions.

GDX-PH \$89

Accessories Used



Electrode Support


ESUP \$10



Stir Station

STIR \$129

Experiment Source



Chemistry with Vernier

Download only: CWV-E \$40

Printed book + download: CWV \$48

Learn more at www.vernier.com/cwv-21

Chemistry with Vernier


Combine *Chemistry with Vernier* with the Starter Package (shown below) to teach students the essentials in chemistry. This lab book contains ready-to-use student experiments and instructor information, including sample data.

Topics include

- Thermochemistry
- Gas laws
- Acid-base reactions
- Equilibrium
- Electrochemistry
- Electrolytes
- States of matter

Learn more at www.vernier.com/cwv

INCLUDES
36
EXPERIMENTS



Download only

CWV-E \$40

Printed book + download

CWV \$48

Chemistry Go Direct Starter Package

This package includes four sensors that all work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3.

- Go Direct Temperature (2)
- Go Direct Gas Pressure
- Go Direct pH

GDP-CH-ST \$316

Learn more at www.vernier.com/gdp-ch-st


Standard package also available
(see page 81)



INVESTIGATION 1

Investigating Food Dyes in Sports Beverages

Use spectroscopy to examine the relationship between % transmittance and concentration of a solution to determine the amount of food dye in a sports drink.



Sensor Used



Go Direct[®] SpectroVis[®] Plus

This spectrophotometer quickly measures a full-wavelength spectrum (380 to 950 nm).

GDx-SVISPL \$399

Recommended Accessories



100 Plastic Cuvettes (Visible Range)

CUV \$19



Cuvette Rack

CUV-RACK \$9

Investigation Source



Vernier Chemistry Investigations for Use with AP^{*} Chemistry

Download only: APCHEM-E \$40


Printed book + download: APCHEM \$48

Learn more at www.vernier.com/apchem-1

INVESTIGATION 8

Determining the Percent Hydrogen Peroxide in a Commercial Product

Test a bottle of commercial hydrogen peroxide and determine the concentration using a potentiometric titration.



Sensors Used



Go Direct ORP

Measure the ability of a solution to act as an oxidizing or reducing agent.

GDx-ORP \$99

Accessory Used



Go Direct Drop Counter

As an alternative to using a buret, the drop counter precisely records the number of drops of titrant added during a titration and then automatically converts it to volume.

GDx-DC \$99



Stir Station

STIR \$129

Investigation Source



Vernier Chemistry Investigations for Use with AP^{*} Chemistry

Download only: APCHEM-E \$40

Printed book + download: APCHEM \$48

Learn more at www.vernier.com/apchem-8

INVESTIGATION 9

Investigating the Components of a Commercial Tablet

A pain medication tablet chips and cracks due to contamination or an incorrect tablet formula. Students use melting point to investigate these two theories.



Sensor Used



Go Direct Melt Station

Accurately determine the melting temperature of solid substances.

GDx-MLT  \$529

Recommended Accessory



Melt Station Capillary Tubes

MLT-TUBE \$19

Investigation Source



Vernier Chemistry Investigations for Use with AP* Chemistry

Download only: APCHEM-E \$40

Printed book + download: APCHEM \$48

Learn more at www.vernier.com/apchem-9

Vernier Chemistry Investigations for Use with AP* Chemistry


This lab book provides AP* Chemistry students with 16 inquiry-based laboratory experiments aligned with the investigations published by the College Board.

Topics include

- Spectroscopy
- Titrations
- Intermolecular forces and properties

Learn more at www.vernier.com/apchem

INCLUDES
16
INVESTIGATIONS




Download only

APCHEM-E \$40

Printed book + download

APCHEM \$48

Chemistry Lab Books with AP* Correlations




Vernier Chemistry Investigations for Use with AP* Chemistry

16 Investigations

Download only: APCHEM-E \$40

Printed book + download: APCHEM \$48



Advanced Chemistry with Vernier

35 Experiments

Download only: CHEM-A-E \$40

Printed book + download: CHEM-A \$48



Investigating Chemistry through Inquiry

25 Investigations

Download only: CHEM-I-E \$40

Printed book + download: CHEM-I \$48

To see all AP correlations, visit www.vernier.com/ap-correlations

* AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

EXPERIMENT 10

The Determination of an Equilibrium Constant

Determine the concentration of ions present in an equilibrium system using spectroscopy. Students calculate the equilibrium constant, K_{eq} , for the reaction.



Sensor Used



Go Direct® SpectroVis® Plus

This spectrophotometer quickly measures a full-wavelength spectrum (380 to 950 nm).

GDX-SVISPL \$399

Recommended Accessories



100 Plastic Cuvettes (Visible Range)

CUV \$19



Cuvette Rack

CUV-RACK \$9

Experiment Source



Advanced Chemistry with Vernier

Download only: CHEM-A-E \$40


Printed book + download: CHEM-A \$48

Learn more at www.vernier.com/chem-a-10

EXPERIMENT 30

Exploring the Properties of Gases

Students conduct a set of experiments, each of which illustrates a gas law such as Boyle's law, shown here. Use the results to derive a single mathematical relationship that relates pressure, volume, temperature, and number of molecules.



Sensors Used



Go Direct Gas Pressure

Explore pressure changes and gas laws with this sensor that measures the absolute pressure of a gas.

GDX-GP \$89



Go Direct Temperature

Students can use this rugged, general-purpose sensor to monitor temperature.

Range: -40 to 125°C

GDX-TMP \$69

Accessories Used



Electrode Support


ESUP \$10

Stir Station

STIR \$129



Experiment Source



Advanced Chemistry with Vernier

Download only: CHEM-A-E \$40


Printed book + download: CHEM-A \$48

Learn more at www.vernier.com/chem-a-30

EXPERIMENT 20

Electrochemistry: Voltaic Cells

Construct voltaic cells to explore oxidation-reduction reactions. Use the measured potentials to identify unknown metal electrodes and create concentration cells for understanding the Nernst equation.



Sensor Used




Go Direct Voltage

This sensor has a wide input voltage and high precision, making it an excellent choice for investigating the basic principles of electrochemical cells.

Range: ± 20 V

GDX-VOLT \$69

Experiment Source



Advanced Chemistry with Vernier

Download only: CHEM-A-E \$40

Printed book + download: CHEM-A \$48

Learn more at www.vernier.com/chem-a-20

Advanced Chemistry with Vernier


The *Advanced Chemistry with Vernier* lab book expands students' skills with experiments appropriate for second year, honors, and AP* Chemistry students.

Topics include

- Redox reactions
- Colligative properties
- Equilibrium

Learn more at www.vernier.com/chem-a

INCLUDES
35
EXPERIMENTS



Download only

CHEM-A-E \$40

Printed book + download

CHEM-A \$48

*AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

Chemistry Go Direct Standard Package

This package includes 8 sensors that all work with our free Vernier Graphical Analysis™ app, as well as Graphical Analysis Pro and LabQuest® 3.

- Go Direct Temperature (2)
- Go Direct Conductivity
- Go Direct Gas Pressure
- Go Direct Colorimeter
- Go Direct pH
- Go Direct Drop Counter
- Go Direct Voltage

GDP-CH-DX \$702

Buy 8 or more packages at \$681 and save \$168


Learn more at www.vernier.com/gdp-ch-dx



INVESTIGATION 8

Evaporation and Intermolecular Attractions

Students study temperature changes caused by the evaporation of different liquids and relate the temperature changes to the strength of intermolecular forces of attraction.



Sensor Used



Go Direct® Temperature

Students can use this rugged, general-purpose sensor to monitor temperature.

Range: -40 to 125°C

GDX-TMP \$69

Investigation Source



Investigating Chemistry through Inquiry

Download only: CHEM-I-E \$40

Printed book + download: CHEM-I \$48

Learn more at www.vernier.com/chem-i-8

Investigating Chemistry through Inquiry

The *Investigating Chemistry through Inquiry* lab book supports both open and guided inquiry experiments. Instructors can help students devise their own researchable questions or choose from a list provided in each experiment.

Topics include

- Chemical kinetics
- Acids and bases
- Thermochemistry

Learn more at www.vernier.com/chem-i

INCLUDES
25
INVESTIGATIONS




Download only

CHEM-I-E \$40

Printed book + download

CHEM-I \$48

Chemistry Lab Books with IB[†] Correlation




Advanced Chemistry with Vernier

Download only: CHEM-A-E \$40

Printed book + download: CHEM-A \$48

35 Experiments



Investigating Chemistry through Inquiry

Download only: CHEM-I-E \$40

Printed book + download: CHEM-I \$48

25 Investigations


To see all IB correlations, visit www.vernier.com/ib-correlations

[†] The IB Diploma Program is an official program of the International Baccalaureate Organization (IBO) which authorizes schools to offer it. The material available here has been developed independently of the IBO and is not endorsed by it.

EXPERIMENT 1

What's the Difference Between Baking Soda and Baking Powder?

Using data-collection technology, students examine the chemical changes that occur when water is added to baking soda and baking powder.



Sensor Used



Go Direct pH

This wireless sensor monitors the pH of aqueous solutions and is perfect for lab and field experiments alike.

GDX-PH \$89



Go Direct CO₂ Gas

Go Direct CO₂ Gas measures gaseous carbon dioxide concentration levels, air temperature, and relative humidity.

GDX-CO2 \$199

Investigation Source



Food Chemistry Experiments

Download only: HSB-FOOD-E \$30

Printed book + download: HSB-FOOD \$38

Learn more at www.vernier.com/hsb-food

NEW Food Chemistry Experiments

This new lab book is filled with experiments that use food as a means to explore crucial chemistry concepts. Students are more likely to engage with science when they see concepts applied to the real world. These experiments use Vernier sensors such as spectrophotometers, temperature probes, and CO₂ gas sensors to investigate complex questions involving food.

Learn more at www.vernier.com/hsb-food

INCLUDES
14
EXPERIMENTS



Download only

HSB-FOOD-E \$30

Printed book + download

HSB-FOOD \$38

Key Products for Food Chemistry Experiments



Go Direct SpectroVis[®] Plus

GDX-SVISPL \$399



Go Direct Polarimeter

GDX-POL \$499



Go Direct Gas Pressure

GDX-GP \$89



Go Direct Conductivity

GDX-CON \$99



Go Direct Temperature

GDX-TMP \$69



Go Direct Ethanol Vapor

GDX-ETOH \$149



Go Direct ORP

GDX-ORP \$99

Chemistry Go Direct Starter Package

4 Sensors • GDP-CH-ST • \$316



This package includes

Go Direct®
Temperature (2)

Go Direct
Gas Pressure

Go Direct
pH


All sensors work with our free Vernier Graphical Analysis™ app,
as well as Graphical Analysis Pro and LabQuest® 3.

Learn more at www.vernier.com/gdp-ch-st

Chemistry Go Direct Standard Package

8 Sensors • GDP-CH-DX • \$702

Buy 8 or more packages at \$681 and save \$168



This package includes

Go Direct
Temperature
(2)

Go Direct
Gas Pressure

Go Direct
pH

Go Direct
Voltage

Go Direct
Conductivity

Go Direct
Colorimeter




Go Direct
Drop Counter

All sensors work with our free Vernier Graphical Analysis app,
as well as Graphical Analysis Pro and LabQuest 3.

Learn more at www.vernier.com/gdp-ch-dx





Featured Products

pH Sensor Comparison

Sensor	Features
Go Direct pH GDX-PH \$89 	Recommended for General Use Go Direct pH is an important and versatile sensor for lab and field activities alike. Conduct acid-base titrations, monitor pH changes during chemical reactions, and investigate household acids and bases. The wireless connection makes it easier to do field-based studies such as testing the pH of surface water.
	Go Direct pH Teacher Pack GDX-PH-TP \$758 Includes 8 Go Direct pH Sensors and a Go Direct Charge Station
Go Direct Tris-Compatible Flat pH GDX-FPH \$115 	Go Direct Tris-Compatible Flat pH is a double-junction electrode for measuring pH in Tris buffers and solutions containing proteins or sulfides. The flat glass shape makes it easy to clean and useful for measuring the pH of semisolids such as soil slurries and certain foods.
Go Direct Glass-Body pH GDX-GPH \$139 	Go Direct Glass-Body pH can be used with non-aqueous solutions and solutions containing solvents, strong acids, and strong bases.

Learn more at www.vernier.com/ph-sensors

Temperature Sensor Comparison

Sensor	Features and Applications
Go Direct Temperature GDX-TMP \$69 Range -40 to 125°C 	Recommended for General Use <ul style="list-style-type: none"> Conduct endothermic and exothermic reactions. Determine the physical properties of water. Measure the energy content of foods. Investigate intermolecular forces.
	Go Direct Temperature Teacher Pack GDX-TMP-TP \$599 Includes 8 Go Direct Temperature Probes and a Go Direct Charge Station
Go Direct Surface Temperature GDX-ST \$79 Range -25 to 125°C 	<ul style="list-style-type: none"> Use this sensor in situations in which low thermal mass or flexibility is required. The exposed thermistor provides an extremely rapid response to temperature changes. Use this sensor in air and water only.
Go Direct Wide-Range Temperature GDX-WRT \$114 Range -20 to 330°C 	<ul style="list-style-type: none"> Determine the melting point of caffeine or the boiling point of different vegetable oils. RTD (Resistance Temperature Detector) technology establishes a $\pm 0.5^\circ\text{C}$ accuracy.
NEW Go Direct Thermocouple GDX-TC \$109 Range (type K) -200 to +1,400°C 	<ul style="list-style-type: none"> Collect reliable data during experiments in which there are extreme temperatures, such as making ice cream with dry ice or testing different elements of a flame. Compatible with Type-K (included), Type-T, and Type-J thermocouple wires

Learn more at www.vernier.com/temperature-sensors

Featured Products

Go Direct Constant Current System

Determine Avogadro's number and perform various electroplating and electrolysis experiments. This system combines a DC power source with a built-in current sensor to eliminate the need for a separate power supply. It can deliver up to 0.6 A at 5 V DC.

GDX-CCS \$74

www.vernier.com/gdx-ccs



Go Direct Melt Station

Teach students the visual detection capillary method of melting point determination with Go Direct® Melt Station. It accurately measures melting temperatures of a solid (up to 260°C), and the real-time graphing provides a unique perspective of the melting process.

GDX-MLT \$529

www.vernier.com/gdx-mlt




Go Direct SpectroVis Plus

Introduce your students to spectroscopy with the affordable Go Direct SpectroVis® Plus Spectrophotometer. With a range of 380 to 950 nm, students can easily collect a full-wavelength spectrum (absorbance, percent transmittance, fluorescence, or intensity), study absorbance vs. concentration (Beer's law), or monitor rates of reaction (kinetics). Collect and analyze data using Vernier Spectral Analysis®, LabQuest® App, or Logger Pro® 3.

GDX-SVISPL \$399

www.vernier.com/gdx-svispl



Absorbance spectra of green food coloring at different concentrations






Pivot Interactives for Chemistry

Pivot Interactives is a powerful supplement to hands-on experimentation, allowing students to vary experimental parameters one at a time to view results from a set of many recordings of the same experiment.

Start a free 30-day trial today at www.pivotinteractives.com



Spectrometer Comparison

Spectrometer	Go Direct SpectroVis Plus	Vernier UV-VIS Spectrophotometer	Vernier Fluorescence/UV-VIS Spectrophotometer
			
Description	The Go Direct SpectroVis Plus Spectrophotometer quickly measures a full-wavelength spectrum. Connect it directly to your device via Bluetooth® wireless technology or via USB.	The Vernier UV-VIS Spectrophotometer generates a full spectrum, Beer's law graph, and kinetics traces of ultraviolet and visible-absorbing samples such as aspirin, DNA, proteins, and NADH.	The Fluorescence/UV-VIS Spectrophotometer measures the fluorescence and absorbance spectra of ultraviolet and visible samples such as quinine sulfate, fluorescein, rhodamine, and DAPI.
Wavelength Range	380 to 950 nm	220 to 850 nm	220 to 850 nm
Light Source	Visible: LED-boosted tungsten Fluorescence: built-in LEDs for excitation at 405 nm and 500 nm	Visible: LED-boosted tungsten UV: Deuterium	Visible: LED-boosted tungsten UV: Deuterium Fluorescence: exchangeable LEDs for excitation at 375 nm, 450 nm, and 525 nm (additional wavelengths sold separately)
Warranty	5 years (1 year on battery, 3 years on lamp, none on consumables)	5 years (1 year on lamp, none on consumables)	5 years (1 year on lamp, none on consumables)
More Information	Innovative use ideas available at www.vernier.com/gdx-svispl	Download free experiments at www.vernier.com/vsp-uv	Download free experiments at www.vernier.com/vsp-fuv
Order Code & Price	GDX-SVISPL \$399	VSP-UV \$2,100	VSP-FUV \$2,899
Optical Fiber Accessory	Vernier Spectrophotometer Optical Fiber This is an optical fiber accessory designed exclusively for emission spectrum experiments with the Vernier-branded spectrophotometers, listed above. It has a wavelength range from 350 to 900 nm. VSP-FIBER \$74		

Learn more at www.vernier.com/spectrometers

Lab Equipment

OHAUS Balances

It is easy to collect mass data from an OHAUS® balance using our popular Logger Pro® 3 software or LabQuest® App. Simply connect a supported balance to the USB port using the OHAUS Scout® USB Cable, start the software, and collect real-time data as if the OHAUS balance were just another Vernier sensor!

OHAUS Scout 120 g

0.001 g precision
OHS-123 ▲ \$660

OHAUS Scout 220 g

0.01 g precision
OHS-222 ▲ \$479

OHAUS Scout 420 g

0.01 g precision
OHS-422 ▲ \$680

All three balances require an OHAUS Scout USB Cable for data collection.

OHAUS Scout USB Cable

OHS-USB \$116

Learn more at www.vernier.com/ohaus



Electrode Support

Our Electrode Support is a great complement to the Vernier Stir Station, as well as a perfect holder for many sensors. It is built to connect to all standard ring stand posts and its large-handled locking nut keeps your sensors firmly in place.

ESUP \$10

Learn more at www.vernier.com/esup




Stir Station

The Stir Station is a high-quality, multi-function magnetic stirrer and ring stand. It includes a Stir Station, Vernier Microstirrer, magnetic stirring bar, AC power adapter, and removable ring stand post. It can be used with AC power (included) or four C batteries (not included).

STIR \$129

Learn more at www.vernier.com/stir



Organic Chemistry

Go Direct Mini GC

Teach students chromatography with an affordable, portable gas chromatograph that detects polar and nonpolar compounds. With the easy-to-use Go Direct® Mini GC™ and the free Vernier Instrumental Analysis™ app, students can separate, analyze, and identify substances contained in a volatile liquid or gaseous sample. Go Direct Mini GC uses Bluetooth® wireless technology or USB to connect to your device.

GDX-GC \$2,499

Learn more at www.vernier.com/gdx-gc

Free Download

Chromatography Experiments with the Go Direct Mini GC e-book

Free with purchase of Go Direct Mini GC




Vernier Instrumental Analysis App

With our free Vernier Instrumental Analysis app, students can collect and analyze data from our Go Direct Mini GC and other advanced instrumentation using computers, Chromebooks, or other mobile devices.

FREE DOWNLOAD

Learn more at www.vernier.com/ia



Organic Chemistry

Polarimeters

Our polarimeters measure chiral properties of optically active samples such as sugars and amino acids. Students no longer have to determine the optical maximum with their eyes but have a graph that shows a clear change in the light's polarization.



NEW Go Direct Polarimeter

GDX-POL \$499



Polarimeter*

CHEM-POL ⚡ \$499

Learn more at www.vernier.com/polarimeters

Melt Stations

Melting point is a physical method of analysis to identify an unknown and purity by its melting temperature. The melt stations accurately measure melting temperatures of a solid (up to 260°C), and the real-time graphing provides a unique perspective of the melting process.



Go Direct Melt Station

GDX-MLT ⚡ \$529



Melt Station*

MLT-BTA ⚡ \$519

Learn more at www.vernier.com/melt-stations

Wide-Range Temperature Probes

The wide-range temperature probes are designed to be used as you would use a thermometer for experiments such as the recrystallization of benzoic acid, simple and fractional distillations, determination of boiling points, the synthesis and analysis of aspirin and other organic compounds, and more.



Go Direct Wide-Range Temperature

GDX-WRT \$114



Wide-Range Temperature Probe*

WRT-BTA \$82

Learn more at www.vernier.com/wr-temp-probes


Organic Chemistry with Vernier

Organic Chemistry with Vernier contains experiments that represent a broad range of topics and techniques taught in most college organic chemistry lab courses. The experiments in this book build upon prior knowledge, laboratory techniques, and skills that students have learned in general chemistry courses.

Topics include

- Distillation
- Chromatography
- Synthesis
- Polarimetry

Learn more at www.vernier.com/chem-o



Download only

CHEM-O-E \$40

Printed book + download


CHEM-O \$48

* requires an interface





⚡ For California Proposition 65 warning, see pp. 142–143.

Featured Products






Go Direct Sensors

Sensor		Order Code	Price
Go Direct® CO ₂ Gas		GDX-CO2	\$199
Go Direct Colorimeter		GDX-COL	\$119
Go Direct Conductivity		GDX-CON	\$99
Go Direct Platinum-Cell Conductivity		GDX-CONPT	\$169
Go Direct Constant Current System		GDX-CCS	\$74
Go Direct Current		GDX-CUR	\$79
Go Direct Drop Counter		GDX-DC	\$99
Go Direct Electrode Amplifier		GDX-EA	\$64
Go Direct Ethanol Vapor		GDX-ETOH	\$149
Go Direct Gas Pressure		GDX-GP	\$89
Go Direct Melt Station		GDX-MLT ▲	\$529
Go Direct ORP		GDX-ORP	\$99

pH Sensors

Go Direct Glass-Body pH		GDX-GPH	\$139
Go Direct pH		GDX-PH	\$89
Go Direct Tris-Compatible Flat pH		GDX-FPH	\$115
Go Direct Radiation Monitor		GDX-RAD	\$179
Go Direct SpectroVis® Plus		GDX-SVISPL	\$399

Temperature Probes

Go Direct Surface Temperature		GDX-ST	\$79
Go Direct Temperature		GDX-TMP	\$69
Go Direct Thermocouple		GDX-TC	\$109
Go Direct Wide-Range Temperature		GDX-WRT	\$114
Go Direct Voltage		GDX-VOLT	\$69

Go Direct Charge Station

Accessory		Order Code	Price
Go Direct Charge Station		GDX-CRG	\$69

See all our products for chemistry at www.vernier.com/chemistry

Looking for Replacement Parts?

Visit www.vernier.com/replacements

LabQuest Sensors

Sensor	Order Code	Price
Colorimeter	COL-BTA	\$119
Conductivity Probes		
Conductivity Probe	CON-BTA	\$99
Platinum-Cell Conductivity Probe	CONPT-BTA	\$149
Current Probes		
Constant Current System	CCS-BTA	\$64
Current Probe	DCP-BTA	\$39
Drop Counter	VDC-BTD	\$99
Electrode Amplifier	EA-BTA	\$49
Gas Pressure Sensors		
Gas Pressure Sensor	GPS-BTA	\$89
Pressure Sensor 400	PS400-BTA	\$189
Instrumentation Amplifier	INA-BTA	\$79
Melt Station	MLT-BTA ⚡	\$519
ORP Sensor	ORP-BTA	\$89
pH Sensors		
Glass-Body pH Electrode BNC (requires Electrode Amplifier)	GPH-BNC	\$85
pH Sensor	PH-BTA	\$88
Tris-Compatible Flat pH Sensor	FPH-BTA	\$104
Polarimeter (Chemical)	CHEM-POL ⚡	\$499
Radiation Monitor	VRM-BTD	\$180

Temperature Probes

Stainless Steel Temperature Probe	TMP-BTA	\$36
Surface Temperature Sensor	STS-BTA	\$25
Thermocouple	TCA-BTA	\$69
Wide-Range Temperature Probe	WRT-BTA	\$82

Voltage Probes

Differential Voltage Probe	DVP-BTA	\$39
Voltage Probe	VP-BTA	\$12

Balances

Sensor	Order Code	Price
OHAUS Scout® (120 g)	OHS-123 ⚡	\$660
OHAUS Scout (220 g)	OHS-222 ⚡	\$479
OHAUS Scout (420 g)	OHS-422 ⚡	\$680

Spectrometers

Spectrometer	Order Code	Price
Go Direct SpectroVis Plus	GDX-SVISPL	\$399
Vernier Emissions Spectrometer	VSP-EM	\$799
Vernier Fluorescence/UV-VIS Spectrophotometer	VSP-FUV	\$2,899
Vernier Spectrometer (Ocean Optics dba Ocean Insight)	V-SPEC	\$1,999
Vernier UV-VIS Spectrophotometer	VSP-UV	\$2,100

Gas Chromatograph

Gas Chromatograph	Order Code	Price
Go Direct Mini GC™	GDX-GC	\$2,499

Lab Equipment and Accessories

Accessory	Order Code	Price
Cuvette Rack	CUV-RACK	\$9
Electrode Support	ESUP	\$10
Melt Station Capillary Tubes	MLT-TUBE	\$19
Plastic Cuvettes (100)	CUV	\$19
Stir Station	STIR	\$129

Lab Books†

Book Title	Order Code	Price
<i>Chemistry with Vernier</i>	CWV	\$48
<i>Advanced Chemistry with Vernier</i>	CHEM-A	\$48
<i>Vernier Chemistry Investigations for Use with AP* Chemistry</i>	APCHEM	\$48
<i>Investigating Chemistry through Inquiry</i>	CHEM-I	\$48
<i>Food Chemistry Experiments</i>	HSB-FOOD	\$38
<i>Organic Chemistry with Vernier</i>	CHEM-O	\$48

† Books listed here include printed book and download; also available as a download only

* AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

See all our products for chemistry at www.vernier.com/chemistry

Physical Science

www.vernier.com/physical-science

From matter and energy to motion and forces, Vernier offers the support you need and the technology your students can use to investigate physical science.



Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training



Physical Science Sets Learning in Motion

Our hands-on physical science investigations help students understand the scientific concepts of real-world phenomena such as energy transfer during phase changes, the cooling effect of evaporation, and principles of simple machines.

Physical Science


Physical Science with Vernier

Physical Science with Vernier contains 40 ready-to-use experiments for physical science. Experiments are included for nine Vernier sensors and cover a variety of topics in chemistry and physics.

Topics include

- Structures and properties of matter
- Forces and interactions
- Waves and electromagnetic radiation
- Chemical reactions

Learn more at www.vernier.com/psv



Download only

PSV-E \$40

Printed book + download

PSV \$48

Go Direct Sensor Carts


With our Go Direct® Sensor Carts, students can explore force, position, velocity, and acceleration directly on their devices via Bluetooth® wireless technology—no wires or additional equipment required. Each cart features built-in sensors to simplify experiment setup.

Go Direct Sensor Cart (Green)

GDX-CART-G \$169

Go Direct Sensor Cart (Yellow)

GDX-CART-Y \$169



www.vernier.com/gdx-cart

EXPERIMENT 23

Reflectivity of Light

After comparing the amount of light reflected from different colors of paper, students apply the results to help answer their questions about planetary albedo.



Sensor Used



Go Direct Light and Color

Students use this sensor to measure the brightness of a light bulb or the reflectance of light off of various objects. They can also measure UV light and relative amounts of red, blue, and green light.

GDX-LC \$79

Experiment Source



Physical Science with Vernier

Download only: PSV-E \$40

Printed book + download: PSV \$48

Learn more at www.vernier.com/psv-23

EXPERIMENT 3

Freezing and Melting of Water

Students measure the temperature of water as it changes from a liquid to a solid. The data are analyzed to make predictions about the freezing patterns of other substances.



Sensor Used



Go Direct® Temperature

This is a rugged, general-purpose sensor that students can use to monitor temperature.

GDX-TMP \$69

Experiment Source



Physical Science with Vernier

Download only: PSV-E \$40

Printed book + download: PSV \$48

Learn more at www.vernier.com/psv-3

EXPERIMENT 21

Pulleys

By comparing the effort force to the resistance force required to lift a mass, students determine the mechanical advantage of different pulley systems.



Sensor Used



Go Direct Force and Acceleration

Students can use this sensor to measure forces of up to 50 N. The included 3-axis accelerometer makes it a versatile sensor for many topics in physical science.

GDX-FOR \$99

Experiment Source



Physical Science with Vernier

Download only: PSV-E \$40

Printed book + download: PSV \$48

Learn more at www.vernier.com/psv-21

Featured Products

Go Direct Sensors

Sensor	Order Code	Price
Go Direct 3-Axis Magnetic Field	GDX-3MG	\$69
Go Direct Acceleration	GDX-ACC	\$99
Carts and Tracks		
Dynamics Cart and Track System with Go Direct Sensor Carts	DTS-GDX	\$535
Go Direct Sensor Cart (Green)	GDX-CART-G	\$169
Go Direct Sensor Cart (Yellow)	GDX-CART-Y	\$169
Go Direct Conductivity	GDX-CON	\$99
Go Direct Current	GDX-CUR	\$79
Go Direct Energy	GDX-NRG	\$89
Go Direct Force and Acceleration	GDX-FOR	\$99
Go Direct Gas Pressure	GDX-GP	\$89
Go Direct Light and Color	GDX-LC	\$79
Go Direct Motion	GDX-MD	\$99
Go Direct pH	GDX-PH	\$89
Go Direct Photogate	GDX-VPG	\$89
Go Direct Sound	GDX-SND	\$89
Go Direct Structures & Materials Tester	GDX-VSMT	\$999
Temperature Probes		
Go Direct Surface Temperature	GDX-ST	\$79
Go Direct Temperature	GDX-TMP	\$69
Go Direct Thermocouple	GDX-TC	\$109
Go Direct Voltage	GDX-VOLT	\$69

Go Direct Charge Station

Accessory	Order Code	Price
Go Direct Charge Station	GDX-CRG	\$69

LabQuest Sensors

Sensor	Order Code	Price
Accelerometers		
3-Axis Accelerometer	3D-BTA	\$99
25-g Accelerometer	ACC-BTA	\$96
Low-g Accelerometer	LGA-BTA	\$89
Conductivity Probe	CON-BTA	\$99
Current Probes		
Current Probe	DCP-BTA	\$39
High Current Sensor	HCS-BTA	\$79
Energy Sensor	VES-BTA	\$88
Force Sensors		
Dual-Range Force Sensor	DFS-BTA	\$109
Force Plate	FP-BTA	\$289
Gas Pressure Sensor	GPS-BTA	\$89
Light Sensor	LS-BTA ⚡	\$59
Magnetic Field Sensor	MG-BTA	\$58
Microphone	MCA-BTA	\$44
Motion Detector	MD-BTD	\$89
pH Sensor	PH-BTA	\$88
Photogate	VPG-BTD	\$49
Sound Level Sensor	SLS-BTA	\$69
Temperature Probes		
Go!Temp® (USB Sensor)	GO-TEMP	\$39
Stainless Steel Temperature Probe	TMP-BTA	\$36
Surface Temperature Sensor	STS-BTA	\$25
Thermocouple	TCA-BTA	\$69

Looking for Replacement Parts?

Visit www.vernier.com/replacements

Voltage Probes		
30-Volt Voltage Probe	30V-BTA	\$49
Differential Voltage Probe	DVP-BTA	\$39
Voltage Probe	VP-BTA	\$12

Accessories & Lab Equipment

Product	Order Code	Price
Balances		
OHAUS Scout® (120 g)	OHS-123 ⚡	\$ 660
OHAUS Scout (220 g)	OHS-222 ⚡	\$ 479
OHAUS Scout (420 g)	OHS-422 ⚡	\$ 680
Electrode Support	ESUP	\$10
pH Storage Solution	PH-SS	\$20
pH Buffer Capsules Kit	PH-BUFCAP	\$29
Stir Station	STIR	\$129
Vernier Circuit Board 2	VCB2 ⚡	\$129

Lab Books

Title	Order Code	Price
<i>Physical Science with Vernier</i>	Printed book + download: PSV	\$48
	Download only: PSV-E	\$40
<i>Chemistry with Vernier</i>	Printed book + download: CWV	\$48
	Download only: CWV-E	\$40
<i>Physics with Vernier</i>	Printed book + download: PWV	\$48
	Download only: PWV-E	\$40

See all our products for physical science at www.vernier.com/physical-science

Physics

www.vernier.com/physics

From kinematics to optics, Vernier technology helps your students connect the dots between the classroom and the real world. Our physics products enable student and educator success so that you can spend less time troubleshooting and more time teaching your students about the scientific principles of the world around them.



Topics

Explore a sampling of our featured experiments by topic to learn how Vernier technology helps your students engage with data-collection technology and deepens their understanding of key physics concepts.

1-D Motion and Force

PAGE 98

2-D Motion and Force

PAGE 106

Electricity and Magnetism

PAGE 108

Thermodynamics

PAGE 110

Waves and Sound



PAGE 112

Light and Optics

PAGE 113

Modern Physics

PAGE 116



A Guide to Vernier Data Collection

GDX

Our Go Direct® technology connects directly to compatible student devices—computers, Chromebooks, LabQuest® 3, and iOS, iPadOS®, and Android™ devices. Its ease of use maximizes valuable lab time so you can focus on teaching.


LQ

With over 80 sensors to choose from, our LabQuest family of sensors offers a wide variety of experiments to integrate into your existing curriculum. Connect LabQuest sensors with an interface to your device, or use LabQuest 3 as a standalone device in the field or lab.

Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training



1-D Motion and Force

Featured Experiments

EXPERIMENT 1

Graph Matching

Kinesthetic experience coupled with real-time graphing helps cement student understanding of the relationships between motion, position vs. time graphs, and velocity vs. time graphs.



Sensor Used

GDX



Go Direct® Motion

Go Direct Motion uses ultrasound to measure the position, velocity, and acceleration of moving objects.

GDX-MD \$99

Can also be done with

LQ

Motion Detector

MD-BTD \$89

Go! Motion® (USB motion detector)

GO-MOT \$129

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-1

EXPERIMENT 12

Static and Kinetic Friction

Make investigating friction easy with a digital force sensor. Students re-create the friction graph from their textbook while determining coefficients of static and kinetic friction.



Sensor Used

GDX



Go Direct Force and Acceleration

Measure forces as small as ± 0.1 N and up to ± 50 N with this sensor that couples a 3-axis accelerometer with a stable and accurate force sensor. Use it to measure pushes and pulls in the classroom or outdoors.

GDX-FOR \$99

Can also be done with

LQ

Dual-Range Force Sensor

DFS-BTA \$109

GDX

Go Direct Sensor Cart (green or yellow)

GDX-CART-G \$169 (green)

GDX-CART-Y \$169 (yellow)

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48


Learn more at www.vernier.com/pwv-12

- GDX** connects directly to devices
- LQ** requires an interface

EXPERIMENT 21

Accelerations in the Real World

In this inquiry activity, students take an acceleration sensor out of the classroom and into different situations, whether it be cars, elevators, amusement parks, or elsewhere.



Sensor Used

GDX



Go Direct Acceleration

Collect acceleration, rotation, and altitude data in the classroom or in the field.

GDX-ACC \$99

Can also be done with

LQ

3-Axis Accelerometer

3D-BTA \$99

GDX

Go Direct Force and Acceleration

GDX-FOR \$99

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-21

EXPERIMENT 14

Pendulum Periods

Take a classic experiment to the next level with precision measurement of pendulum period. Students test three variables to discover which factors influence the period.



Sensor Used

GDX



Go Direct Photogate

This double-gate sensor includes two photogates built into the arms of the sensor. It accurately measures velocity and acceleration.

GDX-VPG \$89

Can also be done with

LQ

Vernier Photogate

VPG-BTD \$49

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-14

Dynamics Cart and Track Systems

One Dynamics System—Three Ways to Collect Data

Depending on your budget and your needs, we offer three ways to collect motion data.

1

Go Direct Sensor Cart GDX

The wireless Go Direct® Sensor Cart includes an optical encoder on a wheel to sense the displacement of the cart, on or off the track. No interface is needed to use this system with our free Vernier Graphical Analysis™ app. Students can perform impulse and momentum experiments with the built-in force sensor, and the 3-axis accelerometer means you can take your Sensor Cart off campus to investigate accelerations on a swing or merry-go-round.




2

The Motion Encoder* LQ

VERNIER EXCLUSIVE

For classrooms already equipped with data-collection interfaces, the Motion Encoder dramatically improves data quality and simplifies experiment setup over the traditional ultrasonic Motion Detector. An optical sensor under the dynamics cart senses the passage of the cart over a striped decal on the track. The displacement information is sent as an encoded IR signal to a receiver at the track's end. This optical-only system provides excellent, repeatable, and noise-resistant data.

* U.S. Patent No. 9,488,503




3

A Traditional Motion Detector GDX LQ

The Motion Detector is the classic method for collecting position data. Use a Motion Detector bracket to measure cart motion for the entire length of the track. You can even use two Motion Detectors at once to study cart collisions.

Unlike the Motion Encoder or Go Direct Sensor Cart, the Motion Detector can be used for dynamics experiments other than cart-on-track experiments. Students can graph their own walking motion, study a simple pendulum, or graph a ball toss with a Motion Detector. If you want to use a Motion Detector for all motion experiments, get the Dynamics Cart and Track System without the Motion Encoder or Go Direct Sensor Cart.



Dynamics Cart and Track System with Go Direct Sensor Cart

BUILT-IN SENSORS = LOWER TOTAL COST

The Dynamics Cart and Track System with Go Direct Sensor Cart includes essential laboratory equipment for teaching dynamics and kinematics. With our Go Direct Sensor Cart, students can explore force, position, velocity, and acceleration directly on their device using Bluetooth® wireless technology. There are no wires to create drag, and no additional equipment is required! Each cart features built-in sensors that simplify experiment setup and make this system the best choice for studying dynamics and kinematics.

with 1.2 m Track DTS-GDX \$535 www.vernier.com/dts-gdx

with 2.2 m Track DTS-GDX-LONG \$639 www.vernier.com/dts-gdx-long



Dynamics Cart and Track System with Motion Encoder


RECOMMENDED OPTION FOR USE WITH LOGGER PRO® 3

The Dynamics Cart and Track System with Motion Encoder includes an optical position sensing system to record cart motion.


with 1.2 m Track DTS-EC \$445 www.vernier.com/dts-ec

with 2.2 m Track DTS-EC-LONG \$549 www.vernier.com/dts-ec-long

Motion encoder data are so pristine that you can usefully graph jerk vs. time.



LQ




Dynamics Cart and Track System

USE WITH SENSORS YOU ALREADY OWN—SENSORS ARE NOT INCLUDED.

The Dynamics Cart and Track System features the Combination Track/Optics Bench, two low-friction plastic carts (one standard and one with an adjustable plunger), and attachment accessories.

with 1.2 m Track DTS \$295 www.vernier.com/dts

with 2.2 m Track DTS-LONG \$399 www.vernier.com/dts-long



Dynamics Cart and Track Systems

EXPERIMENT 4

Determining g on an Incline

Students mimic Galileo's seminal experiment with modern tools using a low-friction setup to determine the acceleration of gravity on Earth.



Watch
a video



Sensor Used



Dynamics Cart and Track System with Go Direct® Sensor Cart

This completely wireless system simplifies experiment setup and allows basic experiments to be conducted with or without the track.

DTS-GDX \$535

Can also be done with

LQ

Dynamics Cart and Track System
with Motion Encoder

DTS-EC \$445

LQ

Motion Detector and Dynamics
Cart and Track System

MD-BTD \$89

DTS \$295

GDX

Go Direct Motion and Dynamics
Cart and Track System

GDX-MD \$99

DTS \$295



Go Direct Sensor Carts

We've added wireless sensors to our popular dynamics cart. Each cart includes an encoder wheel to report position, velocity, and acceleration. Conduct basic physics investigations with or without a track.

Go Direct Sensor Cart
(Green)

GDX-CART-G \$169

Go Direct Sensor Cart
(Yellow)


GDX-CART-Y \$169



INCLUDES
21
EXPERIMENTS

NEW

Sensor Cart Physics



Download only
HSB-SCP-E \$30

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-4a

www.vernier.com/gdx-cart

Dynamics Cart and Track Systems—Featured Kits and Accessories

Fan Cart

The Fan Cart works with a motion detector and the Vernier Dynamics Cart and Track System. Study Newton's second law using variable fan thrust and included mass bars.

CART-F \$109

www.vernier.com/cart-f



Encoder Fan Cart

Use the Encoder Fan Cart with the Motion Encoder System. Study Newton's second law using variable fan thrust and included mass bars.

CART-FEC \$225

www.vernier.com/cart-fec

LQ



Friction Pad DTS

Add a Friction Pad to any of our plastic dynamics carts to study the effect of consistent friction on the motion of the cart.

DTS-PAD \$32

www.vernier.com/dts-pad



Motion Encoder Cart and Receiver

This kit includes a fully assembled Motion Encoder Cart, as well as the Motion Encoder Receiver and Motion Encoder Long Track Strip.

DTS-MEC \$244

www.vernier.com/dts-mec

LQ



Eddy Current Brake

Eddy current brakes are used as a braking system for high-speed trains and roller coasters. Recreate this unusual braking system in your classroom or laboratory by installing our Eddy Current Brake into the end cap of a plastic Vernier dynamics cart. As the cart moves over the track, the magnets in the Eddy Current Brake create an electromagnetic drag on the cart that is proportional to the cart's speed.

DTS-ECB \$19

www.vernier.com/dts-ecb




Bumper and Launcher Kit

With the Bumper and Launcher Kit, students can use the Dynamics Cart and Track System to perform Hooke's law experiments or to study momentum and impulse. The kit includes

- Clay (~20 grams)
- Clay holders (2)
- Dual-magnet bumper
- Force sensor mounting screw
- Hoop bumpers (2)
- Magnetic bumpers (2)
- Rubber bumpers (2)
- Track bracket

BLK \$89

www.vernier.com/blk



Track and Force Sensor not included

Featured Products

Motion Detectors

Go Direct Motion

GDX

Go Direct® Motion uses ultrasound to measure the position, velocity, and acceleration of moving objects. It connects via Bluetooth® wireless technology or via USB to your device.

GDX-MD \$99



Motion Detector

LQ

The Motion Detector uses ultrasound to measure the position of carts, balls, people, and other objects. It can be used with interfaces from the LabQuest® family, LabPro®, and CBL 2.™ It is not supported with Go!Link® or EasyLink®.

MD-BTD \$89



Go! Motion

Go!Motion® is our motion detector that connects directly to a computer or Chromebook™ USB port—eliminating the need for an additional data-collection interface. This USB motion detector works with Logger Pro® 3, Vernier Graphical Analysis™ app, and Graphical Analysis Pro.

GO-MOT \$129



www.vernier.com/motion-detectors

Photogates

Go Direct Photogate

GDX

Go Direct Photogate is a double-gate sensor that includes two photogates built into the arms of the sensor, which accurately measures velocity and acceleration without needing to know anything about the geometry of the object. Go Direct Photogate also includes a single laser gate for use with objects passing outside of the arms of the sensor (required visible light laser not included). The sensor can be used to study free fall, rolling objects, collisions, and pendulums.

GDX-VPG \$89



Photogate

LQ

Study free fall, rolling objects, collisions, and pendulums with the Vernier Photogate. Use the built-in laser detector to create a photogate through which you could drive a truck. It includes an accessory rod for attaching to a ring stand or for adding the Ultra Pulley Attachment (sold separately).

VPG-BTD \$49



Picket Fence

GDX

LQ

PF \$9



Ultra Pulley Attachment

GDX

LQ

SPA \$24



www.vernier.com/photogates

Featured Products

Accelerometers

Go Direct Acceleration

GDX

Collect acceleration, rotation, and altitude data in the classroom or in the field. This 3-axis acceleration sensor has two acceleration ranges plus an altimeter and a 3-axis gyroscope.

Acceleration ranges: $\pm 157 \text{ m/s}^2$, $\pm 1960 \text{ m/s}^2$

Gyroscope: 3 axis, $\pm 35 \text{ rad/s}$

Altimeter: $-1,800$ to $10,000 \text{ m}$

GDX-ACC \$99



Low-g Accelerometer

LQ

Use the Low-g Accelerometer to study the one-dimensional motion of a car (real or toy), pendulum bob, an elevator, or an amusement park ride.

Range: $\pm 50 \text{ m/s}^2$

LGA-BTA \$89



3-Axis Accelerometer

LQ

Range: $\pm 50 \text{ m/s}^2$

3D-BTA \$99



25-g Accelerometer

LQ

Range: $\pm 250 \text{ m/s}^2$

ACC-BTA \$96



www.vernier.com/accelerometers

Force Sensors

Go Direct Force and Acceleration

GDX

Go Direct Force and Acceleration includes a $\pm 50 \text{ N}$ force sensor, a 3-axis accelerometer, and a 3-axis gyroscope. Take it on an amusement park ride, mount it on a dynamics cart, or attach a string and whirl it in a horizontal or vertical circle—in wireless mode, your imagination is the only limiting factor!

Force: $\pm 50 \text{ N}$

Acceleration: 3 axis, $\pm 16 \text{ g}$

Gyroscope: 3 axis, $\pm 35 \text{ rad/s}$

GDX-FOR \$99



Dual-Range Force Sensor

LQ

Using our Dual-Range Force Sensor, students can test Newton's third law of motion, explore Hooke's law, or graph the transition from static friction to kinetic friction.

Ranges: $\pm 10 \text{ N}$, $\pm 50 \text{ N}$

DFS-BTA \$109



Force Plate


LQ

The Force Plate—a force sensor about the size of a bathroom scale—is tough enough to jump on. Two handles are included for pushing or pulling.

Ranges: -850 to $+3500 \text{ N}$

-200 to $+850 \text{ N}$

FP-BTA \$289



www.vernier.com/force-sensors

2-D Motion and Force

Featured Experiments

EXPERIMENT 8B

Projectile Motion

Predict the landing point of a projectile based on the launch velocity and initial height. With precision photogate timing, success depends on student understanding.



Sensor Used



Go Direct Projectile Launcher

Use the Go Direct® Projectile Launcher to investigate important concepts in two-dimensional kinematics. Launch steel balls at angles between 0 and 90 degrees and over distances up to 2.5 m.

GDX-PL \$449

Can also be done with

LQ

Vernier Projectile Launcher

VPL \$389

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-8b

EXPERIMENT 12A

Centripetal Acceleration

Students explore the relationships among force, speed, and radius through reliable data collection using sensors.



Sensors Used



Go Direct Centripetal Force Apparatus

This is an ideal combination to explore rotational dynamics when combined with Go Direct Force and Acceleration (not included).

GDX-CFA \$299



Go Direct Force and Acceleration

This couples a 3-axis accelerometer with a stable and accurate force sensor that measures forces as small as ± 0.1 N and up to ± 50 N. Measure angular rotation using the 3-axis gyroscope.

GDX-FOR \$99

Can also be done with

LQ

Centripetal Force Apparatus

CFA \$549

LQ

Dual-Range Force Sensor

DFS-BTA \$109

LQ

Photogate

VPG-BTD \$49

Experiment Source



Advanced Physics with Vernier—Mechanics

Download only: PHYS-AM-E \$40

Printed book + download: PHYS-AM \$48

Learn more at www.vernier.com/phys-am-12a

GDX

connects directly to devices


LQ

requires an interface

EXPERIMENT 13

Rotational Dynamics

Apply a torque and measure an angular acceleration. Students explore the version of Newton's second law that applies to rotation.



Sensor Used

GDX



Go Direct Rotary Motion

Measure angular displacement, angular velocity, and angular acceleration easily and precisely.

GDX-RMS ▲ \$179

Can also be done with

LQ

Rotary Motion Sensor

RMV-BTD ▲ \$169

Accessories Used



Rotational Motion Accessory Kit

Used with a rotary motion sensor to study the motion of a physical pendulum; the rotational inertia of disks, rings, and point masses; and the conservation of angular momentum

AK-RMV ▲ \$112

Experiment Source



Advanced Physics with Vernier—Mechanics

Download only: PHYS-AM-E \$40

Printed book + download: PHYS-AM \$48

Learn more at www.vernier.com/phys-am-13

Featured Products

Go Direct Acceleration

GDX

Collect acceleration, rotation, and altitude data in the classroom or in the field. This 3-axis acceleration sensor has two acceleration ranges plus an altimeter and a 3-axis gyroscope.

GDX-ACC \$99

www.vernier.com/gdx-acc

Projectile Launcher Accessories

GDX

LQ

Independence of Motion Accessory

The Independence of Motion Accessory enables students to use the Vernier Projectile Launcher to perform the classic experiment where one ball is dropped as another is projected horizontally. The balls strike the floor simultaneously.

IOM-VPL \$59

www.vernier.com/iom-vpl

Time of Flight Pad

The Time of Flight Pad is used with a projectile launcher or photogate (not included) to precisely measure how long a projectile has been in motion.

TOF-VPL \$84

www.vernier.com/tof-vpl

Centripetal Force Apparatus Accessories

Moment of Inertia Kit

CFA-MIK \$179

www.vernier.com/cfa-mik

GDX LQ



Motor Accessory Kit

GDX-CFA-MAK \$189

www.vernier.com/gdx-cfa-mak

GDX



Sensor Bracket

CFA-SBK \$19

www.vernier.com/cfa-sbk

GDX LQ




Electricity and Magnetism

Featured Experiments

EXPERIMENT 6

Electrostatics

Using our Go Direct Static Charge (essentially a digital electroscope), students explore charging by friction, conduction, and induction.



Sensor Used



Go Direct® Static Charge

With Go Direct Static Charge, students can easily measure and analyze static charges. Designed with affordability and ease of use in mind, this sensor ensures enhanced performance so that students can collect accurate data.

GDQ-Q \$99

Accessory Used



Electrostatics Kit

Students use the Electrostatics Kit to perform a range of experiments in electrostatics with the Go Direct Static Charge.

ESK-CRG \$119

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40


Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-6

EXPERIMENT 22

Ohm's Law

Students compare the potential vs. current graphs for resistors and for a light bulb in this exploration of Ohm's law.



Sensors Used



Go Direct Voltage

This sensor combines a wide input voltage range and high precision, making it an excellent choice for investigations of both AC/DC circuits and electromagnetism.

GDQ-VOLT \$69



Go Direct Current

Measure electric currents in circuits with this versatile sensor.

GDQ-CUR \$79

Accessory Used

Vernier Circuit Board 2

VCB2 \$129

Can also be done with

LQ

Differential Voltage Probe

DVP-BTA \$39

LQ

Current Probe

DCP-BTA \$39



Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-22


GDX connects directly to devices

LQ requires an interface

EXPERIMENT 25

Magnetic Field of a Coil

How do different factors affect the magnetic field in the center of a coil of wire? Students investigate the number of turns and the amount of current in a wire coil.



Sensor Used

GDX



Go Direct 3-Axis Magnetic Field

Determine the magnitude and direction of a magnetic field at any point in space with this 3-axis sensor.

GDX-3MG \$69

Can also be
done with

LQ

**Magnetic
Field
Sensor**

MG-BTA
\$58

Accessory Used



Extech® Digital Power Supply

This power supply provides constant current or constant voltage for physics activities that require DC power.

EXPS \$240

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-25

Featured Products

Additional LabQuest Voltage and Current Probes

LQ

Sensor	Range	URL
Current Probe	±0.6 A	www.vernier.com/dcp-bta
High Current Sensor	±10 A	www.vernier.com/hcs-bta
Instrumentation Amplifier	±1 V	www.vernier.com/ina-bta
Differential Voltage Probe	±6 V	www.vernier.com/dvp-bta
Voltage Probe	±10 V	www.vernier.com/vp-bta
30-Volt Voltage Probe	±30 V	www.vernier.com/30v-bta

Power Amplifier

LQ



Use this as a power supply for DC and AC circuit investigations or to drive devices such as speakers, lamps, and small DC motors.

PAMP \$225

High-Voltage Electrostatics Kit

LQ



Investigate the distribution of charge on a sphere, transfer of charge on contact between two spheres, and charging by induction with this kit.

HVEK-CRG \$289

Electrostatic High-Voltage Genecon

LQ



A great addition to the High Voltage Electrostatics Kit, the Electrostatic High-Voltage Genecon generates both positive and negative charges and reliably creates charge differences in high humidity.

HVEK-GEN \$229

Vernier Circuit Board 2



Use this convenient platform to study basic series and parallel circuits as well as RLC circuits. Many components for experimentation are provided.

VCB2 \$129

Optional Breadboard Kit

for the Vernier Circuit Board 2



Install this small breadboard to easily conduct experiments using additional electronic components not permanently mounted on the Vernier Circuit Board 2.

VCB2-OBK \$29

EXPERIMENT 1

Behavior of a Gas

Students collect pressure and temperature data to discover kinetic molecular theory and the iconic expression $PV = nRT$.



Sensors Used

GDx



Go Direct® Gas Pressure

Measure the absolute pressure of a gas.

GDx-GP \$89

GDx



Go Direct Temperature

This is a rugged, general-purpose sensor that students can use to monitor temperature.

GDx-TMP \$69

Can also be done with

LQ

Gas Pressure Sensor

GPS-BTA \$89

LQ

Stainless Steel Temperature Probe

TMP-BTA \$36

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40


Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-1

INNOVATIVE USE

Radiant Energy with FLIR ONE®

Visible light interacts with matter in different ways, depending on the color of the matter. Students use a thermal camera to measure the invisible infrared light that results.



Sensor Used



FLIR ONE Pro Thermal Camera for iOS

Reveal the hidden world of infrared vision. When used with our Vernier Thermal Analysis® Plus app, students can also collect temperature vs. time data for up to four spots or regions, along with a thermal image video.

FLIRPRO-IO5 ⬆ \$399

Can also be done with

FLIR ONE Pro LT

FLIRLT-IO5 ⬆ \$299

FLIR ONE

FLIRONE3-IO5 ⬆ \$199

Software Used



Vernier Thermal Analysis Plus App for FLIR ONE

Students can easily observe temperature changes on the skin, illustrate convection, detect heating due to friction, compare heat conduction in different materials, and analyze the transparency of materials in infrared light.



Experiment Source



FREE DOWNLOAD

www.vernier.com/radiant-energy

Featured Products

FLIR ONE Thermal Cameras

Using a FLIR ONE Thermal Camera, students can observe temperature changes on the skin, illustrate convection, track heating due to friction, compare heat conduction in different materials, analyze the transparency of materials in infrared compared to visible light, and so much more.

FLIR ONE Pro

FLIRPRO-IOS 🚩 \$399



FLIR ONE Pro LT

FLIRLT-IOS 🚩 \$299



FLIR ONE Gen 3

FLIRONE3-IOS 🚩 \$199



www.vernier.com/flir

Vernier Thermal Analysis Plus App

The Vernier Thermal Analysis® Plus app makes it possible to analyze temperatures of up to four spots or regions and collect temperature data as a function of time. Examine the in-app graph, select different points or regions to examine, collect time-lapse videos for longer experiments, or export data to the Logger Pro® 3 or Graphical Analysis™ GW app for further analysis.

www.vernier.com/thermal-analysis



Gas Pressure Sensors

Go Direct Gas Pressure

Range: 0 to 400 kPa

GDX-GP \$89



Gas Pressure Sensor

Range: 0 to 210 kPa

GPS-BTA \$89



www.vernier.com/gas-pressure-sensors

Temperature Probes

Go Direct Surface Temperature

Range: -25 to 125°C

GDX-ST \$79



Surface Temperature Sensor

Range: -25 to 125°C

STS-BTA \$25



www.vernier.com/temperature-sensors

Go Direct Temperature

Range: -40 to 125°C

GDX-TMP \$69



Stainless Steel Temperature Probe

Range: -40 to 135°C

TMP-BTA \$36




Waves and Sound

Featured Experiments

EXPERIMENT 32

Sound Waves and Beats



Compare data from sound waves with sinusoidal functions. What information is contained in each parameter? Students also observe sound wave interference.

Sensor Used

CDX



Go Direct® Sound

Use this sensor to easily capture and evaluate waveforms.

GDY-SND \$89

Can also be done with

LQ Microphone
MCA-BTA \$44

Experiment Source




Physics with Vernier

Download only:
PWV-E \$40
Printed book +
download: PWV \$48

Learn more at www.vernier.com/pwv-32

EXPERIMENT 3

Standing Waves on a String



Students explore waves on a string that is fixed at both ends, create harmonics, and relate string tension and wave speed.

Products Used



Power Amplifier

Drive devices such as speakers, lamps, and small DC motors.

PAMP \$225



Power Amplifier Accessory Speaker

Study mechanical waves on strings and springs.

PAAS-PAMP \$125

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only:
PHYS-ABM-E \$40
Printed book +
download:
PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-3

Featured Products

Microphone



Display and study the waveforms of sounds from voices and musical instruments. This sensor is also appropriate for speed of sound experiments.

MCA-BTA \$44 www.vernier.com/mca-bta

Sound Level Sensor



Use the Sound Level Sensor to easily measure sound level in decibels (dB) in a variety of experiments.

Range: 55 to 110 dB

SLS-BTA \$69 www.vernier.com/sls-bta



EXPERIMENT 29

Light, Brightness, and Distance

Illuminate the inverse square law for light intensity in this experiment, which requires a dark room and a point source of light in addition to a light sensor.



Sensor Used

GDx



Go Direct Light and Color

Measure light intensity in the visible to ultraviolet electromagnetic spectrum. An RGB color sensor detects relative contributions of primary colors in light.

GDx-LC \$79

Can also be done with

LQ

Light Sensor

LS-BTA [▲]
\$59

Accessories Used



Optics Expansion Kit

OEK \$179



Combination 1.2 m Track/Optics Bench

TRACK \$135

Experiment Source



Physics with Vernier

Download only: PWV-E \$40

Printed book + download: PWV \$48

Learn more at www.vernier.com/pwv-29

EXPERIMENT 16

Thin Lenses and Real Images

The number 4 has no symmetry, making it an ideal shape for examining real, inverted images. Students measure object and image distances and sizes to determine focal length and magnification.



Accessories Used



Optics Expansion Kit

Add this kit to your Dynamics Cart and Track System to conduct optics experiments, such as image formation with lenses and light intensity vs. distance. You can even use the kit to build a basic telescope.

OEK \$179



Combination 1.2 m Track/Optics Bench

TRACK \$135

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40

Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-16

EXPERIMENT 15

Curved Mirrors and Images

Students focus real images on a half screen and use parallax to locate a virtual image in this standard optics experiment.



Accessories Used



Optics Expansion Kit

Add this kit to your Dynamics Cart and Track System to conduct optics experiments, such as image formation with lenses and light intensity vs. distance. You can even use the kit to build a basic telescope.

OEK \$179



Mirror Set for Optics Expansion Kit

This set extends the kit so students can easily study image formation by concave and convex mirrors.

M-OEK \$59



Combination 1.2 m Track/Optics Bench

TRACK \$135

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40


Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-15

EXPERIMENT 19

Interference

Explore the wave nature of light with the classic double-slit experiment for light. Students can vary slit width and separation. In addition, they can study single-slit diffraction.



Accessories Used



Diffraction Apparatus

This set extends the kit so students can easily study image formation by concave and convex mirrors.

DAK \$620

Combination 1.2 m Track/Optics Bench

TRACK \$135

Green Diffraction Laser (optional)

Add this to your Diffraction Apparatus to study the effect of wavelength on a diffraction pattern.

GDL-DAK \$250

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40

Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-19

Featured Products

Light Sensors

Go Direct® Light and Color

GDX

This sensor combines the power of visible light, UV, and RGB sensors to measure source emission, transmittance, and reflection of light in the visible light to ultraviolet electromagnetic spectrum.

GDX-LC \$79



Light Sensor

LQ

Investigate polarizers, reflectivity, and solar energy with this sensor that approximates the human eye in spectral response. It's great for inverse square law experiments.

LS-BTA ⚡ \$59



www.vernier.com/light-sensors

Polarizer/Analyzer Set

Using the Polarizer/Analyzer Set, students can study light polarization and do experiments such as Malus's law. The set consists of three adjustable linear polarizers, one of which includes attachment points for either of our rotary motion sensors. It requires components from the Optics Expansion Kit and either a LabQuest® Light Sensor or Go Direct® Light and Color for use.

PAK-OEK \$85

www.vernier.com/pak-oek



Optics Expansion Kit

Use the Optics Expansion Kit with your dynamics track (not included) to conduct optics experiments, such as image formation with lenses and light intensity vs. distance. You can even use the kit to build a basic telescope.

Kit includes

- 3 lenses (100 mm converging lens, 200 mm converging lens, -150 mm diverging lens)
- Screen
- Combination luminous and point light source
- Light Sensor Holder*
- Aperture screen
- Power supply

The Optics Expansion Kit is used in *Physics with Vernier* and *Advanced Physics with Vernier—Beyond Mechanics* experiments.

OEK \$179

Download free sample experiments at www.vernier.com/oek

See website for replacement parts.

* The Light Sensor Holder can be used with any style Vernier light sensor.

Mirror Set

The Mirror Set extends the Optics Expansion Kit so students can easily study image formation by concave and convex mirrors. The set includes a concave mirror, a convex mirror, and a half screen. It requires components from the Optics Expansion Kit for use.

M-OEK \$59

www.vernier.com/m-oek




Light source not included

Combination Dynamics Track and Optical Bench

The Combination Dynamics Track and Optical Bench is aluminum and includes a metric scale. Extremely rigid, this 1.2 (or 2.2) meter track will not sag under use. The track includes two Adjustable Two Foot Levelers.

with 1.2 m Track TRACK \$135 www.vernier.com/track

with 2.2 m Track TRACK-LONG \$239
www.vernier.com/track-long



Color Mixer

The Color Mixer accessory can be used to study the mixing of red, blue, and green light by additive and subtractive mixing. It requires a Combination Track/Optics Bench (not included).

CM-OEK \$175

Download a free sample experiment at www.vernier.com/cm-oek



See website for replacement parts.

EXPERIMENT 21

The Spectrum of Atomic Hydrogen

Compare the spectrum of an incandescent lamp with the few lines of the hydrogen spectrum.



Sensor Used



Vernier Emissions Spectrometer

This emissions spectrometer gives precise measurements over a range of 350–900 nm. Use it to examine spectra of light bulbs, spectrum tubes, or the sun.

VSP-EM \$799

Accessories Used



Spectrum Tube Single Power Supply

These power supplies feature an ultra-safe design for electrifying spectrum tubes.

ST-SPS \$255



Spectrum Tube (Hydrogen)

ST-H \$45



Vernier Emissions Fiber

VSP-EM-FIBER \$88

Experiment Source



Advanced Physics with Vernier—Beyond Mechanics

Download only: PHYS-ABM-E \$40


Printed book + download: PHYS-ABM \$48

Learn more at www.vernier.com/phys-abm-21

EXPERIMENT 2

Distance and Radiation

Students use a gamma emitter and radiation monitor to determine the relationship between radiation counts and distance. This is a great follow-up to our Light, Brightness, and Distance experiment (see page 113)!



Sensor Used

GDX



Go Direct® Radiation Monitor

Use this sensor to detect alpha, beta, gamma, and X-ray radiation.

GDX-RAD \$179

Can also be done with

LQ

Vernier Radiation Monitor

VRM-BTD \$180

Experiment Source



Nuclear Radiation with Vernier

FREE DOWNLOAD www.vernier.com/nrv

Featured Products

Vernier Emissions Spectrometer

The Vernier Emissions Spectrometer gives precise measurements over a range of 350–900 nm. Use it with or without an optical fiber (not included) to examine spectra of light bulbs, spectrum tubes, or the sun.

VSP-EM \$799

www.vernier.com/vsp-em



Vernier Emissions Fiber

VSP-EM-FIBER \$88

www.vernier.com/vsp-em-fiber



Spectrum Tube Power Supplies

Spectrum Tube Single Power Supply

These power supplies feature an ultra-safe design for electrifying spectrum tubes.

ST-SPS ⚡ \$255

www.vernier.com/st-sps



Spectrum Tube Carousel Power Supply

These power supplies hold eight gas spectrum tubes.

ST-CAR ⚡ \$319


www.vernier.com/st-car



Vernier Spectral Analysis App

Our free Vernier Spectral Analysis® app makes it easy to incorporate spectroscopy into your physics lab. Using the app, students can analyze spectra from diverse sources such as spectrum tubes, light bulbs, and the sun.

www.vernier.com/spectral-analysis



Spectrum Tubes

Spectrum Tubes

Spectrum Tubes are permanently enclosed in protective plastic carriers, with no exposed high voltage.

All Spectrum Tubes are sold separately:

Hydrogen	ST-H ⚡	\$45
Nitrogen	ST-N ⚡	\$45
Helium	ST-HE ⚡	\$45
Neon	ST-NE ⚡	\$45
Carbon Dioxide	ST-CO2 ⚡	\$45
Air	ST-AIR ⚡	\$45
Argon	ST-AR ⚡	\$45

www.vernier.com/spectrum-tubes

Spectrum Tubes carry a two-year warranty (hydrogen tube: two years or 40 hours, whichever comes first; all other tubes: two years or 100 hours, whichever comes first).

Radiation Monitors

Vernier Radiation Monitor

LQ

The Vernier Radiation Monitor detects alpha, beta, gamma, and X-ray radiation and can be used for experiments in nuclear counting statistics, shielding, and decay rate measurements.

VRM-BTD \$180



Go Direct Radiation Monitor

GDX

Explore radiation statistics, measure the rate of nuclear decay, and monitor radon progeny. Go Direct Radiation Monitor detects alpha, beta, gamma, and X-ray radiation, and it includes LED and audible indicators.

GDX-RAD \$179




www.vernier.com/radiation-monitors

Nuclear Radiation with Vernier

This free e-book includes six experiments for data collection with a radiation monitor:

- Distance and Radiation
- Counting Statistics
- Lifetime Measurement
- Background Radiation Sources
- Radiation Shielding
- Alpha, Beta, and Gamma



FREE DOWNLOAD

www.vernier.com/nrv

Lab Books

NEW Vernier Video Analysis: Motion and Sports

This new e-book features 12 investigations using the Vernier Video Analysis™ app covering common concepts such as velocity and acceleration, as well as analysis of sports activities.

Download only: HSB-VVAMS-E \$25

NEW Sensor Cart Physics

Students use the Vernier Go Direct® Sensor Cart to complete the 21 investigations in this new e-book—providing a stimulating structure to explore introductory through AP* physics concepts.

Download only: HSB-SCP-E \$30

Physics with Vernier

This book features 35 experiments in mechanics, sound, light, electricity, and magnetism, using Vernier motion detectors, force sensors, light sensors, and more.

Download only: PWV-E \$40

Printed book + download: PWV \$48

Advanced Physics with Vernier—Mechanics

Advanced Physics with Vernier—Mechanics and *Advanced Physics with Vernier—Beyond Mechanics* is a two-volume set of experiments for more in-depth introductory physics courses, such as college physics, AP* Physics, and IB† Physics.

Download only: PHYS-AM-E \$40

Download only: PHYS-ABM-E \$40

Printed book + download: PHYS-AM \$48

Printed book + download: PHYS-ABM \$48

Advanced Physics with Vernier—Beyond Mechanics

Physics Explorations and Projects

Physics Explorations and Projects is a collection of investigations aligned to the NGSS. These investigations invite students to explore phenomena without extensive instructions. The guided-inquiry format involves students having some choice in what they measure and analyze.

Download only: PEP-E \$40

Printed book + download: PEP \$48

* AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

† The IB Diploma Program is an official program of the International Baccalaureate Organization (IBO) which authorizes schools to offer it. The material available here has been developed independently of the IBO and is not endorsed by it.

Learn more at www.vernier.com/lab-books

Pivot Interactives



Students overlay measurement tools onto high-quality videos to make measurements, such as in this activity where students calculate torque.



Deepen Student Understanding with Pivot Interactives

Pivot Interactives provides students with instant access to a robust collection of web-based interactive video exercises.

Each activity consists of student-controlled videos that allow variation of experimental parameters one at a time. Each video exercise challenges students to answer open-ended questions, collect their own data, and develop a mathematical model that describes the relationship between the variables.

Subscriptions start at \$5 per student (10 student minimum).

Features

- Classroom-ready experiments with teacher guides and grading/feedback tools
- Libraries (or matrices) of videos for each topic in introductory physics
- Web-based access on computers, Chromebooks, and mobile devices



See Pivot Interactives in Action




Watch a video

Start a free 30-day trial today at www.pivotinteractives.com

Vernier Video Analysis



Students can capture their own videos for analysis using mobile devices.



Investigate projectile motion

Study Motion Everywhere

The Vernier Video Analysis app brings video analysis to your students in an easy-to-use, streamlined application.

Benefits

Vernier Video Analysis makes studying motion easy and accessible. Students can use it to analyze their own recorded videos as the subject of their scientific investigations. This app brings video analysis to all your students regardless of device—it even works with Chromebooks!

Features

- Vernier Video Analysis is compatible with multiple devices and platforms: macOS®, iPadOS®, iOS, Windows® 10, Chrome OS™, and Android™.
- Students can use prepared videos, found videos, or collect their own videos for analysis.
- The app makes it possible to do experiments that cannot be done with sensors, such as analyzing the motion of a basketball in flight—objects can be tracked automatically by the app.
- Analysis is easy with multiple graphing options, so students are able to think critically about the collected data—they can even analyze the motion of multiple objects in a single video.
- With this app, you can apply vectors and vector components over the video after tracking a moving object, illuminating changes in position, velocity, and acceleration.
- When multiple objects have been marked, just enter their masses and the app can automatically calculate and display the center of mass location.
- Annual site-licensing makes purchasing and renewing quick and easy.

NEW

Vernier Video Analysis: Motion and Sports

Vernier Video Analysis: Motion and Sports features 12 investigations using Vernier Video Analysis. In addition to traditional physics concepts such as velocity and acceleration, its investigations of sports activities expand learning opportunities and further connect the study of motion to students' daily lives.

Download only
HSB-VVAMS-E \$25




Free 30-Day Trial

Get a 30-day free trial and learn about site license options at www.vernier.com/video-analysis

Physics Go Direct Package GDX

12 Products • GDP-PHY-DX • \$883
Buy 8 or more packages at \$857 and save \$208



This package includes


Go Direct® Motion	Go Direct Force and Acceleration	Go Direct Voltage	Go Direct Current (×2)
Go Direct Photogate	Ultra Pulley Attachment	Picket Fence	Go Direct Acceleration
Go Direct Sound	Go Direct Light and Color	Go Direct 3-Axis Magnetic Field	

All sensors work with our free Vernier Graphical Analysis™ app,
as well as Graphical Analysis Pro and LabQuest® 3.

Learn more at www.vernier.com/gdp-phy-dx

LabQuest 3 Physics Standard Package LQ

13 Products • LQ3-PHY-DX • \$1,082
Buy 8 or more packages at \$1,050 and save \$256



This package includes

Vernier LabQuest 3 Interface	Motion Detector	Go Direct Force and Acceleration	Differential Voltage Probe
Current Probe (×2)	Go Direct Photogate	Ultra Pulley Attachment	Picket Fence
Go Direct Acceleration	Go Direct Sound	Light Sensor	Go Direct 3-Axis Magnetic Field

All sensors work with our free Vernier Graphical Analysis app,
as well as Graphical Analysis Pro and LabQuest 3.

Learn more at www.vernier.com/lq3-phy-dx

More packages available online at www.vernier.com/physics-packages

Featured Products

Go Direct Sensors

Sensor	Order Code	Price
Go Direct 3-Axis Magnetic Field	GDX-3MG	\$69
Go Direct Acceleration	GDX-ACC	\$99
Carts and Tracks		
Dynamics Cart and Track System with Go Direct Sensor Carts	DTS-GDX	\$535
Go Direct Sensor Cart (Green)	GDX-CART-G	\$169
Go Direct Sensor Cart (Yellow)	GDX-CART-Y	\$169
Go Direct Centripetal Force Apparatus	GDX-CFA	\$299
Go Direct Current	GDX-CUR	\$79
Go Direct Force and Acceleration	GDX-FOR	\$99
Go Direct Gas Pressure	GDX-GP	\$89
Go Direct Light and Color	GDX-LC	\$79
Go Direct Motion	GDX-MD	\$99
Go Direct Photogate	GDX-VPG	\$89
Go Direct Projectile Launcher	GDX-PL	\$449
Go Direct Radiation Monitor	GDX-RAD	\$179
Go Direct Rotary Motion	GDX-RMS [⚡]	\$179
Go Direct Sound	GDX-SND	\$89
Go Direct Static Charge	GDX-Q	\$99
Temperature Probes		
Go Direct Surface Temperature	GDX-ST	\$79
Go Direct Temperature	GDX-TMP	\$69
Go Direct Voltage	GDX-VOLT	\$69

Go Direct Charge Station

Sensor	Order Code	Price
Go Direct Charge Station	GDX-CRG	\$69

LabQuest Sensors

Sensor	Order Code	Price
Accelerometers		
3-Axis Accelerometer	3D-BTA	\$99
25-g Accelerometer	ACC-BTA	\$96
Low-g Accelerometer	LGA-BTA	\$89
Carts and Tracks		
Dynamics Cart and Track System with Motion Encoder	DTS-EC	\$445
Encoder Fan Cart	CART-FEC	\$225
Current Sensors		
Current Probe	DCP-BTA	\$39
High Current Sensor	HCS-BTA	\$79
Electricity and Magnetism Sensors		
Charge Sensor	CRG-BTA	\$79
Magnetic Field Sensor	MG-BTA	\$58
Force Sensors		
Dual-Range Force Sensor	DFS-BTA	\$109
Force Plate	FP-BTA	\$289
Gas Pressure Sensor	GPS-BTA	\$89
Light Sensors		
Diffraction Apparatus	DAK	\$620
Light Sensor	LS-BTA [⚡]	\$59
Motion Detectors		
Go!Motion® (USB sensor)	GO-MOT	\$129

Looking for Replacement Parts?

Visit www.vernier.com/replacements

Motion Detector	MD-BTD	\$89
Photogate	VPG-BTD	\$49
Power Amplifier	PAMP	\$225
Projectiles		
Projectile Launcher	VPL	\$389
Time of Flight Pad	TOF-VPL	\$84
Radiation Monitor	VRM-BTD	\$180
Rotary Motion Sensor	RMV-BTD [⚡]	\$169
Sound Sensors		
Microphone	MCA-BTA	\$44
Sound Level Sensor	SLS-BTA	\$69
Temperature Probes		
Stainless Steel Temperature Probe	TMP-BTA	\$36
Surface Temperature Sensor	STS-BTA	\$25
Voltage Probes		
30-Volt Voltage Probe	30V-BTA	\$49
Differential Voltage Probe	DVP-BTA	\$39
Instrumentation Amplifier	INA-BTA	\$79
Voltage Probe	VP-BTA	\$12

Emissions Spectrometer

Spectrometer	Order Code	Price
Vernier Emissions Spectrometer	VSP-EM	\$799

Infrared Cameras

Camera	URL
FLIR ONE® Thermal Cameras	www.vernier.com/flir

See all our products for physics at www.vernier.com/physics



Engineering, Coding, and Robotics

www.vernier.com/engineering

Encourage curiosity, build confidence, and spark an interest in STEM careers in your students. Vernier solutions give your students practical ways to learn engineering design principles, integrate sensor data into computer science concepts, and learn coding with robotics.

Topics

Explore a sampling of our featured experiments and investigations by topic to learn how Vernier technology helps your students engage with data-collection technology and deepens their understanding of key engineering, computer science, and STEM concepts.

Engineering


PAGE 124



Bridge and Structure Testing



Renewable Energy



Arduino®

Coding with Sensors

PAGE 128



Block-Based



Python®



JavaScript™

Robotics

PAGE 129



LEGO® Robotics



Makeblock®

Professional Development

We are here to help. Our virtual professional development workshops, webinars, and personalized online training options offer innovative ways to engage students with STEM in a traditional classroom or virtual environment.

www.vernier.com/training

Our solutions help your students understand the engineering design process, critical thinking, and teamwork. Your students learn to build and design bridges, wind turbines, and more. Plus, our world-class technical support ensures success in the classroom.

Coding introduces problem solving, nurtures creativity, increases critical thinking, and builds confidence. We have added coding support to our Go Direct® sensors so that your students can develop computational thinking as they learn to code.

When your students design robots and develop code, they express ideas in new ways. With robotics, your students learn skills extending beyond the screen as they program robots to interact with the physical world.

Bridge and Structure Testing


FEATURED ACTIVITY

Bridge Competition

In this team competition, students use the engineering design process to design a bridge with the highest efficiency, following a set of constraints and design requirements.



Equipment Used



Go Direct® Structures & Materials Tester

Use our new Go Direct Structures & Materials Tester to evaluate the strength of model bridges and engineered structures by measuring the applied load. Utilizing both load and displacement sensors, your students can evaluate the properties of materials.

Benefits

- Force and displacement sensors connect via Bluetooth® wireless technology or via USB
- Uses our free Vernier Graphical Analysis™ app or Graphical Analysis Pro to collect and analyze data
- Exact force and displacement for bends and breaks
- Accurate positioning for center and off-center loading
- Easy loading for different sizes and shapes
- Includes *Materials Testing: Beams to Bridges* e-book

GDX-VSMT \$999

Activity Source



Materials Testing: Beams to Bridges with Go Direct Structures & Materials Tester

GDXVSMT-BB-E \$20*

*Free with purchase of Go Direct Structures & Materials Tester

Learn more at www.vernier.com/gdxvsmt-bb-e

Materials Testing: Beams to Bridges with Go Direct Structures & Materials Tester

With the activities in this e-book, students use the Go Direct Structures & Materials Tester to investigate materials and structures.

Topics include

- Beams: Investigate the relationship between dimensions and flexibility.
- Trusses: Explore why trusses fail and how to compensate for weaknesses.
- Bridges: Use the engineering design process to build and test bridges.

www.vernier.com/gdxvsmt-bb-e


Truss Tester Accessory

The Truss Tester Accessory attaches to the Go Direct Structures & Materials Tester, holds a single truss upright, and allows the load to be applied in a variety of locations.

VSMT-TRUSS \$128

www.vernier.com/vsmt-truss

INCLUDES
5
ACTIVITIES



GDXVSMT-BB-E \$20†

†Free with purchase of Go Direct Structures & Materials Tester



PLTW Engineering

PLTW Engineering (9–12) empowers students to step into the role of an engineer and adopt a problem-solving mindset, inspiring students to believe in their own potential and see themselves in a career that improves communities.

Learn more at www.vernier.com/pltw



Renewable Energy

FEATURED EXPERIMENT

Project: Maximum Energy Output

Challenge your students to design their own wind turbines following the provided design requirements, constraints, and deliverables.



Sensor Used



Go Direct Energy

Use Go Direct Energy with our free Vernier Graphical Analysis app or Graphical Analysis Pro to determine the power output of a renewable energy system. Connect a source, such as KidWind solar panels or wind turbines, and students can quantitatively evaluate the effects of their design changes.

GDX-NRG \$89

Accessory Used



Vernier Variable Load

The Vernier Variable Load provides a range of resistive loads for projects with wind turbines or solar panels. This load is used in our *Renewable Energy with Vernier* lab book.

VES-VL \$64

Experiment Source



Renewable Energy with Vernier

Download only: REV-E \$40

Printed book + download: REV \$48


Learn more at www.vernier.com/rev-15

Renewable Energy with Vernier

The *Renewable Energy with Vernier* lab book features 26 experiments in wind and solar energy. The book contains a combination of explorations, classic experiments, inquiry investigations, engineering projects, and more.

Learn more at www.vernier.com/rev

INCLUDES
26
EXPERIMENTS



Download only

REV-E \$40

Download + printed book

REV \$48

Additional Products

KidWind Advanced Wind Experiment Kit

Discover advanced aspects of wind turbine technology. Test different blade designs, gear ratios, generators, and devices to measure electrical and weightlifting power.

KW-AWX \$154

More KidWind renewable energy products can be found at

www.vernier.com/kidwind



Engineering

Arduino

FEATURED PROJECT

Functions

This activity uses Arduino® to introduce students to the concept of functions. Students explore how functions can make their Arduino code more efficient and easier to understand. Students also learn formatting for creating and calling a function and how to distinguish between local and global variables.



Products Used



Gas Pressure Sensor

Use the Gas Pressure Sensor with an Arduino microcontroller to introduce the basics of sensor technology.

GPS-BTA \$89



Vernier Arduino® Interface Shield

The Vernier Arduino Interface Shield provides a convenient way to make connections from Arduino microcontrollers to Vernier sensors.

BT-ARD \$29



SparkFun® RedBoard with Cable

The SparkFun RedBoard is an Arduino-compatible board, which is perfect for use with the Vernier Arduino Interface Shield.

ARD-RED \$25

Project Source



Vernier Coding Activities with Arduino: Analog Sensors

VCA-AS-E \$20*

*Free with the purchase of the Vernier Coding with Arduino—Analog Sensor Package or the Vernier Arduino Interface Shield

Learn more at www.vernier.com/arduino

NEW Vernier Coding Activities with Arduino: Analog Sensors

The activities in this e-book provide an introduction to coding and sensor technology using Vernier sensors and Arduino microcontrollers. Teaching students about microcontrollers and sensors opens the door for them to explore how technology and coding affect the world beyond the screen. This e-book is available for individual purchase or is free with the purchase of the Vernier Interface Shield. It is also included with the purchase of the Vernier Coding with Arduino—Analog Sensor Package.

VCA-AS-E \$20†


†Free with purchase of the Vernier Coding with Arduino—Analog Sensor Package or the Vernier Arduino Interface Shield

www.vernier.com/arduino

Vernier Coding Activities with Arduino®

Analog Sensors

VCA-AS-E



INCLUDES
8
ACTIVITIES

NEW Vernier Coding with Arduino—Analog Sensor Package


This package has all the equipment and activities you need to get students started using Vernier sensors with Arduino microcontrollers. The package includes the new Vernier Coding Activities with Arduino: Analog Sensors e-book at no additional cost.

This package includes

- Gas Pressure Sensor
- Vernier Arduino Interface Shield
- SparkFun® RedBoard with Cable
- Vernier Coding Activities with Arduino: Analog Sensors

VCA-AS-PKG \$143

Learn more at www.vernier.com/vca-as-pkg













Featured Products

Bridge and Structure Testing

Product	Order Code	Price
Go Direct® Structures & Materials Tester	GDX-VSMT	\$999
Truss Tester Accessory	VSMT-TRUSS	\$128
Materials Testing: Beams to Bridges with the Go Direct Structures & Materials Tester lab book	GDXVSMT-BB-E	\$20 [†]

Arduino

Product	Order Code	Price
SparkFun RedBoard with Cable	ARD-RED	\$25
Vernier Arduino Interface Shield	BT-ARD	\$29
Vernier Coding Activities with Arduino: Analog Sensors lab book	VCA-AS-E	\$20 [†]


Anemometer		ANM-BTA	\$89
Dual-Range Force Sensor		DFS-BTA	\$109
Gas Pressure Sensor		GPS-BTA	\$89
Light Sensor		LS-BTA [‡]	\$59
Motion Detector		MD-BTD	\$89
pH Sensor		PH-BTA	\$88
Photogate		VPG-BTD	\$49
Soil Moisture Sensor		SMS-BTA	\$109
Stainless Steel Temperature Probe		TMP-BTA	\$36
Surface Temperature Sensor		STS-BTA	\$25

Renewable Energy

Product	Order Code	Price
Go Direct Energy	GDX-NRG	\$89
Vernier Variable Load	VES-VL	\$64
KidWind Advanced Wind Experiment Kit	KW-AWX	\$154
KidWind Balsa Blade Sheets	KW-BBS10	\$12
KidWind Wind Turbine Generator with Wires	KW-GEN	\$7
KidWind Tower and Base Set	KW-TBS	\$24
KidWind Basic Turbine Building Parts	KW-BTPART	\$16
Renewable Energy with Vernier lab book	Printed book + download: REV	\$48
	Download only: REV-E	\$40

See all of our products for engineering at www.vernier.com/engineering

Coding with Go Direct Sensors



Coding with Go Direct® Sensors

Vernier offers a range of coding solutions—from entry-level to advanced instrument-control programming. With Vernier technology and an appropriate coding application, your students can create code to visualize scientific data, incorporate sensor input, and create sensor-controlled projects.

Learn more at www.vernier.com/hs-engineering



Block-Based Coding

Scratch

Block-based programming is ideal for students new to coding. With Scratch, students can develop their coding skills with fun hands-on projects. Block-based coding in Scratch helps students get started making natural connections between their digital and physical worlds.

Learn more at www.vernier.com/scratch



Connecting to Python®

With our Python getting started guide and examples, you can connect Vernier Go Direct sensors to your Python project. Your students can write Python programs to visualize Go Direct sensor data or integrate that data into a larger Python project.



Using JavaScript™

Use JavaScript to integrate Go Direct sensor data into your custom web applications. Integrate coding, sensor data collection, and web design by combining the Vernier Go Direct library with other libraries including Chart.js, Desmos.js, and p5.js.





Computer Science

PLTW Computer Science (9–12) engages students in real-world activities and projects that challenge them to apply computational thinking and logic to solve big problems.

Learn more at www.vernier.com/pltw

Robotics



Vernier Robotics

When your students learn to program robots, they learn to organize, express, and share their ideas in a whole new way. With robotics in the classroom, your students learn coding skills that extend beyond the screen as they program robots to interact with the physical world.

We recognize that educators partner with dependable providers that they know and love. We strive to do the same, which is why we work with LEGO® Education and Makeblock. Boost your students' understanding of robotics concepts with downloadable e-books that incorporate problem-solving, engineering design, and critical thinking skills.

LEGO® Education

LEGO® MINDSTORMS® Education EV3 is a hands-on cross-curricular robotics STEM solution that engages students by providing the resources to design, build, and program their creations while helping them develop essential 21st-century skills such as creativity, critical thinking, collaboration, and communication.

See page 130.



Makeblock®

Help your students learn how to organize, express, and share their ideas in a whole new way through coding. With Makeblock robots coupled with exclusive STEM activities from Vernier, your students learn coding skills as they program robots to interact with the physical world.

See page 131.



Learn more at www.vernier.com/coding-robotics

LEGO® Education

LEGO® MINDSTORMS® Education EV3 Core Set with Charger*

LEGO® MINDSTORMS® Education EV3 Core Set is a hands-on cross-curricular STEM solution that engages students by providing the resources to design, build, and program their creations. A Core Set supports two students as they practice collaboration, communication, and critical thinking. The software is Windows®, macOS®, Chrome OS™, iPadOS®, and iOS compatible.

The LEGO® MINDSTORMS® Education EV3 kit includes 541 elements, including an EV3 brick, interactive servo-motors, gears, sensors, and wheels, that can be used for teaching science, technology, engineering, math, and computer science.

LEGO-EV3-CORE \$439.90

www.vernier.com/lego-ev3-core



LEGO® MINDSTORMS® Education EV3 Expansion Set*

The Expansion Set contains a wide range of structural and mechanical elements to augment the LEGO® MINDSTORMS® Education EV3 Core Set. Students can deepen their experience with the additional building programs and instructions.

LEGO-EV3-EXP \$119.95

www.vernier.com/lego-ev3-exp



NXT Sensor Adapter for EV3 and NXT

The Vernier NXT Sensor Adapter allows certain Vernier LabQuest® sensors to work on the LEGO® MINDSTORMS® EV3 and LEGO® MINDSTORMS® NXT robotics systems. Enhance your robots with sensors for measuring everything from temperature to force, light level, UV level, pH, and more.

BTA-NXT \$39

www.vernier.com/bta-nxt



* Vernier Software & Technology is an authorized sales agent of LEGO® Brand Retail Inc. d/b/a LEGO Education North America. LEGO® Education North America will operationally fulfill all LEGO® Education products ordered through Vernier.

Robotics Makeblock®

ACTIVITY C7

Security and Emergency Systems

Students program and troubleshoot their block-based code in order to create a car alarm, a security system, and an emergency vehicle (with warning lights and sirens) for their mBot robot.



Accessory Used



This activity can also be done with

mBot Explorer

MBOT-S \$83.99

mBot by Makeblock

mBot provides students with a fun and tactile way to learn entry-level coding with simple Scratch-based software. Included with your purchase is our *Coding with mBot: Self-Driving Vehicles* e-book.

MBOT-P (pink) or MBOT-B (blue) \$73.49 each

Experiment Source



Coding with mBot: Self-Driving Vehicles

MBOT-MSDV-E \$20*

*Free with purchase of mBot from Vernier

Learn more at www.vernier.com/mbot-msdv-e-c7

mBot STEM Classroom Kit and Coding with mBot: Life Hacks

With the mBot™ STEM Classroom Kit, a complete robotics kit, students can create anything—from a simple mBot to a complex robotics system. Your purchase includes our comprehensive STEM activities e-book with lessons that have students solve a number of practical problems using robotics and coding.

The kit includes

- Blue mBot robot
- Perception Gizmos Add-on Pack
- Variety Gizmos Add-on Pack
- *Coding with mBot: Life Hacks* e-book

MBOT-SKIT \$168 www.vernier.com/mbot-skit



mBot Ranger by Makeblock

mBot Ranger is a STEM robot kit that can be constructed into three unique designs, like an off-road tank or a spinning raptor, for a wider range of learning.

MB-RANGER \$157.49

www.vernier.com/mb-ranger



Makeblock Accessories

Products	Order Code	Price
Add-on Packs for mBot		
mBot Servo Pack Add-on Pack	MBOT-SERVO	\$27
mBot Interactive Light and Sound Add-on Pack	MBOT-LS	\$27
Perception Gizmos Add-on Pack	MBOT-PER	\$52.49
Variety Gizmos Add-on Pack	MBOT-VAR	\$41.99
mBot Six-Legged Robot Add-on Pack	MBOT-6LR	\$27
Makeblock Bluetooth® Dongle	MB-BLE	\$15.99
mBot 3.7 V LiPo Battery	MBOT-BAT	\$9.99
Me 7-Segment Serial Display	MBOT-DSPL	\$9
Me LED Matrix 8 × 16	MBOT-MTRX	\$13

TI-Nspire™ CX II Handheld

TI-Nspire CX II handheld is the latest in learning technology from Texas Instruments. The handheld includes an easy-glide touchpad that works like a computer with a mouse.

Recommended for algebra, geometry, trigonometry, and precalculus

Includes TI-Nspire CX II handheld, rechargeable battery, slide cover, and unit-to-computer USB connectivity and charging cable


TI-NSCX2 \$136

TI-Nspire CX II Teacher Pack

Includes 10 TI-Nspire CX II EZ-Spot handhelds with the words "School Property" on the keypad, 10 rechargeable batteries, and a 10-unit docking station

TI-NSCX2-TPK \$1,479

Learn more at www.vernier.com/ti-nsx2



CBR 2™

The CBR 2 connects directly to a TI calculator. This motion detector collects distance, velocity, and acceleration data.

CBR2 \$99

Learn more at www.vernier.com/cbr2



TI-Nspire CX II CAS Handheld

TI-Nspire CX II CAS handheld has all the features of the TI-Nspire CX II handheld plus a built-in Computer Algebra System (CAS) for factoring and expanding expressions, solving for common denominator, and performing symbolic calculations.

Recommended for geometry, trigonometry, precalculus, and calculus

Includes TI-Nspire CX II CAS handheld, rechargeable battery, slide cover, and unit-to-computer USB connectivity and charging cable


TI-NSXCAS2 \$139

TI-Nspire CX II CAS Teacher Pack

Includes 10 TI-Nspire CX II CAS handhelds, 10 rechargeable batteries, and a 10-unit docking station

TI-NSXCAS2-TPK \$1,509

Learn more at www.vernier.com/ti-nsxcas2



Vernier EasyTemp®

EasyTemp is a temperature probe designed for use with TI-84 Plus calculators and TI-Nspire handhelds.

Range: -20 to 115°C

EZ-TMP \$38

Learn more at www.vernier.com/ez-tmp



Vernier EasyLink®

EasyLink is a single-channel sensor interface that plugs into the USB port of a TI-84 Plus calculator or TI-Nspire handheld. It supports any one of over 60 Vernier sensors.

EZ-LINK \$67

Learn more at www.vernier.com/ez-link



TI-84 Plus CE

The TI-84 Plus CE has a full-color, high-resolution, backlit screen, making it easy to read. The calculator comes with a rechargeable battery, so there is never a need to buy AAA batteries.

- Supported USB sensors: CBR 2 and EasyTemp
- Supported interface: EasyLink

Includes TI-84 Plus CE calculator, rechargeable battery, unit-to-computer connectivity and charging cable, slide cover, and AC wall adapter

TI-84PCE \$129

TI-84 Plus CE Teacher Pack

Includes 10 TI-84 Plus CE EZ-Spot calculators, 10 rechargeable batteries, and a 10-unit charging station

TI-84PCE-TPK \$1,345

Learn more at www.vernier.com/ti-84pce



TI-84 Plus

The TI-84 Plus is a lower-price alternative to the TI-84 Plus CE calculator. The TI-84 Plus supports data collection with 78 Vernier sensors, including microphones, photogates, and drop counters, when used with a CBL 2™ sensor interface.

- Supported USB sensors: CBR 2 and EasyTemp
- Support interfaces: EasyLink and CBL 2
- Collect data from multiple sensors simultaneously with CBL 2.


TI-84PL \$108

TI-84 EZ-Spot Teacher Pack

Includes 10 TI-84 Plus EZ-Spot calculators and 40 AAA batteries

TI-84PL-TPK \$1,080

Learn more at www.vernier.com/ti-84pl



Calculator Products

Product	Order Code	Price
Books	<i>Real-World Math with Vernier</i> (download only)	RWV-E \$40
	TI-84 Plus CE	TI-84PCE \$129
	TI-84 Plus CE Teacher Pack (10 EZ-Spot calculators & charging station)	TI-84PCE-TPK \$1,345
	TI-84 Plus Calculator	TI-84PL \$108
Calculators	TI-84 Plus EZ-Spot Teacher Pack (10 EZ-Spot calculators)	TI-84PL-TPK \$1,080
	TI-Nspire CX II Handheld	TI-NSCX2 \$136
	TI-Nspire CX II Teacher Pack (10 EZ-Spot handhelds & docking station)	TI-NSCX2-TPK \$1,479
	TI-Nspire CX II CAS Handheld	TI-NSCXCAS2 \$139
	TI-Nspire CX II CAS Teacher Pack (10 handhelds & docking station)	TI-NSCXCAS2-TPK \$1,509
	TI-84 Plus CE Charging Station	TI-84PCE-CS \$70
Charging/ Docking Station	TI-Nspire CX Docking Station	TI-NSCX-DS \$120
	CBL 2*	CBL2 \$166
Data Collection	EasyLink	EZ-LINK \$67
	EasyTemp	EZ-TMP \$38
	CBR 2	CBR2 \$99
Emulator/ Computer Software	TI-SmartView™ Emulator software for TI-84	
	TI-Nspire CX Student Software	www.vernier.com/ti-software
	TI-Nspire CX Premium Teacher Software	
Miscellaneous Accessories	Easy to Go! USB Adapter	MINI-USB \$17
	Go! to Easy USB Adapter	USB-MINI \$9
TI Navigator System	30-User TI-Nspire CX Navigator System	TI-NAV-CX30 \$2,025
	10-User TI-Nspire CX Navigator System	TI-NAV-CX10 \$1,160
	10-User TI-Nspire CX Navigator Add-on†	TI-NAV-10ADDON \$756

TI products purchased in the United States are covered by a one-year warranty based on the date of purchase. Units are warranted against defective materials or workmanship.

*Cannot be used with TI-84 Plus CE or TI-Nspire handhelds

†Requires purchase of a Navigator system

Sensors & Accessories

The Vernier Sensor Advantage

Outstanding Performance

With 40 years of experience developing technology for education, we design our sensors for active, hands-on experiments. Vernier sensors are rugged, classroom-proven technology that are well supported and easy to use. The sensors provide consistent, high-quality results for the demands of the classroom.

Connect & Collect

Simply connect, and you're ready to collect. All Vernier sensors on the following pages are automatically detected and set up for data collection when used with Vernier software.

Go Direct Sensors

Our Go Direct® sensors connect directly to a computer, Chromebook™, or a mobile device via Bluetooth® wireless technology or USB connection. Most sensors include a rechargeable battery to power the sensor when used wirelessly.

LabQuest Sensors

Our LabQuest® sensors require an interface from the LabQuest family, such as LabQuest 3, LabQuest Stream®, or LabQuest Mini. The interface sends information from the sensor to the data-collection and analysis software on a device such as a computer, Chromebook, or mobile device.

For more information on sensor compatibility, visit www.vernier.com/sensors

Generous Warranty

Buy with confidence. Most Vernier sensors are covered by a 5-year limited warranty. During the warranty period, Vernier will repair or replace the item if there is a defect in materials or workmanship. Outside the warranty, Vernier will attempt to repair most products, often at no charge.

Go Direct Sensors

Sensor	Order Code	Price
Go Direct 3-Axis Magnetic Field	GDX-3MG	\$69
Go Direct Acceleration	GDX-ACC	\$99
Go Direct Blood Pressure	GDX-BP	\$105
Carts and Tracks		
Dynamics Cart and Track System with Go Direct Sensor Carts	DTS-GDX	\$535
Go Direct Sensor Cart (Green)	GDX-CART-G	\$169
Go Direct Sensor Cart (Yellow)	GDX-CART-Y	\$169
Go Direct Centripetal Force Apparatus (requires Go Direct Force and Acceleration)	GDX-CFA	\$299
Go Direct CO ₂ Gas	GDX-CO2	\$199
Go Direct Colorimeter	GDX-COL	\$119
Conductivity Probes		
Go Direct Conductivity	GDX-CON	\$99
Go Direct Platinum-Cell Conductivity	GDX-CONPT	\$169
Go Direct Constant Current System	GDX-CCS	\$74
Go Direct Current	GDX-CUR	\$79
Go Direct Drop Counter	GDX-DC	\$99
Go Direct EKG	GDX-EKG	\$159
Go Direct Electrode Amplifier	GDX-EA	\$64
Go Direct Energy	GDX-NRG	\$89
Go Direct Ethanol Vapor	GDX-ETOH	\$149
Go Direct Force and Acceleration	GDX-FOR	\$99
Go Direct Gas Pressure	GDX-GP	\$89
Go Direct Hand Dynamometer	GDX-HD	\$109
Heart Rate Monitors		
Go Wireless Exercise Heart Rate	GW-EHR	\$79
Go Wireless Heart Rate	GW-HR	\$89
Go Direct Ion-Selective Electrode Amplifier	GDX-ISEA	\$69

** Ion-Selective Electrodes require excellent chemical technique and careful calibration to obtain accurate results; they are not recommended for elementary or middle school students.*

Ion-Selective Electrodes (ISE)*

Go Direct Ammonium ISE	GDX-NH4	\$249
Go Direct Calcium ISE	GDX-CA	\$249
Go Direct Chloride ISE	GDX-CL	\$249
Go Direct Nitrate ISE	GDX-NO3	\$249
Go Direct Potassium ISE	GDX-K	\$249
Go Direct Light and Color	GDX-LC	\$79
Go Direct Melt Station	GDX-MLT ▲	\$529
Go Direct Motion	GDX-MD	\$99
Go Direct Mini GC	GDX-GC	\$2,499
Go Direct O ₂ Gas	GDX-O2	\$189
Go Direct Optical Dissolved Oxygen	GDX-ODO	\$298
Go Direct ORP	GDX-ORP	\$99


pH Sensors





Go Direct Glass-Body pH	GDX-GPH	\$139
Go Direct pH	GDX-PH	\$89
Go Direct Tris-Compatible Flat pH	GDX-FPH	\$115
Go Direct Photogate	GDX-VPG	\$89
Go Direct Polarimeter	GDX-POL	\$499
Go Direct Projectile Launcher	GDX-PL	\$449
Go Direct Radiation Monitor	GDX-RAD	\$179
Go Direct Respiration Belt	GDX-RB	\$99
Go Direct Rotary Motion	GDX-RMS ▲	\$179
Go Direct Sound	GDX-SND	\$89
Go Direct SpectroVis Plus	GDX-SVISPL	\$399
Go Direct Spirometer	GDX-SPR	\$199
Go Direct Static Charge	GDX-Q	\$99
Go Direct Structures & Materials Tester	GDX-VSMT	\$999

Temperature Probes

Go Direct Surface Temperature	GDX-ST	\$79
Go Direct Temperature	GDX-TMP	\$69
Go Direct Thermocouple	GDX-TC	\$109
Go Direct Wide-Range Temperature	GDX-WRT	\$114
Go Direct Voltage	GDX-VOLT	\$69
Go Direct Weather	GDX-WTHR	\$99

LabQuest Sensors

Sensor	Order Code	Price
Accelerometers		
3-Axis Accelerometer	3D-BTA	\$99
25-g Accelerometer	ACC-BTA	\$96
Low-g Accelerometer	LGA-BTA	\$89
Anemometer	ANM-BTA	\$89
Barometer	BAR-BTA	\$71
Blood Pressure Sensor	BPS-BTA	\$109
Charge Sensor	CRG-BTA	\$79
CO ₂ Gas Sensor	CO2-BTA	\$269
Colorimeter	COL-BTA	\$119
Conductivity Probes		
Conductivity Probe	CON-BTA	\$99
Platinum-Cell Conductivity Probe	CONPT-BTA	\$149
Constant Current System	CCS-BTA	\$64
Current Probes		
Current Probe	DCP-BTA	\$39
High Current Sensor	HCS-BTA	\$79
Diffraction Apparatus	DAK	\$620
Digital Control Unit	DCU-BTD	\$61
Drop Counter	VDC-BTD	\$99
EKG Sensor	EKG-BTA	\$158
Electrode Amplifier	EA-BTA	\$49
Energy Sensor	VES-BTA	\$88
Ethanol Sensor	ETH-BTA	\$119
Flow Rate Sensor	FLO-BTA	\$129
Force Sensors		
Dual-Range Force Sensor	DFS-BTA	\$109
Force Plate	FP-BTA	\$289
Gas Pressure Sensors		
Gas Pressure Sensor	GPS-BTA	\$89
Pressure Sensor 400	PS400-BTA	\$189
Goniometer	GNM-BTA 	\$159
Hand Dynamometer	HD-BTA	\$110

Heart Rate Monitors		
Exercise Heart Rate Monitor	EHR-BTA	\$99
Hand-Grip Heart Rate Monitor	HGH-BTA	\$119
Instrumentation Amplifier	INA-BTA	\$79
Ion-Selective Electrodes (ISE)*		
Ammonium ISE	NH4-BTA	\$219
Calcium ISE	CA-BTA	\$219
Chloride ISE	CL-BTA	\$219
Nitrate ISE	NO3-BTA	\$219
Potassium ISE	K-BTA	\$219
Light Sensor	LS-BTA 	\$59
Magnetic Field Sensor	MG-BTA	\$58
Melt Station	MLT-BTA 	\$519
Microphone	MCA-BTA	\$44
Motion Detectors		
Dynamics Cart and Track System with Motion Encoder	DTS-EC	\$445
Motion Detector	MD-BTD	\$89
O ₂ Gas Sensor	O2-BTA	\$199
Optical DO Probe	ODO-BTA	\$299
ORP Sensor	ORP-BTA	\$89
PAR Sensor	PAR-BTA	\$229
pH Sensors		
Glass-Body pH Electrode BNC (requires Electrode Amplifier)	GPH-BNC	\$85
pH Sensor	PH-BTA	\$88
Tris-Compatible Flat pH Sensor	FPH-BTA	\$104
Photogate	VPG-BTD	\$49
Polarimeter (Chemical)	CHEM-POL 	\$499
Power Amplifier	PAMP	\$225
Projectile Launcher	VPL	\$389
Pyranometer	PYR-BTA	\$229
Qubit Sensors	www.vernier.com/qubit	
Radiation Monitor	VRM-BTD	\$180
Relative Humidity Sensor	RH-BTA	\$69
Respiration Monitor Belt (requires Gas Pressure Sensor)	RMB	\$63
Rotary Motion Sensor	RMV-BTD 	\$169

Salinity Sensor	SAL-BTA	\$119
Soil Moisture Sensor	SMS-BTA	\$109
Sound Level Sensor	SLS-BTA	\$69
Spirometer	SPR-BTA	\$219
Temperature Probes		
Extra-Long Temperature Probe	TPL-BTA	\$99
Stainless Steel Temperature Probe	TMP-BTA	\$36
Surface Temperature Sensor	STS-BTA	\$25
Thermocouple	TCA-BTA	\$69
Wide-Range Temperature Probe	WRT-BTA	\$82
Turbidity Sensor	TRB-BTA	\$112
UV Sensors		
UVA Sensor	UVA-BTA	\$109
UVB Sensor	UVB-BTA	\$110
Voltage Probes		
30-Volt Voltage Probe	30V-BTA	\$49
Differential Voltage Probe	DVP-BTA	\$39
Voltage Probe	VP-BTA	\$12

USB-Only Sensors

Sensor	Order Code	Price
Go!Motion	GO-MOT	\$129
Go!Temp	GO-TEMP	\$39
OHAUS® Balances	www.vernier.com/ohaus	
Spectrometers		
Go Direct SpectroVis® Plus (USB and Wireless)	GDX-SVISPL	\$399
Vernier Emissions Spectrometer	VSP-EM	\$799
Vernier Flash Photolysis Spectrometer	VSP-FP	\$4,999
Vernier Fluorescence/UV-VIS Spectrophotometer	VSP-FUV	\$2,899
Vernier Spectrometer (Ocean Optics dba Ocean Insight)	V-SPEC	\$1,999
Vernier UV-VIS Spectrophotometer	VSP-UV	\$2,100

Accessories & Replacement Parts





Sensors

Part Name	Order Code	Price
Blood Pressure Sensors		
Small Blood Pressure Cuff	CUFF-SM	\$32
Standard Blood Pressure Cuff	CUFF-STD	\$30
Large Blood Pressure Cuff	CUFF-LG	\$35
CO ₂ and/or O ₂ Gas Sensors		
250 mL Nalgene® Bottle (1 opening)	CO2-BTL	\$5
BioChamber 250 (250 mL) (2 openings)	BC-250 ⚡	\$8
BioChamber 2000 (2000 mL) (2 openings)	BC-2000 ⚡	\$22
Colorimeters		
Cuvette Lids (pkg. of 100)	CUV-LID	\$9
Cuvette Rack	CUV-RACK	\$9
Plastic Cuvettes (Visible Range) (pkg. of 100)	CUV	\$19
Conductivity Probes		
Conductivity Low Standard (500 mL)	CON-LST	\$20
Conductivity Middle Standard (500 mL)	CON-MST	\$20
Conductivity High Standard (500 mL)	CON-HST	\$20
Dissolved Oxygen Probe (Go Direct®, order code GDX-ODO)		
Go Direct Optical Dissolved Oxygen Replacement Cap	GDX-ODO-CAP	\$69
Dissolved Oxygen Probe (Optical, order code ODO-BTA)		
Optical DO Probe Metal Guard	ODO-GRD	\$49
Optical DO Probe Replacement Cap	ODO-CAP	\$54
Dissolved Oxygen Probe (Non-optical, order code DO-BTA)		
DO Calibration Solution (60 mL)	DO-CAL	\$5
DO Filling Solution (130 mL)	FS	\$6
DO Polishing Strips	PS	\$4
DO Probe Membrane Cap	MEM	\$13
Drop Counters		
Microstirrer	MSTIR	\$9
Reagent Reservoir, 2 Valves, and Tip	VDC-RR	\$10
Stopper Stem	PS-STEM	\$1
Plastic 2-Way Valve	PS-2WAY	\$2
EKG Sensors		
EKG Electrodes (100)	ELEC	\$15
Electrode Amplifier (Go Direct, order code GDX-EA)		
Go Direct pH Electrode BNC	GDX-PH-BNC	\$40
Go Direct Glass-Body pH Electrode BNC	GDX-GPH-BNC	\$84



Go Direct Flat pH Electrode BNC	GDX-FPH-BNC	\$73
Go Direct ORP Electrode BNC	GDX-ORP-BNC	\$49
Electrode Amplifier (LabQuest, order code EA-BTA)		
pH Electrode BNC	PH-BNC	\$41
Glass-Body pH Electrode BNC	GPH-BNC	\$85
Flat pH Electrode BNC	FPH-BNC	\$74
ORP Electrode BNC	ORP-BNC	\$48
Energy Sensors		
Vernier Resistor Board	VES-RB	\$18
Vernier Variable Load	VES-VL	\$64
Ethanol Sensors		
Ethanol Cap Assemblies (pkg. of 3)	ETH-CAPS	\$10
Ethanol Stopper	ETH-STOP	\$4
Ethanol Tape	ETH-TAPE	\$3
Force Sensors		
Reflex Hammer Accessory Kit	RFX-ACC	\$29
Replacement Accessory Rod	ACC-ROD	\$4
Springs Set	SPRINGS	\$18
Dual-Range Force Sensor Replacement Parts Kit	DFS-RPK	\$24
Bumper Launcher Kit	BLK	\$89
Hoop Bumpers for Bumper and Launcher Kit	HOOPS-BLK	\$15
Gas Chromatographs		
GC Septa (pkg. of 4)	GC-SEP	\$29
GC Syringe, 1 µL Hamilton	GC-SYR-MIC	\$90
Gas Pressure Sensors		
Gas Pressure Sensor Bulb (1)	GPS-BULB1	\$6
Gas Pressure Sensor Bulb (set of 4)	GPS-BULB4	\$21
Pressure Sensor Accessories Kit	PS-ACC	\$12
#11-Hole Rubber Stopper	PS-STOP1	\$1
#5 2-Hole Rubber Stopper	PS-STOP5	\$1.50
Luer-Lock Connector	PS-LUER	\$1
Plastic 2-Way Valve	PS-2WAY	\$2
Plastic Tubing	PS-TUBING	\$1
Plastic Tubing Clamps (pkg. of 100)	PTC	\$49
Stopper Stem	PS-STEM	\$1
Syringe (20 mL, plastic)	PS-SYR	\$2
Syringe (20 mL, plastic) (pkg. of 10)	PS-SYR10	\$18

Heart Rate Sensors		
Heart Rate Hand Grips	HR-GRIP	\$31
Exercise Heart Rate Strap	HR-STRAP	\$21
Polar Transmitter Module	HR-TRANS	\$58
Ion-Selective Electrodes		
ISE Ammonium Replacement Module†	NH4-MOD	\$79
ISE Calcium Replacement Module†	CA-MOD	\$79
ISE Nitrate Replacement Module†	NO3-MOD	\$79
ISE Potassium Replacement Module†	K-MOD	\$79
ISE Ammonium Low Standard (500 mL)	NH4-LST	\$20
ISE Ammonium High Standard (500 mL)	NH4-HST	\$20
ISE Calcium Low Standard (500 mL)	CA-LST	\$20
ISE Calcium High Standard (500 mL)	CA-HST	\$20
ISE Chloride Low Standard (500 mL)	CL-LST	\$20
ISE Chloride High Standard (500 mL)	CL-HST	\$20
ISE Nitrate Low Standard (500 mL)	NO3-LST	\$20
ISE Nitrate High Standard (500 mL)	NO3-HST	\$20
ISE Potassium Low Standard (500 mL)	K-LST	\$20
ISE Potassium High Standard (500 mL)	K-HST	\$20
Melt Stations		
Melt Station Capillary Tubes (pkg. of 100)	MLT-TUBE	\$19
Motion Detectors		
Go!Motion to Computer Cable	GMC-USB	\$5
Motion Detector Cable	MDC-BTD	\$5
Motion Detector Clamp	MD-CLAMP	\$15
pH and ORP Sensors		
Microstirrer	MSTIR	\$9
pH Buffer Capsules (10 each of pH 4, 7, 10)	PH-BUFCAP	\$29
pH Storage Bottles (pkg. of 5)	BTL	\$10
pH Storage Solution (500 mL)	PH-SS	\$20
Photogates		
Cart Picket Fence	PF-CART	\$6
Go Direct Photogate Timing Cable	VPG-CB-GDX ⚡	\$8
Go Direct Time of Flight Pad Cable	TOF-CB-GDX ⚡	\$8
Laser Pointer	LASER	\$19
Laser Pointer Stand	STAND	\$14
Photogate Bar Tape Kit	TAPE-VPG ⚡	\$17
Picket Fence	PF	\$9


† ISE modules have a life expectancy of 1 to 2 years. We recommend that you do not purchase ISE replacement modules too far in advance of their expected time of use; degradation occurs while replacement modules are stored on the shelf.

Pulley Bracket	B-SPA	\$12
Ultra Pulley Attachment	SPA	\$24
Polarimeters (Chemical)		
Polarimeter Sample Cells (pkg. of 4)	CELLS-POL	\$66
Power Amplifier		
Accessory Speaker	PAAS-PAMP 	\$125
Projectile Launchers		
Goggles (set of 2)	GGL-VPL	\$6
Time of Flight Pad	TOF-VPL	\$84
Steel Balls (set of 6)	STB-VPL	\$7
Projectile Stop	PS-VPL	\$40
Independence of Motion Accessory	IOM-VPL	\$59
Wax Tape (300 ft.)	WXT-VPL	\$18
Rotary Motion Sensors		
Rotational Motion Accessory Kit	AK-RMV 	\$112
Rotary Motion Motor Kit	MK-RMV 	\$12
Rotary Motion Sensor Replacement Pulley	RMV-PULLEY	\$5
Rotary Motion Sensor Replacement Parts Kit	RMV-RPK	\$25
Salinity Sensors		
Salinity Standard (500 mL)	SAL-ST	\$20
Spectrophotometers/Spectrometers		
Cuvette Lids (pkg. of 100)	CUV-LID	\$9
Cuvette Rack	CUV-RACK	\$9
Plastic Cuvettes (visible) (pkg. of 100)	CUV	\$19
Plastic Cuvettes (UV-VIS) (pkg. of 100)	CUV-UV 	\$153
Quartz Cuvettes (pkg. of 2)	CUV-QUARTZ	\$199
Fluorescence/UV Quartz Cuvette (pkg. of 1)	CUV-QUARTZ-FUV	\$179
Spectrophotometer Optical Fiber (for GDX-SVISPL, VSP-UV, VSP-FUV)	VSP-FIBER	\$74
Vernier Emissions Fiber (for VSP-EM)	VSP-EM-FIBER	\$88
Spirometers		
Disposable Bacterial Filter (pkg. of 10)	SPR-FIL10	\$45
Disposable Bacterial Filter (pkg. of 30)	SPR-FIL30	\$119
Disposable Mouthpiece (pkg. of 30)	SPR-MP30	\$15
Disposable Mouthpiece (pkg. of 100)	SPR-MP100	\$36
Noseclip (pkg. of 10)	SPR-NOSE10	\$10
Noseclip (pkg. of 30)	SPR-NOSE30	\$25
O ₂ Gas Sensor to Spirometer Adapter	O2-SPR	\$8
Structures & Materials Testers		
Truss Tester Accessory	VSMT-TRUSS	\$128
Turbidity Sensor (order code TRB-BTA)		
Turbidity Accessories Replacement Kit	TRB-ACC	\$39
Turbidity Bottles (pkg. of 6)	TRB-BOT	\$29

Voltage and Current Probes

Inductor	IND	\$40
Miniature Alligator Clips for Vernier Circuit Board	VCB-GATOR	\$15
Optional Breadboard Kit for the Vernier Circuit Board 2	VCB2-OB BK	\$29
Replacement Lamps for Vernier Circuit Board	VCB-BULB	\$12
Resistivity Rods	RRS 	\$54
Vernier Circuit Board 2	VCB2 	\$129

Dynamics Cart and Track System

Part Name	Order Code	Price
For any Cart and Track System		
Adjustable Two Foot Leveler	AL-VDS	\$10
Adjustable End Stop	AS-VDS	\$8
Anti-Roll Pegs	VDS-ARPI0	\$3
Axles and Wheels for Cart	WHEELS-VDS	\$15
Cart Picket Fence	PF-CART	\$6
Cart—Plunger Cart (plastic)	DTS-CART-P	\$79
Cart—Standard Cart (plastic)	DTS-CART-S	\$68
Motion Detector Bracket	DTS-MDB	\$11
Optics Accessories	page 115	
Photogate Bracket	PGB-VDS	\$5
Pulley Bracket	B-SPA	\$12
Vernier Dynamics System Replacement Parts Kit	VDS-RPK 	\$25

For Dynamics Cart and Track Systems Only (Plastic Carts)

DFS/Accelerometer Fasteners	DTS-ACC	\$9
Eddy Current Brake	DTS-ECB	\$19
Friction Pad DTS (for plastic carts)	DTS-PAD	\$32
Mass DTS (hexagonal bars)	DTS-MASS	\$16
Motion Detector Reflector Flag	DTS-FLAG	\$9

For Vernier Dynamics Systems Only (Metal Carts)

Friction Pad (for metal carts)	PAD-VDS	\$35
Mass for Dynamics Carts (500 g block)	MASS	\$12

Go Direct

Part Name	Order Code	Price
Go Direct Charge Station	GDX-CRG	\$69
Go Direct Sensor Clamp	GDX-CLAMP	\$12
Go Direct USB Radio	GDX-RADIO	\$29
Vernier Micro USB Cable	CB-USB-MICRO	\$5
Vernier USB Type C to Micro USB Cable	CB-USB-C-MICRO	\$9

LabQuest 3, LabQuest 2, and Original LabQuest

Part Name	Order Code	Price
Vernier Mini USB Cable	CB-USB-MINI	\$5
Vernier USB Type C to Mini USB Cable	CB-USB-C-MINI	\$9

For LabQuest® 3 Only

LabQuest 3 Battery	LQ3-BAT	\$34
LabQuest 3 Lanyard	LQ3-LAN	\$9
LabQuest 3 Charging Station	LQ3-CRG	\$129
LabQuest 3 Power Supply	LQ3-PS	\$14
LabQuest 3 Stand	LQ3-STN	\$5

For LabQuest 2 and Original LabQuest Only

LabQuest Charge Station	LQ2-CRG	\$129
LabQuest Power Supply	LQ-PS	\$11
LabQuest Tether (pkg. of 5)	LQ-TETH-5	\$5
LabQuest Lanyard	LQ-LAN	\$5
LabQuest Battery Boost 3	LQ-BOOST3	\$119
LabQuest SD Card	LQ-SD	\$12

For LabQuest 2 Only

LabQuest 2 Lab Armor	LQ2-ARMOR	\$15
LabQuest 2 Stand	LQ2-STN	\$5
LabQuest 2 Battery	LQ2-BAT	\$19
LabQuest 2 Stylus (pkg. of 5)	LQ2-STYL-5	\$5

For Original LabQuest Only

Original LabQuest Battery	LQ-BAT	\$19
Original LabQuest Stylus (pkg. of 5)	LQ-STYL-5	\$5

Cables/Adapters/Power Supplies

Part Name	Order Code	Price
BTA/BTD Cables and Adapters		
Analog Bare Wire Cable	CB-BTA	\$5
Digital Bare Wire Cable	CB-BTD	\$5
Analog Breadboard Cable	BB-BTA	\$12
Digital Breadboard Cable	BB-BTD	\$11
Analog ProtoBoard Adapter	BTA-ELV	\$10
Digital ProtoBoard Adapter	BTD-ELV	\$12
Analog Sensor Extension Cable (2 m)	EXT-BTA	\$12
Digital Sensor Extension Cable (2 m)	EXT-BTD	\$12

For LabPro®

AC Adapter (for LabPro, CBL 2, or DCU)	IPS	\$12
LabPro USB Cable	CB-USB	\$5

Additional Replacement Parts Available Online
 Visit www.vernier.com/replacements

Index

A

Accelerometers
 3-Axis Accelerometer 105
 25-g Accelerometer 105
 Go Direct Acceleration 105
 Go Direct Force and Acceleration 105
 Low-g Accelerometer 105
 Accessories and replacement parts 136–137
Advanced Biology with Vernier
www.vernier.com/bio-a
Advanced Chemistry with Vernier 81
Advanced Physics with Vernier—Beyond Mechanics 118
Advanced Physics with Vernier—Mechanics 118
Agricultural Science with Vernier
www.vernier.com/awv
 Ammonium ion-selective electrodes
 Ammonium ISE www.vernier.com/nh4-bta
 Go Direct Ammonium ISE
www.vernier.com/gdx-nh4
 Anemometer www.vernier.com/anm-bta
 Arduino® products 126

B

Balances 88
 Barometer www.vernier.com/bar-bta
 BioChamber 250 www.vernier.com/bc-250
 BioChamber 2000 www.vernier.com/bc-2000
Biology with Vernier 45
 Bio-Rad® 55
 Biotechnology 54–55
 Blood pressure sensors
 Blood Pressure Sensor
www.vernier.com/bps-bta
 Go Direct Blood Pressure 50
 BNC electrodes 136
 BlueView Transilluminator 54
 Bumper and Launcher Kit 103

C

Cables 137
 Calcium ion-selective electrodes
 Calcium ISE www.vernier.com/ca-bta
 Go Direct Calcium ISE
www.vernier.com/gdx-ca

Calculators 132–133
 Calibration standards 136–137
 Canadian sales 141
 CASE 51
 CBL 2™ www.vernier.com/cbl2
 CBR 2™ 132
 Celestron® Digital Microscope Imagers 55
 Centripetal force apparatuses
 Centripetal Force Apparatus
www.vernier.com/cfa
 Go Direct Centripetal Force Apparatus 106
 Charge sensors
 Charge Sensor 108
 Go Direct Static Charge
www.vernier.com/gdx-q
 Charging stations
 Go Direct 137
 LabQuest 34
 TI-84 Plus CE 133
 TI-Nspire™ CX 133
 Chemical polarimeters 89
Chemistry with Vernier 77
 Chloride ion-selective electrodes
 Chloride ISE www.vernier.com/cl-bta
 Go Direct Chloride ISE
www.vernier.com/gdx-cl
 CO₂ gas sensors
 CO₂ Gas Sensor www.vernier.com/co2-bta
 Go Direct CO₂ Gas 20, 44
 Coding 13, 24, 128
Coding with Codey Rocky: Mission to Mars 13
Coding with mBot: Self Driving Vehicles 24
 Color Mixer Kit 115
 Colorimeters
 Colorimeter www.vernier.com/col-bta
 Go Direct Colorimeter 81
 Conductivity probes
 Conductivity Probe
www.vernier.com/con-bta
 Go Direct Conductivity 62, 83
 Platinum-Cell Conductivity Probe
www.vernier.com/conpt-bta
 Constant current systems
 Constant Current System
www.vernier.com/ccs-bta
 Go Direct Constant Current System 86
 Current sensors
 Current Probe www.vernier.com/dcp-bta
 Go Direct Current 108
 High Current Sensor
www.vernier.com/hcs-bta
 Cuvette Rack 137
 Cuvettes 137

D

Davis® weather stations 63
 Differential Voltage Probe
www.vernier.com/dvp-bta
 Diffraction Apparatus 114
 Digital curriculum 40–41
 Digital microscopes 55
 Dissolved oxygen probes
 Go Direct Optical Dissolved Oxygen 61
 Optical DO Probe
www.vernier.com/odo-bta
 Drop counters
 Drop Counter www.vernier.com/vdc-btd
 Go Direct Drop Counter 78
 Dual-Range Force Sensor 105
 Dynamics systems and accessories 100–103

E

Earth Science with Vernier 72
 EasyLink 132
 EasyTemp 132
 Eddy Current Brake 103
 EKG electrodes 136
 EKG sensors
 EKG Sensor www.vernier.com/ekg-bta
 Go Direct EKG 48
 Electrode amplifiers
 Electrode Amplifier
www.vernier.com/ea-bta
 Go Direct Electrode Amplifier
www.vernier.com/gdx-ea
 Ion-Selective Electrode Amplifier
www.vernier.com/gdx-isea
 Electrode Support 88
 Electronic lab books (e-books) 29
 Electrostatics kits 108–109
Elementary Science with Vernier 11
 ELVIS protoboard adapters
www.vernier.com/protoboard-adapters
 Emissions Spectrometer 117
 Energy sensors
 Energy Sensor www.vernier.com/ves-bta
 Go Direct Energy 66
Engineering Projects with NI LabVIEW™ and Vernier www.vernier.com/epv-e
 Equipment return 141
 Ethanol sensors
 Ethanol Sensor www.vernier.com/eth-bta
 Go Direct Ethanol Vapor
www.vernier.com/gdx-eto

Exercise heart rate monitors
 Exercise Heart Rate Monitor
www.vernier.com/ehr-bta
 Go Wireless Exercise Heart Rate
www.vernier.com/gw-ehr
Exploring Motion and Force with Go Direct Sensor Cart 22
Exploring Earth and Space Science 23
Exploring Life Science 21
Exploring Physical Science 22
 Exttech® Power Supply 109
 Extra-Long Temperature Probe
www.vernier.com/tpl-bta

F

Fan carts 103
 FLIR ONE® Thermal Cameras 111
 Flow Rate Sensor www.vernier.com/flo-bta
 Fluorescence UV/VIS Spectrophotometer 87
 Force sensors
 Dual-Range Force Sensor 105
 Force Plate 105
 Go Direct Force and Acceleration 9, 105
Forensics with Vernier www.vernier.com/fwv
 Friction Pad 103

G

Gas chromatograph 88
 Gas pressure sensors
 Gas Pressure Sensor 111
 Go Direct Gas Pressure 8, 76
 Pressure Sensor 400
www.vernier.com/ps400-bta
 Glass-Body pH Electrode BNC
www.vernier.com/gph-bnc
 GLOBE® 63
 Go Direct Charge Station 137
 Go Direct Centripetal Force Apparatus 106
 Go Direct sensors
 Go Direct 3-Axis Magnetic Field 10, 109
 Go Direct Acceleration 105
 Go Direct Ammonium Ion-Selective Electrode www.vernier.com/gdx-nh4
 Go Direct Blood Pressure 50
 Go Direct Calcium Ion-Selective Electrode
www.vernier.com/gdx-ca
 Go Direct Chloride Ion-Selective Electrode
www.vernier.com/gdx-cl
 Go Direct CO₂ Gas 20, 44
 Go Direct Colorimeter 81

H

Go Direct Conductivity 62, 83
 Go Direct Constant Current System 86
 Go Direct Current 108
 Go Direct Drop Counter 78
 Go Direct EKG 48
 Go Direct Electrode Amplifier
www.vernier.com/gdx-ea
 Go Direct Energy 66
 Go Direct Ethanol Vapor
www.vernier.com/gdx-etoh
 Go Direct Force and Acceleration 9, 105
 Go Direct Gas Pressure 8, 76
 Go Direct Glass-Body pH 85
 Go Direct Hand Dynamometer 49
 Go Direct Ion-Selective Electrode Amplifier
www.vernier.com/gdx-isea
 Go Direct Light and Color 9, 20, 93
 Go Direct Mini GC 88
 Go Direct Melt Station 89
 Go Direct Motion 8, 104
 Go Direct Nitrate Ion-Selective Electrode
www.vernier.com/gdx-no3
 Go Direct O₂ Gas 46
 Go Direct Optical Dissolved Oxygen 61
 Go Direct ORP 78
 Go Direct pH 85
 Go Direct Photogate 104
 Go Direct Polarimeter 89
 Go Direct Platinum-Cell Conductivity
www.vernier.com/gdx-conpt
 Go Direct Projectile Launcher 106
 Go Direct Radiation Monitor 117
 Go Direct Respiration Belt 50
 Go Direct Rotary Motion 107
 Go Direct Sensor Carts
 (Green and Yellow) 102
 Go Direct Sound 112
 Go Direct SpectroVis Plus
 Spectrophotometer 53, 86
 Go Direct Spirometer 50
 Go Direct Static Charge
www.vernier.com/gdx-q
 Go Direct Structures & Materials
 Tester 124
 Go Direct Surface Temperature 85
 Go Direct Tris-Compatible Flat pH 60
 Go Direct Temperature 7, 21, 85
 Go Direct Voltage 10, 108
 Go Direct Weather 63
 Go Direct Wide-Range Temperature 89
 Go Direct Sensor Clamp 62
 Go!Link www.vernier.com/go-link
 Go!Motion 104
 Go!Temp www.vernier.com/go-temp
 Goniometer www.vernier.com/gnm-bta
 Go Wireless Exercise Heart Rate
www.vernier.com/gw-ehr
 Go Wireless Heart Rate 48
 Graphical Analysis app 36
 Graphical Analysis Pro app 37
 Green Diffraction Laser 114

Hand dynamometers
 Go Direct Hand Dynamometer 49
 Hand Dynamometer
www.vernier.com/hd-bta
 Heart rate monitors
 Exercise Heart Rate Monitor
www.vernier.com/ehr-bta
 Go Wireless Exercise Heart Rate
www.vernier.com/gw-ehr
 Go Wireless Heart Rate 48
 Hand-Grip Heart Rate Monitor
www.vernier.com/hgh-bta
 High Current Sensor www.vernier.com/hcs-bta
 High-Voltage Electrostatics Kit 109
Human Physiology Experiments
www.vernier.com/hsb-hp
Human Physiology with Vernier
www.vernier.com/hp-a

I

Independence of Motion Accessory 107
 Instrumental Analysis app 88
 Instrumentation Amplifier
www.vernier.com/ina-bta
 Interfaces for LabQuest sensors
 CBL 2 www.vernier.com/cbl2
 EasyLink 132
 Go!Link www.vernier.com/go-link
 LabQuest 3 32–33
 LabQuest Mini 35
 LabQuest Stream 35
 NXT/EV3 Adapter 130
 SensorDAQ www.vernier.com/sdaq
 International sales 141
Investigating Biology through Inquiry 52
Investigating Chemistry through Inquiry 82
Investigating Environmental Science through Inquiry 61
Investigating Force 9
Investigating Gas Pressure 8
Investigating Light 9
Investigating Magnetism 10
Investigating Motion 8
Investigating Solar Energy 12
Investigating Temperature 7
Investigating Voltage 10
Investigating Wind Energy 12
 Ion-Selective Electrodes (ISE)
www.vernier.com/ise
 ISE standards 136

J

JavaScript™ 128

K

KidWind Challenge 66
 KidWind products 64–67

L

LabQuest 3 32–33
 LabQuest accessories 34
 LabQuest Mini 35
 LabQuest Stream 35
 LabQuest Viewer 34
 LEGO® MINDSTORMS® robotics 24, 130
 Light sensors
 Go Direct Light and Color 9, 20, 93
 Light Sensor 115
 Logger Pro 3 software 38

M

Magnetic field sensors
 Go Direct 3-Axis Magnetic Field 10, 109
 Magnetic Field Sensor
www.vernier.com/mg-bta
 Makeblock® products
 Codey Rocky™ 13
 Makeblock Bluetooth® Dongle 131
 mBot™ (blue and pink) 24, 131
 mBot Explorer 131
 mBot STEM Classroom Kit 131
 mBot Ranger 131
Materials Testing: Beams to Bridges with Go Direct VSMT 124
 Melt stations
 Go Direct Melt Station 89
 Melt Station 89
 Microscopes (Digital) 55
 Microphone sensors
 Go Direct Sound 112
 Microphone 112
Middle School Explorations: Chemical Reactions 21
Middle School Science with Vernier
www.vernier.com/msv
 Mini GC 88
 Mirror Set 115
 Moment of Inertia Kit 107
 Motion detectors
 CBR 2 132
 Go Direct Motion 8, 104
 Go!Motion 104
 Motion Detector 104
 Motion Encoder
 Cart and Receiver 103
 Dynamics Cart and Track Systems 100–101
 Fan Cart 103
 MyDAQ Adapter www.vernier.com/bt-mdaq

N

Nitrate ion-selective electrodes
 Go Direct Nitrate ISE
www.vernier.com/gdx-no3
 Nitrate ISE www.vernier.com/no3-bta
 NXT/EV3 Adapter 130

O

O₂ gas sensors
 Go Direct O₂ Gas 46
 O₂ Gas Sensor www.vernier.com/o2-bta
 OHAUS® balances 88
 OpenSciEd 19–21
 Optical DO probes
 Go Direct Optical Dissolved Oxygen 61
 Optical DO Probe
www.vernier.com/odo-bta
 Optical fibers 137
 Optics accessories 114–115
Organic Chemistry with Vernier 89
 ORP sensors
 Go Direct ORP 78
 ORP Sensor www.vernier.com/orp-bta

P

Packages www.vernier.com/packages
 Elementary packages 11–12
 Middle school packages 22–23, 25
 High school packages 45, 47, 49, 61, 84, 120
 PAR Sensor www.vernier.com/par-bta
 pH Buffer Capsules 136
 pH sensors
 Glass-Body pH Electrode BNC
www.vernier.com/gph-bnc
 Go Direct Glass-Body pH 85
 Go Direct pH 85
 Go Direct Tris-Compatible Flat pH 60
 pH Sensor www.vernier.com/ph-bta
 Tris-Compatible Flat pH Sensor
www.vernier.com/fph-bta
 pH Storage Solution 136
 Photogates
 Go Direct Photogate 104
 Photogate 104
Physical Science with Vernier 93
Physics Explorations and Projects 118
Physics with Vernier 118
Physics with Video Analysis www.vernier.com/pva
 Picket Fence 104
 Pivot Interactives 41
 Platinum-Cell Conductivity Probe
www.vernier.com/conpt-bta
 PLTW 49, 124, 128

Polarimeters (Chemical)
 Go Direct Polarimeter 89
 Polarimeter (Chemical) 89
 Polarizer/Analyzer Set 115
 Potassium ion-selective electrodes
 Go Direct Potassium ISE
www.vernier.com/gdx-k
 Potassium ISE www.vernier.com/k-bta
 Power Amplifier 109
 Power Amplifier Accessory Speaker 112
 Power (AC) adapters 137
 Pressure sensors
 Go Direct Gas Pressure 8, 76
 Gas Pressure Sensor
www.vernier.com/gps-bta
 Pressure Sensor 400
www.vernier.com/ps400-bta
 Primary Productivity Kit www.vernier.com/ppk
 Professional development
www.vernier.com/training
 Projectile launchers
 Go Direct Projectile Launcher 106
 Projectile Launcher www.vernier.com/vpl
 Prop 65 (California) 142–143
 ProScope kits www.vernier.com/proscope
 Protoboard adapters
www.vernier.com/protoboard-adapters
 Pyranometer www.vernier.com/pyr-bta
 Python® 128

Q

Qubit Systems sensors www.vernier.com/qubit

R

Radiation monitors
 Go Direct Radiation Monitor 117
 Vernier Radiation Monitor 117
Real-World Math with Vernier
www.vernier.com/rwv
 Reflex Hammer Accessory Kit 50
 Relative Humidity Sensor
www.vernier.com/rh-bta
 Renewable energy products 64–67
Renewable Energy with Vernier 66
 Respiration monitors
 Go Direct Respiration Belt 50
 Respiration Monitor Belt
www.vernier.com/rmb
 Returns 141
 Robotics 13, 24, 129–131

Rotary motion sensors
 Go Direct Rotary Motion 107
 Rotary Motion Sensor
www.vernier.com/rmv-btd
 Rotary Motion Motor Kit
www.vernier.com/mk-rmv
 Rotational Motion Accessory Kit 107

S

Salinity Sensor www.vernier.com/sal-bta
 SAM Labs 24
 Scratch 13, 24
 Sensor carts 102
 Site license policy 141
 Software
 Graphical Analysis app 36
 Graphical Analysis Pro app 37
 Instrumental Analysis app 88
 LabQuest App 33
 LabQuest Viewer 34
 Logger Pro 3 38
 Spectral Analysis app 39
 Thermal Analysis Plus app 111
 TI-Nspire™ software
www.vernier.com/ti-software
 TI-SmartView™ www.vernier.com/ti-sv
 Video Analysis app 40
 Soil Moisture Sensor www.vernier.com/sms-bta
 Solar Energy Exploration Kit 67
Solar Energy Explorations 25
 Solar panel 65
 Solar Thermal Exploration Kit 65
 Sound level sensors
 Go Direct Sound 112
 Sound Level Sensor 112
 SparkFun® RedBoard www.vernier.com/ard-red
 Spectral Analysis app 39
 Spectrometers/Spectrophotometers
 Emissions Spectrometer 117
 Fluorescence/UV-VIS Spectrometer 87
 Go Direct SpectroVis Plus 53, 86
 Vernier Spectrometer
www.vernier.com/v-spec
 UV-VIS Spectrometer 87
 Spectrum Tube Power Supplies 117
 Spectrum tubes 117
 Spirometer accessories 137
 Spirometers
 Go Direct Spirometer 50
 Spirometer www.vernier.com/spr-bta
 Stainless Steel Temperature Probe 111
 Static Genecon 109
 Stir Station 88
 Structures & Materials Tester 124
 Surface temperature sensors
 Go Direct Surface Temperature 85
 Surface Temperature Sensor 111

T

Technical specifications
 LabQuest 3 www.vernier.com/labq3
 LabQuest Mini www.vernier.com/lq-mini
 LabQuest Stream
www.vernier.com/lq-stream
 Vernier sensors www.vernier.com/sensors
 Temperature probes
 EasyTemp 132
 Extra-Long Temperature Probe
www.vernier.com/tpl-bta
 Go!Temp www.vernier.com/go-temp
 Go Direct Surface Temperature 85
 Go Direct Temperature 7, 21, 85
 Go Direct Wide-Range Temperature 89
 Stainless Steel Temperature Probe 111
 Surface Temperature Sensor 111
 Thermocouple www.vernier.com/tca-bta
 Wide-Range Temperature Probe 89
 Texas Instruments products 132–133
 Thermocouple www.vernier.com/tca-bta
 Time of Flight Pad 107
 Track/optics bench 115
 Transilluminator 54
 Tris-Compatible pH sensors
 Go Direct Tris-Compatible Flat pH 60
 Tris-Compatible Flat pH Sensor
www.vernier.com/fph-bta
 Truss Tester Accessory 124
 Turbidity Sensor www.vernier.com/trb-bta

U

Ultra Pulley Attachment 104
 Ultraviolet light sensors
 Go Direct Light and Color 9, 20, 93
 UVA Sensor www.vernier.com/uva-bta
 UVB Sensor www.vernier.com/uvb-bta
 USB cables 137
 USB digital microscopes 55
 UV/VIS Spectrophotometer 87

V

Vernier Chemistry Investigations for Use with AP® Chemistry 79
 Vernier Circuit Board 2 109
 Optional Breadboard Kit 109
 Vernier dynamics cart and track systems 100–101
 Vernier Emissions Spectrometer 117
 Vernier Energy Sensor www.vernier.com/ves-bta
 Vernier Fluorescence/UV-VIS Spectrometer 87
 Vernier Radiation Monitor 117
 Vernier Resistor Board 136
 Vernier Spectrometer www.vernier.com/v-spec
 Vernier UV-VIS Spectrophotometer 87
 Vernier Variable Load 125
Vernier Video Analysis: Motion and Sports 118
 Video Analysis app 40
 Voltage probes
 30-Volt Voltage Probe
www.vernier.com/30v-bta
 Differential Voltage Probe
www.vernier.com/dvp-bta
 Go Direct Voltage 10, 108
 Instrumentation Amplifier
www.vernier.com/ina-bta
 Voltage Probe www.vernier.com/vp-bta

W

Warranty information 141
 Water Depth Sampler www.vernier.com/wds
 Water quality bottles www.vernier.com/wq-bot
Water Quality with Vernier 62
 Weather sensor 63
 Weather stations 63
 White paper 1
 Wide-range temperature probes
 Go Direct Wide-Range Temperature 89
 Wide-Range Temperature Probe 89
Wind Energy Explorations 25

* AP and Advanced Placement Program are registered trademarks of the College Entrance Examination Board, which was not involved in the production of and does not endorse this product.

Satisfaction Guarantee

Vernier has been selling science education software and data-collection hardware since 1981. We pride ourselves on the quality and affordability of our products and our service to our customers. If at any time you are unhappy with any of our products or service, please get in touch.

Vernier Software & Technology

13979 SW Millikan Way
Beaverton, OR 97005-2886
www.vernier.com • info@vernier.com
Toll Free: 888-VERNIER (888-837-6437)
Fax: 503-277-2440

Product Usage

Vernier products are designed for educational use. Our products are not designed nor are they recommended for any industrial, medical, or commercial process, such as life support, patient diagnosis, control of a manufacturing process, or industrial testing of any kind. We design our products with the specifications and features that educators and students need to be successful. In our effort to keep our products affordable and easy to use, we may not meet the specifications or include the features that an industrial scientist or medical professional might want.

Equipment Return

Any product that does not meet your needs may be returned within 30 days for a full refund. Equipment returned after 30 days may be subject to a restocking fee.

A Return Merchandise Authorization, available from Vernier, is required for any product return. Equipment returned for exchange or credit must be in new condition and in its original packaging.

Prices and Shipping

Prices are effective January 1, 2021 and supersede previously published prices. Prices are in US dollars and are f.o.b. shipping point. Shipping charges may vary, depending on method of shipment. Increased shipping charges for heavier or bulkier items may apply due to weight or dimensions. Applicable sales tax may be charged. Prices are for US educational institutions only and are subject to change without notice.

International Sales

All Vernier orders for use outside of the US and Canada are handled by us and the worldwide network of Vernier dealers. Contact us for more information.

Sales of Vernier products in Canada are handled by

Vernier Canada

7030 Woodbine Ave. Suite 500
Markham, Ontario L3R 6G2
Canada
www.verniercanada.ca • info@verniercanada.ca
Phone: (800) 376-4210
Local: (705) 915-3656

Preview Policy

Most Vernier products are available for a 30-day preview (or longer, if requested) to US educational institutions.

Warranties

Most Vernier-branded products carry a five-year limited warranty. During the warranty period, Vernier will repair or replace the item if there is a defect in materials or workmanship. Outside the warranty, Vernier will attempt to repair most products. The Vernier warranty covers products when used by educational institutions only. Products manufactured by anyone other than Vernier are subject to the conditions of the warranty supplied by the manufacturer.

Additional exclusions and limitations can be found at www.vernier.com/warranty

Privacy Policy

Vernier Software & Technology does not sell, lease, or loan our mailing list or portions thereof to anyone at any time. We do not store credit card information on our online store or in our accounting system. For more information on our privacy policy, see www.vernier.com/privacy-policy

If you wish to be removed from our mailing list, simply write to us at updates@vernier.com, and we will remove you immediately.

Software Licenses

We have a very generous site license policy for our software. The purchase of one copy of Logger Pro 3 or LabQuest Viewer computer software entitles you to install it on every computer in your school or, for post-secondary institutions, department. Installation to local machines over a network is allowed. Purchasers are also permitted to distribute Logger Pro 3 to their students and instructors for home use. The license is limited to a single campus if your institution has multiple campuses.

Vernier Graphical Analysis, Vernier Spectral Analysis, and Vernier Instrumental Analysis are available as free downloads from our website or distributed through the appropriate web store. Vernier Graphical Analysis Pro is available as a subscription service. Vernier Video Analysis is available as a subscription service and is distributed as a progressive web app. Video Physics and Thermal Analysis Plus are available for purchase through the App Store. Apps for iOS, iPadOS, Android, and Chrome are distributed through their respective stores. Terms and licensing are thus determined entirely by these stores.

Other Software

Software from Pivot Interactives, Texas Instruments, Davis Instruments, and Bodelin Technologies are licensed under separate agreements by their respective companies.

Trademarks

Logger Pro 3, LabQuest, LabQuest Stream, SpectroVis, SensorDAQ, Vernier and caliper design, Go Direct, Go Wireless, Go!, Go! Link, Go!Temp, Go!Motion, LabQuest Viewer, Vernier Spectral Analysis, Vernier Thermal Analysis, Vernier EasyLink, and Vernier EasyTemp are our registered trademarks. Vernier Software & Technology, www.vernier.com, BlueView, Video Physics, Vernier Graphical Analysis, Vernier Graphical Analysis Pro, Vernier Video Analysis, and Vernier Instrumental Analysis are our trademarks or trade dress.

Apple, the Apple logo, iPhone, iPadOS and macOS are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Arduino® and  are trademarks of Arduino SA.

CBL 2, CBR 2, TI Navigator, SmartView, and TI-Nspire are trademarks of Texas Instruments.

National Instruments, NI, and LabVIEW are trademarks or trade names of National Instruments Corporation.

LEGO, the LEGO logo, MINDSTORMS and the MINDSTORMS EV3 logo are trademarks and/or copyrights of the LEGO Group. ©2020 The LEGO Group. All rights reserved.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Vernier Software & Technology is under license.

mBot, Codey Rocky, and mBlock are trademarks of Makeblock.

All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.

Technical Support

We are readily available to help you with individual questions about our software and hardware—simply email support@vernier.com or call us at our toll-free number: 888-VERNIER (888-837-6437).

Our email newsletter, *The Caliper*, makes it easy to access new ideas, learn about new products, and get inspired by fellow educators. Sign up at www.vernier.com/newsletter

How to Order



www.vernier.com



888-VERNIER
(888-837-6437)



orders@vernier.com

California Proposition 65 Warning

🚩 **PROP 65**—For more information, go to [P65Warnings.ca.gov](https://www.P65Warnings.ca.gov)

Vernier Products Affected	WARNING
Adjustable End Stop	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
BioChamber 2000	⚠️ WARNING: This product can expose you to chemicals, including bisphenol A (BPA), which are known to the State of California to cause cancer, and methyl isobutyl ketone (MIBK), which are known to the State of California to cause birth defects or other reproductive harm.
BioChamber 250	⚠️ WARNING: This product can expose you to chemicals, including methyl isobutyl ketone (MIBK), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Biology Go Direct Standard Package	⚠️ WARNING: This product can expose you to chemicals, including bisphenol A (BPA), which are known to the State of California to cause cancer, and methyl isobutyl ketone (MIBK), which are known to the State of California to cause birth defects or other reproductive harm.
BlueView Transilluminator	⚠️ WARNING: This product can expose you to chemicals, including ethyl acrylate, which are known to the State of California to cause cancer.
Celestron® Digital Imager 5MP	⚠️ WARNING: Cancer and Reproductive Harm— www.P65Warnings.ca.gov
Celestron Digital Microscope Imager	⚠️ WARNING: Cancer and Reproductive Harm— www.P65Warnings.ca.gov
FLIR ONE® Gen III Camera (iOS)	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
FLIR ONE Pro Camera (iOS)	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
FLIR ONE Pro LT	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Go Direct® Melt Station	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Go Direct Photogate Timing Cable	⚠️ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Go Direct Rotary Motion Sensor	⚠️ WARNING: This product can expose you to chemicals, including chromium, which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Go Direct Time of Flight Pad Cable	⚠️ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Go To Easy Adapter	⚠️ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Human Physiology Go Direct Standard Package	⚠️ WARNING: This product can expose you to chemicals, including methyl isobutyl ketone (MIBK), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Light Sensor	⚠️ WARNING: This product can expose you to chemicals, including antimony, which are known to the State of California to cause cancer.
Melt Station	⚠️ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.

P65Warnings.ca.gov

⚡ **PROP 65**—For more information, go to [P65Warnings.ca.gov](https://www.p65warnings.ca.gov)

Vernier Products Affected	WARNING
OHAUS Scout® 120 g	⚠ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
OHAUS Scout 220 g	⚠ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
OHAUS Scout 420 g	⚠ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Photogate Bar Tape Kit	⚠ WARNING: This product can expose you to chemicals, including formaldehyde, which are known to the State of California to cause cancer.
Plastic Cuvettes (UV-VIS)	⚠ WARNING: This product can expose you to chemicals, including Di(2-ethylhexyl) phthalate (DEHP), which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Polarimeter (Chemical)	⚠ WARNING: This product can expose you to chemicals, including chromium, which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Power Amp Accessory Speaker	⚠ WARNING: This product can expose you to chemicals, including chromium, which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Resistivity Rod Set	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Rotary Motion Motor Kit	⚠ WARNING: This product can expose you to chemicals, including formaldehyde, which are known to the State of California to cause cancer.
Rotational Motion Accessory Kit	⚠ WARNING: This product can expose you to chemicals, including chromium, which are known to the State of California to cause cancer and birth defects or other reproductive harm.
Spectrum Tube Carousel Power Supply	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Spectrum Tube Single Power Supply	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Spectrum Tubes (Air, Argon, Carbon Dioxide, Hydrogen, Helium, Neon, Nitrogen)	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Vernier Circuit Board 2	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Vernier Dynamics System Replacement Parts Kit	⚠ WARNING: This product can expose you to chemicals, including nickel, which are known to the State of California to cause cancer.
Vernier Rotary Motion Sensor	⚠ WARNING: This product can expose you to chemicals, including chromium, which are known to the State of California to cause cancer and birth defects or other reproductive harm.

Education is in our company DNA.

For over four decades, the people of Vernier Software & Technology have been pioneering technologies and sharing our passion for STEM education to give teachers and students around the world more enriching and relevant classroom experiences.

Vernier technology
is used in

150

Countries

22,138

K-12 schools across the US



I taught chemistry and physics for 34 years and used Vernier technology for over 30 years. In this trying time, I enjoy helping teachers find ways to bring authentic experiences to their students while working remotely.

*Nüsret Hisim,
Tech Support*



As a Master Recycler and a member of the Vernier Green Team, it is great to know our products will inspire tomorrow's leaders in STEM with today's technology.

*Natalie Tunison,
Warehouse*



As a Quality Assurance engineer, I make sure we release software with high quality so teachers have a flawless experience in the classroom.

*Vidya Selvamani,
Software Development*




In supporting the content and curriculum needs of our customers, I strive to demystify complex concepts with consistent, clear instructions that are easy to follow.

*Chad Krieger,
Instructional Content & Curriculum*

We're Scientists,
Engineers,
and Educators.
We're Your
People.





Vernier Software & Technology

13979 SW Millikan Way
Beaverton, OR 97005-2886

888-VERNIER
(888-837-6437)
fax 503-277-2440

www.vernier.com

info@vernier.com



Recipient not at your school?
Please send updates to
updates@vernier.com

 Celebrating
40 years

