

# Go Direct<sup>®</sup> Sensor Cart Charge Station

## Order Code GDX-CART-CRG

The Go Direct Charge Station is used to charge Go Direct Sensor Carts (not included).

For additional product information, visit [www.vernier.com/gdx-cart-crg](http://www.vernier.com/gdx-cart-crg)

### What's Included

- 4-slot Charge Station
- Charge Station Cable



### Using the Product

Perform the following steps to charge your Go Direct Sensor Carts in the charge station:

1. Connect the power cable to the charge station and to an outlet.
2. Place the Go Direct Sensor Cart in the slot so that its USB power port is easily accessible.
3. Remove the charge station connector cord from its storage hole and attach the connector to the USB port on the Go Direct Sensor Cart.

### Sensor Cart Charging Indicator Lights

- The LED next to the battery icon is orange while the sensor cart is charging.
- The LED next to the battery icon is green when the sensor cart is fully charged.

### Precautions and Operational Tips

- You can charge Go Direct Sensor Carts for any amount of time, even for a few minutes between classes.
- The charge station is not water resistant and should not be exposed to a wet environment.
- Go Direct Sensor Carts can be stored in the charger with the charge station connected to AC. Go Direct Sensor Cart batteries cannot be overcharged.
- When not using a charging slot, place the connector cord in the storage hole. This protects the cord from damage.

### Troubleshooting

Try these troubleshooting tips if your charge station is not properly charging your Go Direct Sensor Carts.

- Verify the wall outlet is providing power.
- Verify the Charge Station Cable is connected to the charge station base.
- If one port is not working, try other ports.
- If the charging LED does not come on, try charging the cart with the USB cable that came with the Go Direct Sensor Cart to verify the LED comes on using the USB cable.

### Repair Information

If you have followed the troubleshooting steps and are still having trouble with your Go Direct Sensor Cart Charge Station, contact Vernier Technical Support at [support@vernier.com](mailto:support@vernier.com) or call 888-837-6437. Support specialists will work with you to determine if the unit needs to be sent in for repair. At that time, a Return Merchandise Authorization (RMA) number will be issued and instructions will be communicated on how to return the unit for repair.

### Accessories/Replacements

#### Replacement Parts

Charge Station Cable

Charge Station Cable Euro

#### Order Code

CRG-CB

CRG-CB-EURO

### Disposal Instruction

When disposing of this electronic product, do not treat it as household waste. Its disposal is subject to regulations that vary by country and region. This item should be given to an applicable collection point for the recycling of electrical and electronic equipment. By ensuring that this product is disposed of correctly, you help prevent potential negative consequences on human health or on the environment. The recycling of materials will help to conserve natural resources. For more detail information about recycling this product, contact your local city office or your disposal service.



The symbol, show here, indicates that this product must not be disposed of in a standard waste container.

### Warranty

Vernier warrants this product to be free from defects in materials and workmanship for a period of five years from the date of shipment to the customer. This warranty does not cover damage to the product caused by abuse or improper use.



**Vernier Software & Technology**

13979 S. W. Millikan Way • Beaverton, OR 97005-2886

Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440

info@vernier.com • www.vernier.com

Rev. 3/30/2021

Go Direct and other marks shown are our trademarks or registered trademarks in the United States.

All other marks not owned by us that appear herein are the property of their respective owners, who may or may not be affiliated with, connected to, or sponsored by us.