

Friction Pad DTS

(Order Code DTS-PAD)



The Friction Pad DTS provides adjustable friction for use with the Vernier Go Direct® Sensor Cart, Standard Cart, Plunger Cart, or Motion Encoder Cart. It allows users to apply a consistent frictional force when performing dynamics experiments. The friction level can be adjusted by turning the control knob.


For best results, use the Friction Pad DTS with the Vernier Dynamics Cart and Track System. The Friction Pad DTS is not compatible with older metal carts; use the Cart Friction Pad (order code PAD-VDS) for those models.

Assembly

Remove both collision tabs from one end of the cart. Insert the Friction Pad DTS into the exposed slots. Use the thumb screw on the Friction Pad DTS to adjust the friction.

Sample Experiments

- Use a Motion Detector to study the accelerated motion of a Vernier Dynamics Cart rolling up and down an inclined plane. Set the Friction Pad so that a small friction is applied and repeat the experiment to observe the change in motion. Repeat the experiment with larger frictional forces.
- Use a Motion Detector to study the motion of a Dynamics Cart as it is pushed across a horizontal surface. Attach the Friction Pad and repeat the experiment. How did friction change the motion?



Cart moving away from the Motion Detector and slowing down

Warranty

Vernier warrants this product to be free from defects in materials and workmanship for a period of five years from the date of shipment to the customer. This warranty does not cover damage to the product caused by abuse or improper use.



Vernier Science Education
13979 SW Millikan Way • Beaverton, OR 97005-2886
Toll Free (888) 837-6437 • (503) 277-2299 • Fax (503) 277-2440
info@vernier.com • www.vernier.com

Vernier Science Education
13979 S. W. Millikan Way • Beaverton, OR 97005-2886
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440
info@vernier.com • www.vernier.com

Rev. 10/9/2025



Printed on recycled paper