

Periodic Trends

In 1870, Dmitri Mendeleev first proposed a new way of studying and organizing the then known 63 elements. The modern form of the table has been modified and improved many times since Mendeleev's tables. Pioneers like Moseley (1913) and Seaborg (1941) have made the properties of the elements much simpler to study and understand.

In this activity, you are going to discover some of the trends of the properties that exist on the modern periodic table. To do this, you are going to graph a property to determine the trend then share your trends with the class.

OBJECTIVES

In this experiment, you will

- define some of the properties of an element
- graph the properties of the elements as a function of their atomic number
- determine the general trends found on the periodic table
- understand the power of the periodic table

MATERIALS

LabQuest
LabQuest App

Periodic Table App

PROCEDURE

1. See your teacher to determine which property you will study.
2. Use the Periodic Table Application found on your LabQuest to complete the data table.

Set up your LabQuest for data entry

3. With no sensors connected to your LabQuest, select New from the File menu.
4. Tap the Table tab to display the data table.
5. Set up the **Atomic Number** column
 - a. Tap on the X column and change the column name to "Atomic Number."
 - b. Select the Generate Values check box. Modify the Start value to **1**, the Stop value to **100**, and the Increment value to **1**.
 - c. Adjust the Displayed Precision to **0** Decimal Places.
 - d. Select OK.
6. Tap on the Y column and change the column name to "Group" then select OK.

Periodic Trends

7. Set up the properties column.
 - a. From the Table menu, select New Manual Column.
 - b. Change the column name to match the property you are studying.
 - c. Enter the appropriate units.
 - d. Adjust the Displayed Precision to match your data.
 - e. Select OK.
8. Tap the first cell in the Group column and enter the data from your data table.

Viewing the data

9. Tap the Graph tab to view the graph.
10. From the Graph menu, choose Show Graph and select Graph 1.
11. Tap on the vertical axis label and select your property data. Use this graph to investigate trends in the data related to the different periods of the periodic table.
12. Sketch a graph of your data. Be sure to label the axes.



13. Tap on the horizontal axis label and select the Group data. This will display the data as if sorted by group. For example, the Alkali metals (group 1A) will graph between 1 and 2 and the Noble gases (group 0) will graph between 18 and 19.
14. From the Graph menu, choose Graph Options, uncheck the Connect Point check box, and tap OK.
15. To investigate trends in the data based on the Groups, tap and drag a region of the graph to highlight one of the Groups you wish to study (for example for the Halogens, highlight the region between 17 and 18.)
16. From the Graph menu, choose Zoom In to view the data. (To return to the original view, from the Graph menu select either AutoScale Once or Zoom Out.)
17. Sketch several graphs that show the trend of your property based on Group.



DATA

Property _____ Units _____

| Atomic Symbol | Atomic Number | Group (coded) | Value | Atomic Symbol | Atomic Number | Group (coded) | Value |
|---------------|---------------|---------------|-------|---------------|---------------|---------------|-------|
| | 1 | 1.1 | | | 26 | 8.4 | |
| | 2 | 18.1 | | | 27 | 9.4 | |
| | 3 | 1.2 | | | 28 | 10.4 | |
| | 4 | 2.2 | | | 29 | 11.4 | |
| | 5 | 13.2 | | | 30 | 12.4 | |
| | 6 | 14.2 | | | 31 | 13.4 | |
| | 7 | 15.2 | | | 32 | 14.4 | |
| | 8 | 16.2 | | | 33 | 15.4 | |
| | 9 | 17.2 | | | 34 | 16.4 | |
| | 10 | 18.2 | | | 35 | 17.4 | |
| | 11 | 1.3 | | | 36 | 18.4 | |
| | 12 | 2.3 | | | 37 | 1.5 | |
| | 13 | 13.3 | | | 38 | 2.5 | |
| | 14 | 14.3 | | | 39 | 3.5 | |
| | 15 | 15.3 | | | 40 | 4.5 | |
| | 16 | 16.3 | | | 41 | 5.5 | |
| | 17 | 17.3 | | | 42 | 6.5 | |
| | 18 | 18.3 | | | 43 | 7.5 | |
| | 19 | 1.4 | | | 44 | 8.5 | |
| | 20 | 2.4 | | | 45 | 9.5 | |
| | 21 | 3.4 | | | 46 | 10.5 | |
| | 22 | 4.4 | | | 47 | 11.5 | |
| | 23 | 5.4 | | | 48 | 12.5 | |
| | 24 | 6.4 | | | 49 | 13.5 | |
| | 25 | 7.4 | | | 50 | 14.5 | |

DATA (continued)

Property _____ Units _____

| Atomic Symbol | Atomic Number | Group (coded) | Value | Atomic Symbol | Atomic Number | Group (coded) | Value |
|---------------|---------------|---------------|-------|---------------|---------------|---------------|-------|
| | 51 | 15.5 | | | 76 | 8.6 | |
| | 52 | 16.5 | | | 77 | 9.6 | |
| | 53 | 17.5 | | | 78 | 10.6 | |
| | 54 | 18.5 | | | 79 | 11.6 | |
| | 55 | 1.6 | | | 80 | 12.6 | |
| | 56 | 2.6 | | | 81 | 13.6 | |
| | 57 | 20.601 | | | 82 | 14.6 | |
| | 58 | 20.602 | | | 83 | 15.6 | |
| | 59 | 20.603 | | | 84 | 16.6 | |
| | 60 | 20.604 | | | 85 | 17.6 | |
| | 61 | 20.605 | | | 86 | 18.6 | |
| | 62 | 20.606 | | | 87 | 1.7 | |
| | 63 | 20.607 | | | 88 | 2.7 | |
| | 64 | 20.608 | | | 89 | 20.701 | |
| | 65 | 20.609 | | | 90 | 20.702 | |
| | 66 | 20.610 | | | 91 | 20.703 | |
| | 67 | 20.611 | | | 92 | 20.704 | |
| | 68 | 20.612 | | | 93 | 20.705 | |
| | 69 | 20.613 | | | 94 | 20.706 | |
| | 70 | 20.614 | | | 95 | 20.707 | |
| | 71 | 3.6 | | | 96 | 20.708 | |
| | 72 | 4.6 | | | 97 | 20.709 | |
| | 73 | 5.6 | | | 98 | 20.710 | |
| | 74 | 6.6 | | | 99 | 20.711 | |
| | 75 | 7.6 | | | 100 | 20.712 | |

- Periodic Table of the Elements**
- | | | | | | | | | | | | | | | | | | | |
|-----------|-----------|---|-----------|-----------|-----------|-----------|-----------|-----------|-----------|---|------------|------------|------------|------------|------------|------------|------------|----------|
| H | | | | | | | | | | | | | | | | | He | |
| 3 Li | 4 Be | <ul style="list-style-type: none"> hydrogen alkali metals alkali earth metals transition metals | | | | | | | | <ul style="list-style-type: none"> poor metals metalloids noble gases rare earth metals | | | | | | | | 10 Ne |
| 11 Na | 12 Mg | | | | | | | | | | | | | | | | | 18 Ar |
| 19 K | 20 Ca | 21 Sc | 22 Ti | 23 V | 24 Cr | 25 Mn | 26 Fe | 27 Co | 28 Ni | 29 Cu | 30 Zn | 31 Ga | 32 Ge | 33 As | 34 Se | 35 Br | 36 Kr | |
| 37 Rb | 38 Sr | 39 Y | 40 Zr | 41 Nb | 42 Mo | 43 Tc | 44 Ru | 45 Rh | 46 Pd | 47 Ag | 48 Cd | 49 In | 50 Sn | 51 Sb | 52 Te | 53 I | 54 Xe | |
| 55 Cs | 56 Ba | 57 La | 58 Ce | 59 Pr | 60 Nd | 61 Pm | 62 Sm | 63 Eu | 64 Gd | 65 Tb | 66 Dy | 67 Ho | 68 Er | 69 Tm | 70 Yb | 71 Lu | 72 Hf | |
| 73 Ta | 74 W | 75 Re | 76 Os | 77 Ir | 78 Pt | 79 Au | 80 Hg | 81 Tl | 82 Pb | 83 Bi | 84 Po | 85 At | 86 Rn | 87 Fr | 88 Ra | 89 Ac | 90 Th | |
| 91 Pa | 92 U | 93 Np | 94 Pu | 95 Am | 96 Cm | 97 Bk | 98 Cf | 99 Es | 100 Fm | 101 Md | 102 No | 103 Lr | 104 Rf | 105 Db | 106 Sg | 107 Bh | 108 Hs | |
| 109 Mt | 110 Ds | 111 Rg | 112 Cn | 113 Nh | 114 Fl | 115 Mc | 116 Lv | 117 Ts | 118 Og | 119 Uue | 120 Uub | 121 Uut | 122 Uuq | 123 Uup | 124 Uuh | 125 Uus | 126 Uuq | |

- ## QUESTIONS

- 5

Periodic Trends

2. Describe, in terms of atomic structure, how each of the general trends from Question 1 might occur.
 - a.
 - b.
 - c.
 - d.
 - e.
3. In general, what trends do you notice as you enter the transitional elements? Based on atomic structure, how might this occur?
4. In general, what trends do you notice as you enter the lanthanide and actinide series of elements? Based on atomic structure, how might this occur?

TEACHER INFORMATION

Periodic Trends

- This activity uses the LabQuest Periodic Table Application as the source of data. You can use other sources of data, including the internet. Be prepared to discuss differences in the data from different sources.
- For some situations, you may want to have the students investigate only the first 36 or 54 elements.
- You can use the Periodic Trends Summary Sheet as a way to have your students record the entire class data.
- You can copy the cards below, cut them, and distribute them as a means of assigning physical properties to investigate. You may choose to do any combination of these, some give better results than others. Since several of the ideas below are similar, you may choose to do only some of these.

PHYSICAL PROPERTIES

| | | | |
|-------------------------------|----------------------|-------------------------|------------------------|
| FIRST IONIZATION ENERGY | VALENCE ELECTRONS | DENSITY | MELTING POINT |
| BOILING POINT | ATOMIC RADIUS | COVALENT RADIUS | ATOMIC VOLUME |
| SPECIFIC HEAT | HEAT OF FUSION | HEAT OF VAPORIZATION | ELECTRO- NEGATIVITY |

ANSWERS TO QUESTIONS

Atomic Radius, Atomic Volume, and Covalent Radius

Atomic Radius is the radius of the atom measured in picometers and is the basis for the calculated **Atomic Volume** measure in cm^3/mol . **Covalent Radius** is the measurement of the part of an atom that forms a covalent bond. Covalent Radius is measured in picometers

Atomic Radius and Volume decrease as you go across a period.

Atomic Radius and Volume increase as you move down a group.

Periodic Table of the Elements

Legend:

- hydrogen
- alkali metals
- alkali earth metals
- transition metals
- poor metals
- metals
- noble gases
- rare earth metals

Electronegativity and First Ionization Energy

Electronegativity is an arbitrary scale from 0 – 4 and is used to measure how much an atom wants one more electron from another atom. **First Ionization Energy** is the amount of energy required to remove an electron from an atom.

Electronegativity and first ionization energy increase as you move across a period.

Electronegativity and first ionization energy decrease as you move down a group.

Periodic Table of the Elements

Legend:

- hydrogen
- alkali metals
- alkali earth metals
- transition metals
- poor metals
- metals
- noble gases
- rare earth metals

Heat of Fusion, Heat of Vaporization, and Specific Heat

Heat of Fusion is the amount of heat necessary to change one mole of solid into one mole of liquid at the melting point. **Heat of Vaporization** is the amount of heat necessary to change one mole of liquid into one mole of vapor at the melting point. Both Heat of Fusion and Heat of Vaporization are measured in kJ/mol . **Specific Heat** is the amount of heat necessary to raise 1 kg of a substance 1 degree Kelvin. Specific Heat is measured in $\text{J/g}\cdot\text{K}$.

Heat of fusion, heat of vaporization, and specific heat all decrease as you increase atomic number within a period.

Heat of fusion, heat of vaporization, and specific heat all decrease as you increase atomic number within a group.

Periodic Table of the Elements

Legend:


- hydrogen
- alkali metals
- alkali earth metals
- transition metals
- poor metals
- metals
- noble gases
- rare earth metals

Density

Density is a calculated value of mass divided by volume. Density is measured in g/cm^3 .

Within a period the density is greatest towards the middle of the period.

Density increases going down a group.




Melting Point and Boiling Point

Melting Point and Boiling Point are the temperatures where changes in phase occur. Melting point is the temperature where a solid changes to a liquid. The boiling point is the temperature at which a liquid changes to a gas. Melting point and boiling point are measured in degrees Kelvin.

Melting and boiling points tend to increase as you move toward the center of the period.

The melting and boiling points for metals generally decreases as you go down a group. Melting and boiling points for non-metals generally increases as you go down a group.




Valence Electrons

Valence Electrons are the electrons in the outermost electron shell. This is a count and does not have any associated units.

The number of valence electrons generally increases as you move across a period. Exceptions include the transition metals that generally have two valence electrons.

The number of valence electrons is constant for elements in the same group.



Periodic Trends Summary Sheet

Property _____

1. Define your property and how it is determined.
2. On the representation of the periodic table below, draw arrows to represent an **increasing** trend.

| Periodic Table of the Elements | | | | | | | | | | | | | | | | | | H | | | | | | | | | |
|--------------------------------|----|--|--|--|--|--|--|--|--|----|----|----|-----|-----|-----|-----|-----|-----|-----|----|----|----|----|----|----|----|----|
| 1 | | | | | | | | | | | | | | | | | | 3 | 4 | | | | | | | | |
| Li | Be | | | | | | | | | | | | | | | | | B | C | N | O | F | Ne | | | | |
| | | | | | | | | | | Na | Mg | | | | | | | Al | Si | P | S | Cl | Ar | | | | |
| 11 | 12 | | | | | | | | | K | Ca | Sc | Ti | V | Cr | Mn | Fe | Co | Ni | Cu | Zn | Ga | Ge | As | Se | Br | Kr |
| 19 | 20 | | | | | | | | | Rb | Sr | Y | Zr | Nb | Mo | Tc | Ru | Rh | Pd | Ag | Cd | In | Sn | Sb | Te | I | Xe |
| 27 | 28 | | | | | | | | | Cs | Ba | La | Hf | Ta | W | Re | Os | Ir | Pt | Au | Hg | Tl | Pb | Bi | Po | At | Rn |
| 35 | 36 | | | | | | | | | Fr | Ra | Ac | Unq | Unp | Unh | Uns | Uno | Une | Unr | | | | | | | | |
| | | | | | | | | | | | | Ce | Pr | Nd | Pm | Sm | Eu | Gd | Tb | Dy | Ho | Er | Tm | Yb | Lu | | |
| | | | | | | | | | | | | Th | Pa | U | Np | Pu | Am | Cm | Bk | Cf | Es | Fm | Md | No | Lr | | |

Property _____

1. Define your property and how it is determined.
2. On the representation of the periodic table below, draw arrows to represent an **increasing** trend.

| Periodic Table of the Elements | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------------------------------|--|--|--|--|--|--|--|--|--|----|----|----|-----|-----|-----|-----|-----|-----|-----|----|----|----|----|----|----|----|----|
| 1 | | | | | | | | | | | | | | | | | | 3 | 4 | | | | | | | | |
| H | | | | | | | | | | | | | | | | | | Li | Be | | | | | | | | |
| | | | | | | | | | | B | C | N | O | F | Ne | | | | | | | | | | | | |
| 3 | | | | | | | | | | Na | Mg | | | | | | | Al | Si | P | S | Cl | Ar | | | | |
| 11 | | | | | | | | | | K | Ca | Sc | Ti | V | Cr | Mn | Fe | Co | Ni | Cu | Zn | Ga | Ge | As | Se | Br | Kr |
| 19 | | | | | | | | | | Rb | Sr | Y | Zr | Nb | Mo | Tc | Ru | Rh | Pd | Ag | Cd | In | Sn | Sb | Te | I | Xe |
| 27 | | | | | | | | | | Cs | Ba | La | Hf | Ta | W | Re | Os | Ir | Pt | Au | Hg | Tl | Pb | Bi | Po | At | Rn |
| 35 | | | | | | | | | | Fr | Ra | Ac | Unq | Unp | Unh | Uns | Uno | Une | Unr | | | | | | | | |
| | | | | | | | | | | | | Ce | Pr | Nd | Pm | Sm | Eu | Gd | Tb | Dy | Ho | Er | Tm | Yb | Lu | | |
| | | | | | | | | | | | | Th | Pa | U | Np | Pu | Am | Cm | Bk | Cf | Es | Fm | Md | No | Lr | | |

Property _____

1. Define your property and how it is determined.
2. On the representation of the periodic table below, draw arrows to represent an **increasing** trend.

| Periodic Table of the Elements | | | | | | | | | | | | | | | | | | H | | | | | | | | | |
|--------------------------------|----|--|--|--|--|--|--|--|--|----|----|----|-----|-----|-----|-----|-----|-----|-----|----|----|----|----|----|----|----|----|
| 1 | 2 | | | | | | | | | | | | | | | | | 3 | 4 | | | | | | | | |
| H | He | | | | | | | | | | | | | | | | | Li | Be | | | | | | | | |
| | | | | | | | | | | B | C | N | O | F | Ne | | | | | | | | | | | | |
| 3 | 4 | | | | | | | | | Na | Mg | | | | | | | Al | Si | P | S | Cl | Ar | | | | |
| 11 | 12 | | | | | | | | | K | Ca | Sc | Ti | V | Cr | Mn | Fe | Co | Ni | Cu | Zn | Ga | Ge | As | Se | Br | Kr |
| 19 | 20 | | | | | | | | | Rb | Sr | Y | Zr | Nb | Mo | Tc | Ru | Rh | Pd | Ag | Cd | In | Sn | Sb | Te | I | Xe |
| 27 | 28 | | | | | | | | | Cs | Ba | La | Hf | Ta | W | Re | Os | Ir | Pt | Au | Hg | Tl | Pb | Bi | Po | At | Rn |
| 35 | 36 | | | | | | | | | Fr | Ra | Ac | Unq | Unp | Unh | Uns | Uno | Une | Unr | | | | | | | | |
| | | | | | | | | | | | | Ce | Pr | Nd | Pm | Sm | Eu | Gd | Tb | Dy | Ho | Er | Tm | Yb | Lu | | |
| | | | | | | | | | | | | Th | Pa | U | Np | Pu | Am | Cm | Bk | Cf | Es | Fm | Md | No | Lr | | |

