


---

# LabQuest Stand

(Order Code LQ-STN)




The portability of the Vernier LabQuest<sup>®</sup> allows it to be used in a variety of locations, e.g., at a workbench or at a desk. This stand allows you to set the original LabQuest at different angles, which can improve the viewing angle. The stand also lifts the LabQuest off the workspace, which reduces the possibility of damage due to spills.



Set up the stand for use:

1. Pull the arm down and bend it up with your thumb.
2. Push the end of the arm through one of the holes in the back panel. To change the desired angle, select a different hole in the back panel.



Measure. Analyze. Learn.<sup>™</sup>  
**Vernier Software & Technology**

13979 S.W. Millikan Way • Beaverton, OR 97005-2886  
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440  
info@vernier.com • www.vernier.com